

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [MEDca](#) /

› [MEDca Digital Personal Sound and Voice Amplifier Instruction Manual \(Model AVA-68V4\)](#)

MEDca AVA-68V4

MEDca Digital Personal Sound and Voice Amplifier Instruction Manual

Model: AVA-68V4

INTRODUCTION

The MEDca Digital Personal Sound and Voice Amplifier is designed to enhance sound and voice clarity, making it easier to hear in various environments. This device is suitable for individuals who require sound amplification for activities such as watching television, participating in conversations, or observing nature. Its compact design and versatile earbud options provide a comfortable and effective listening experience.

PACKAGE CONTENTS

Please verify that all the following items are included in your package:

- MEDca Digital Personal Sound and Voice Amplifier Unit
- Wired Headphones
- Single Earbud Cable
- Double Earbud Cable
- Assorted Ear Tips
- 1 AAA Battery



Image: The complete package contents, showing the amplifier, two types of earbud cables, multiple ear tips, and an AAA battery.



Image: The MEDca Digital Personal Sound and Voice Amplifier unit with earbuds attached.

SETUP

Battery Installation

1. Locate the battery compartment on the amplifier unit.
2. Open the battery compartment cover.
3. Insert the provided AAA battery, ensuring the polarity (+/-) matches the markings inside the compartment.
4. Close the battery compartment cover securely.

Earbud Connection and Selection

The amplifier supports both single (mono) and double (stereo) earbud configurations. Choose the cable that best suits your hearing needs.

1. Select either the single or double earbud cable.
2. Connect the chosen earbud cable to the amplifier's audio jack.
3. Choose the most comfortable ear tips from the assorted sizes provided and attach them to the

earbuds.

4. Place the earbuds securely in your ears.



Image: This diagram shows how the removable cable allows for either single or double earbud use.

OPERATING INSTRUCTIONS

Power On/Off and Volume Control

1. To power on the device, slide the ON/OFF switch to the 'ON' position.
2. Adjust the volume using the control dial located on the amplifier. Start with a low volume and gradually increase it to a comfortable listening level.
3. To power off, slide the ON/OFF switch to the 'OFF' position.

Pitch Control

The amplifier features a pitch control switch (H/N) to adjust the sound frequency response. 'H' typically enhances higher frequencies, while 'N' provides a more natural sound profile. Experiment with both settings to find the one that offers the clearest sound for your environment.

Usage Tips

- Position the amplifier's swivel microphone towards the sound source for optimal amplification.
- The device includes a belt clip for convenient portability, allowing hands-free use.
- The amplifier is designed to pick up ambient sounds, providing a comprehensive listening experience.



Image: Detailed view of the amplifier's controls, including the volume dial, power switch, and pitch control.

MAINTENANCE

Cleaning

- Wipe the amplifier unit with a soft, dry cloth. Do not use liquid cleaners or immerse the device in water.
- Clean ear tips regularly with a mild soap and water solution, then rinse thoroughly and allow to air dry completely before reattaching.

Battery Replacement

Replace the AAA battery when sound quality diminishes or the device stops functioning. Follow the battery installation steps outlined in the Setup section.

Storage

Store the amplifier in a cool, dry place away from direct sunlight and extreme temperatures. If storing for an extended period, remove the battery to prevent leakage.

TROUBLESHOOTING

Issue	Possible Cause	Solution
No Sound	Device is off; low battery; earbuds not connected; volume too low.	Ensure device is ON; replace battery; check earbud connection; increase volume.
Weak Sound	Low battery; ear tips not sealed properly; volume too low.	Replace battery; try different ear tip size for better seal; increase volume.
Static or Feedback	Earbuds not inserted correctly; volume too high; proximity to microphone.	Re-insert earbuds for a snug fit; lower volume; ensure microphone is not too close to earbuds.
Short Battery Life	Old or low-quality battery; continuous high volume use.	Use fresh, high-quality AAA batteries; reduce volume when possible.

SPECIFICATIONS

Feature	Detail
Model Number	AVA-68V4
Maximum Output Power	130+/dB
Maximum Gain	68 dB
Harmonic Distortion	<5%
Frequency Range	450~3500Hz
Input Noise	<30dB
Battery Type	1 x AAA (1.5V)
Battery Life	Up to 30 hours (current <4mA)
Product Dimensions	4.06 x 3.9 x 2.13 inches
Weight	3.52 ounces
Material	Plastic

EACH SET INCLUDES



ADAPTER



AAA BATTERY



Image: Information regarding the AAA battery and its extended battery life.

WARRANTY AND SUPPORT

MEDca is committed to customer satisfaction. If you encounter any issues with your Digital Personal Sound and Voice Amplifier, please contact MEDca customer support for assistance. While specific warranty periods may vary, MEDca strives to ensure your needs are met.

For further information or support, please visit the official MEDca store or contact their customer service department.

Related Documents - AVA-68V4

 <p>Guía de inicialización de archivo de vídeo Avigilon</p>  <p>AVA-H501-2021E, AVA-H501-2021E, AVA-H501-2021E, AVA-H501-2021E, AVA-H501-2021E, AVA-H501-2021E, AVA-H501-2021E, AVA-H501-2021E</p>	<p>Guía de inicialización del archivo de vídeo Avigilon</p> <p>Esta guía proporciona instrucciones detalladas para la configuración e inicialización del archivo de vídeo Avigilon, incluyendo la conexión de NVRs, partición de almacenamiento, y configuración de SNMP y adaptadores de red para la monitorización del estado.</p>
<p>INSTRUCCIONES DE FUNCIONAMIENTO ORIGINALES Les abstenemos y conserve estas instrucciones MANUAL DE INSTRUCCIONES ORIGINAL Lire ces avertissements et conserve ces instructions</p>  <p>Emisor térmico Emisor térmico AVA HORIZONTAL Modèles: AVA HORIZONTAL 500 AVA HORIZONTAL 1000 AVA HORIZONTAL 1500 AVA HORIZONTAL 2000</p>	<p>Manual de Usuario Calefactor Eléctrico AVA HORIZONTAL - Instrucciones y Funcionamiento</p> <p>Este manual proporciona las instrucciones originales de funcionamiento para la serie AVA HORIZONTAL de emisores térmicos. Incluye advertencias de seguridad, guía de instalación, operación detallada y consejos de mantenimiento para los modelos AVA HORIZONTAL 500, 1000, 1500 y 2000.</p>
 	<p>BOSS AVA-1204 & AVA-1404 7-Band Graphic Equalizer User Manual</p> <p>User manual for the BOSS AVA-1204 and AVA-1404 7-Band Graphic Equalizers, detailing features, installation, electrical connections, and specifications.</p>
<p>AVA BRAIN</p> <p>IMPORTANT PRODUCT INFORMATION</p>	<p>AVA Innovations AG Product Information, Terms & Conditions, and Safety Guidelines</p> <p>This document provides important product information, terms and conditions of use, regulatory compliance details, and safety instructions for AVA products, including the AVA Nano Brain. It is available in multiple languages.</p>
<p>HJM NOTICE D'INSTRUCTION D'ORIGINE À lire attentivement et à conserver</p>  <p>Radiateur à Inertie Céramique Références: AVA V WIFI Modèles: AVA V 1200 WIFI, AVA V 1500 WIFI & AVA V 2000 WIFI</p>	<p>HJM AVA V WIFI Ceramic Inertia Radiator - User Manual and Installation Guide</p> <p>Comprehensive user manual and installation guide for the HJM AVA V WIFI Ceramic Inertia Radiator. Covers safety, installation, electrical connection, app control via WiFi, thermostat functions, troubleshooting, maintenance, and compliance information.</p>
	<p>SVAKOM Ava Neo Interactive Thrusting Vibrator User Manual</p> <p>User manual for the SVAKOM Ava Neo interactive thrusting vibrator, detailing features, functions, safety, compliance, and maintenance.</p>

