

Manuals.plus /

- › RENPHO /
- › RENPHO Eyeris 1 Eye Massager Instruction Manual

## RENPHO RF-EM001

# RENPHO Eyeris 1 Eye Massager

## INSTRUCTION MANUAL

Model: RF-EM001 | Brand: RENPHO

### 1. Introduction

The RENPHO Eyeris 1 Eye Massager is designed to provide a relaxing and soothing experience for your eyes. It combines comfortable heating, oscillating pressure, and rhythmic percussion massaging to help reduce eye strain and enhance sleep quality.

Key features include:

- **Comfortable Heating:** Built-in heating pads provide a temperature between 104°F and 113°F (40°C-45°C) to soothe tired eyes.
- **Relaxing Massage:** Utilizes oscillating pressure and rhythmic percussion for comprehensive eye relaxation.
- **Portable Design:** 180° foldable for easy transport, making it ideal for office, travel, or home use.
- **Adjustable Headband:** Ensures a comfortable fit for various head sizes.

### 2. Safety Information

Please read all safety instructions carefully before using the device.

- This product is not suitable for individuals who have undergone eye operations, have retina conditions, cataracts, glaucoma, or similar eye conditions.
- The device produces normal operating noise below 45dB during use.
- If compression feels too strong, adjust to a lower compression setting.
- Ensure you are not allergic to protein leather before use.
- Always use a 5V/1A adapter (not a PD adapter) to charge the device.

### 3. Setup

#### 3.1 Unboxing and Initial Inspection

Carefully remove the eye massager and its accessories from the packaging. Inspect all components for any signs of damage.



Figure 1: RENPHO Eyeris 1 Eye Massager in its packaging.

#### 3.2 Charging the Device

Before first use, fully charge the eye massager. Connect the provided charging cable to the device's charging port and plug it into a 5V/1A power adapter. The indicator light will show charging status.

# 104°F – 113°F Heat Comfort

Alleviate eye fatigue and maximize eye cream absorption



Figure 2: Folded eye massager showing the charging port.

### 3.3 Adjusting the Headband

Place the eye massager over your eyes. Adjust the elastic headband to ensure a snug yet comfortable fit. The headband is adjustable to accommodate various head sizes.

## 4. Operating Instructions

## 4.1 Power On/Off and Mode Selection

The device features intuitive touch controls on the side.

- **Power On/Off:** Press and hold the power button ( ) to turn the device on or off.
- **Mode Change:** Briefly press the power button ( ) to cycle through the 5 built-in modes:
  1. Compression + Heat
  2. Compression Only
  3. Compression + Heat + Continuous Vibration
  4. Intermittent Vibration Only
  5. Heat Only

# Customize Your Favorite Music

Connect your smartphone and stream your own music



Figure 3: Control buttons and available modes.

## 4.2 Compression and Volume Control

- **Compression Adjustment:** Press and hold the compression button ( ) to switch between Voice Mode and Silent Mode. Briefly press the compression button to adjust compression intensity: Soft, Strong, or No Pressure.
- **Volume Control:** Press and hold the music button (🎵) to skip to the next track. Briefly press the music button to adjust the volume (Volume + / -).

## 4.3 Bluetooth Music Playback

The RENPHO Eyeris 1 allows you to connect your smartphone via Bluetooth to stream your own music during your massage session.

# Compression Relief

Combine heat and vibration massage for better relaxation



Figure 4: Bluetooth music integration.

## 4.4 Silent Mode

For a disturbance-free experience, activate the Silent Mode. This will turn off voice prompts and music.

# Disturbance-Free Silent Mode

Press and hold to select voice / silent mode



Figure 5: Activating Silent Mode.

## 5. Maintenance

### 5.1 Cleaning

The interior material of the eye massager is made of soft, skin-friendly PU leather, which is easy to clean. If makeup or other substances get on the surface, gently wipe it off with a damp cloth. Do not use harsh chemicals or abrasive cleaners.

## 5.2 Storage

Thanks to its 180° foldable design, the eye massager can be easily folded into a compact size for convenient storage. Store it in a cool, dry place away from direct sunlight and extreme temperatures.

## 6. Troubleshooting

Problem	Possible Cause	Solution
Device not turning on	Low battery	Charge the device fully using a 5V/1A adapter.
Compression feels too strong	High compression setting	Adjust to a lower compression setting by briefly pressing the compression button.
Unusual noise during operation	Normal operating noise	The product will cause normal noise <45dB during use. This is not a defect.
Bluetooth not connecting	Device not in pairing mode or phone issue	Ensure the eye massager is in Bluetooth pairing mode. Restart your phone's Bluetooth and try again.

## 7. Specifications

- **Model Number:** RF-EM001
- **Product Dimensions:** 7.87 x 5.51 x 4.33 inches
- **Item Weight:** 1.46 Pounds
- **Batteries:** 1 Lithium Ion battery (included)
- **Power Source:** Battery Powered
- **Material:** Durable ABS + Soft Skin-friendly PU
- **Manufacturer:** RENPHO

## 8. Warranty and Support

For any questions or support regarding your RENPHO Eyeris 1 Eye Massager, our dedicated team is available 24 hours a day to assist you. Please refer to the contact information provided in your product packaging or visit the official RENPHO website for further assistance.