



**B07TK7R2TP**

# Universal Tip Rod for Dental Curing Light

Instruction Manual

## PRODUCT OVERVIEW

---

This universal tip rod is designed for use with Woodpecker LED dental curing light lamps. It serves as a light guide, directing the curing light precisely to the desired area during dental procedures. Its robust construction ensures durability and efficient light transmission.



*Image: The universal tip rod, featuring a black curved body, a silver base with grooves for attachment, and a red cap at the light-emitting end.*

## SETUP AND INSTALLATION

---

Follow these steps to properly attach the tip rod to your compatible dental curing light unit:

1. **Inspect the Tip Rod:** Before first use, carefully inspect the tip rod for any visible damage, cracks, or debris. Ensure the light-emitting surface is clean and clear.
2. **Prepare the Curing Light:** Ensure your Woodpecker LED dental curing light is turned off and disconnected from its power source for safety.
3. **Attach the Tip Rod:** Align the base of the universal tip rod with the attachment port on your curing light unit. Gently insert and twist (if applicable, refer to your curing light's manual) until it is securely seated. Do not force the attachment.
4. **Verify Secure Connection:** Once attached, gently tug on the tip rod to ensure it is firmly in place and will not detach during operation.

## OPERATING INSTRUCTIONS

---

This tip rod is an accessory for dental curing lights. Its operation is directly linked to the curing light unit itself.

Always refer to the specific operating manual of your Woodpecker LED dental curing light for detailed instructions on power settings, curing times, and safety precautions.

- **Positioning:** Position the tip of the rod as close as possible to the restorative material without touching it, ensuring the light is directed accurately.
- **Eye Protection:** Always use appropriate eye protection (orange filter or safety glasses) when operating the curing light to protect against harmful blue light.
- **Curing Protocol:** Follow the manufacturer's recommendations for curing times and light intensity based on the specific dental material being used.

## MAINTENANCE AND CARE

---

Proper maintenance ensures the longevity and optimal performance of your tip rod. Always follow standard dental sterilization protocols.

- **Cleaning:** After each use, wipe the tip rod with a soft cloth dampened with a suitable disinfectant solution. Ensure no residue remains.
- **Sterilization:** The tip rod is designed for sterilization. Consult your dental curing light's manual and local regulations for appropriate sterilization methods (e.g., autoclaving). Ensure the tip rod is completely dry before sterilization.
- **Storage:** Store the tip rod in a clean, dry, and protected environment when not in use to prevent damage and contamination.
- **Inspection:** Regularly inspect the tip rod for any signs of wear, cracks, or cloudiness on the light-emitting surface, which can impair light transmission. Replace if damage is observed.

## TROUBLESHOOTING

---

If you encounter issues with the tip rod, consider the following common problems and solutions:

Problem	Possible Cause	Solution
Poor light transmission	Dirty or scratched tip surface; internal damage.	Clean the tip thoroughly. If scratched or damaged, replace the tip rod.
Tip rod loose or detaching	Improper attachment; worn attachment mechanism on curing light or tip rod.	Ensure the tip rod is fully inserted and secured. If the issue persists, inspect both the tip rod and the curing light's attachment port for wear.
Light not emitting from tip	Curing light malfunction; tip rod not properly seated.	Verify the curing light is functioning correctly with another tip (if available). Re-seat the tip rod firmly. If the problem continues, consult the curing light's manual or a qualified technician.

## SPECIFICATIONS

---

- **Product Type:** Universal Tip Rod for Dental Curing Light
- **Compatibility:** Designed for Woodpecker LED Dental Curing Light Lamps
- **Manufacturer:** Generic
- **ASIN:** B07TK7R2TP
- **Net Quantity:** 1 Count
- **First Available:** 30 January 2020

## WARRANTY AND SUPPORT

---

For information regarding warranty coverage, returns, or technical support, please refer to the purchase documentation or contact the seller directly. Keep your proof of purchase for any warranty claims.