

[Manuals.plus](#) /

> [Good-Lite](#) /

> Good-Lite Flip-Up Occluder Glasses Fun Frames User Manual

Good-Lite flip-up-occluder-glasses-fun-47e4f9c03c

Good-Lite Flip-Up Occluder Glasses Fun Frames User Manual

Model: flip-up-occluder-glasses-fun-47e4f9c03c

INTRODUCTION

The Good-Lite Flip-Up Occluder Glasses Fun Frames are designed for pediatric patients requiring vision therapy. These glasses feature expansive circular lenses and a flip-up mechanism for simple occlusion. Constructed from durable TPE (thermoplastic elastomers), they are flexible, virtually indestructible, and latex-free. The attached rubber 'feet' are designed to prevent slipping, ensuring a secure and comfortable fit.

SETUP

The Flip-Up Occluder Glasses are ready for use directly out of the packaging. No assembly is required. Ensure the glasses fit comfortably on the user's face, with the ear temples resting securely behind the ears.



Figure 1: Front view of the Flip-Up Occluder Glasses, showing the red and yellow frame with clear lenses.

The flexible frame material allows for minor adjustments to ensure a snug fit without causing discomfort. The rubber 'feet' on the temple tips are designed to enhance stability and prevent the glasses from sliding down.

OPERATING INSTRUCTIONS

These glasses are designed for easy occlusion of either eye as needed for vision therapy.

1. **To Occlude an Eye:** Gently push the desired lens frame upwards until it clicks into its flipped-up position. This action covers the eye, providing occlusion.
2. **To Unocclude an Eye:** Gently pull the flipped-up lens frame downwards until it clicks back into its original position, revealing the eye.



Figure 2: Side view of the glasses, illustrating the design that allows lenses to flip up for occlusion.

The flip-up mechanism is designed for durability and ease of use, allowing for quick transitions between occluded

and unoccluded states.

MAINTENANCE AND CARE

Proper care will extend the life and maintain the clarity of your Flip-Up Occluder Glasses.

- **Cleaning:** Clean the glasses with a damp cloth and mild soap. Gently wipe the lenses and frame.
- **Avoid Harsh Chemicals:** Do not use harsh chemicals, abrasive materials, or alcohol-based cleaners, as these can damage the lens coating or frame material.
- **Storage:** Store the glasses in a cool, dry place when not in use. Avoid leaving them in direct sunlight or extreme temperatures, which could affect the material's integrity.
- **Handling:** While the glasses are designed to be robust, avoid excessive force when flipping the occluders or bending the frame to prevent potential damage.



Figure 3: Back view of the glasses, highlighting the rubber 'feet' on the temple tips for secure fit.

TROUBLESHOOTING

If you encounter any issues with your Flip-Up Occluder Glasses, consider the following:

- **Glasses slipping:** Ensure the glasses are properly positioned on the nose bridge and that the temple tips with rubber 'feet' are securely behind the ears. The flexible frame can be gently adjusted for a better fit.
- **Occluder not flipping smoothly:** Check for any debris or obstruction in the hinge mechanism. Clean gently with a damp cloth. Avoid forcing the mechanism.
- **Lenses appear cloudy or scratched:** Refer to the cleaning instructions in the Maintenance section. If scratches are significant and affect vision, replacement may be necessary.

PRODUCT SPECIFICATIONS

Brand	Good-Lite
Model Number	flip-up-occluder-glasses-fun-47e4f9c03c
Department	Children
Age Range Description	Kid
Item Dimensions (L x W x H)	5.25 x 4.81 x 1.06 inches
Material	TPE (thermoplastic elastomers)
Color	Red, Yellow
Lens Material Type	Plastic
Lens Coating Description	Scratch-Resistant Coating
Frame Material Type	Plastic
Frame Type	Full Rim
Item Shape	Oval
UPC	632016061258

WARRANTY INFORMATION

Specific warranty details for the Good-Lite Flip-Up Occluder Glasses Fun Frames are not provided in this manual. For information regarding product warranty, please contact Good-Lite customer support directly.

CUSTOMER SUPPORT

For further assistance, questions, or to report any issues with your Good-Lite Flip-Up Occluder Glasses, please contact the manufacturer:

- **Manufacturer:** Good-Lite
- **Website:** Please refer to the product packaging or the official Good-Lite website for contact information.