

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [OasisSpace](#) /
- › [OasisSpace Heavy Duty Bariatric Rollator Walker User Manual](#)

OasisSpace B07WR8Y6F4

OasisSpace Heavy Duty Bariatric Rollator Walker User Manual

1. INTRODUCTION

This manual provides essential information for the safe and effective use of your OasisSpace Heavy Duty Bariatric Rollator Walker. Designed for individuals requiring robust mobility support, this rollator combines durability, comfort, and ease of use. Please read these instructions thoroughly before assembly and operation to ensure proper function and user safety.



Figure 1: The OasisSpace Heavy Duty Bariatric Rollator Walker, featuring a blue frame, black padded seat, and backrest.

2. SETUP AND ASSEMBLY

The OasisSpace Rollator Walker is designed for quick and tool-free assembly, typically taking around 12 minutes. Follow these steps to set up your walker:

1. **Unpack Components:** Carefully remove all parts from the packaging. Ensure all components are present and undamaged.
2. **Unfold the Frame:** Gently pull the sides of the walker frame apart until it locks into its open position.
3. **Install Wheels:** Attach the front and rear wheels to their designated slots on the frame. Ensure they click securely into place.
4. **Adjust Handle Height:** Locate the adjustment pins on the handle tubes. Press the pins and slide the handles up or down to your desired height (adjustable from 37 to 41 inches). Ensure both handles are at

the same height and the pins are fully engaged.

5. **Adjust Seat Height:** The rollator legs are adjustable, allowing the seat height to be set between 20.5 and 23.6 inches. Adjust the leg height using the pins to achieve a comfortable seating position where your feet can touch the floor.
6. **Secure Backrest:** Ensure the padded backrest is securely attached to the frame.
7. **Attach Storage Bag:** Place the large storage bag under the seat, securing it to the frame as intended.

Two Adjustment Modes.

Convenient for life



Figure 2: Illustration of the handle and seat adjustment mechanisms for personalized fit.

Setup Video Guide

Video 1: A detailed guide demonstrating the assembly and initial setup of the OasisSpace Heavy Duty Rollator Walker.

3. OPERATING INSTRUCTIONS

Proper operation ensures safety and maximizes the benefits of your rollator walker.

- **Walking:** Push the rollator forward, maintaining a comfortable distance. Keep your body upright and centered between the handles for optimal balance.
- **Using Brakes:** Squeeze the handbrakes to slow down or stop. For temporary parking, push the brake levers down until they lock into place. To release, pull the levers up. Always engage the brakes when

sitting or standing.

- **Sitting:** Ensure the rollator is on a flat, stable surface. Engage both handbrakes by pushing the levers down. Turn around and slowly lower yourself onto the padded seat, using the handles for support. The padded backrest provides comfort.
- **Folding for Storage/Transport:** To fold the rollator, lift the padded seat upwards. The frame will collapse inward. Secure it in the folded position if a latch is present. This compact design facilitates easy storage in a car trunk or small spaces.
- **Utilizing Storage:** The large storage bag under the seat is ideal for carrying personal items, groceries, or other necessities.



Figure 3: Demonstrating proper posture and use of the rollator walker during outdoor activity.

Perfect & secure Seat

Support up to 450lbs



Figure 4: The rollator's padded seat and backrest provide a secure and comfortable resting spot.

Operation Video Guide

Video 2: Demonstrates the ease of use and maneuverability of the rollator in various indoor and outdoor settings.

4. KEY FEATURES

- **High Weight Capacity:** Supports users up to 500 lbs, thanks to its heavy-duty aluminum frame.
- **Adjustable Height:** Handles adjust from 37 to 41 inches, and seat height is adjustable from 20.5 to 23.6 inches, accommodating a wide range of user heights.
- **Advanced Braking System:** Features adjustable locking handbrakes for immediate stopping power and secure parking.
- **Enhanced Maneuverability:** Equipped with 8-inch wheels that provide smooth navigation over various terrains, including rough surfaces, grass, and boardwalks. The 25.8-inch widest distance allows passage through most standard doorways.
- **Comfortable Seating:** A large (15" x 18"), soft padded seat with memory foam and a thick (1.5-inch) breathable backrest offer superior comfort during rests.
- **Convenient Storage:** Includes a spacious storage bag located discreetly under the seat for personal

belongings.

- **Foldable Design:** Easily folds for compact storage and transport, making it convenient for travel.



Figure 5: Close-up of the thickened aluminum frame, emphasizing its robust construction for heavy-duty use.



Figure 6: Dimensions of the breathable padded seat and wide backrest for user comfort.

Fits For Large People.

Not for standard door



Figure 7: The rollator's dimensions, highlighting its suitability for larger individuals and its foldable design for transport.

5. MAINTENANCE

Regular maintenance ensures the longevity and safe operation of your rollator walker.

- **Cleaning:** Wipe down the frame and seat regularly with a damp cloth and mild soap. Avoid abrasive cleaners.
- **Brake Inspection:** Periodically check the brake cables and levers for proper function. Ensure brakes engage smoothly and hold securely. Adjust if necessary according to manufacturer guidelines.
- **Wheel Check:** Inspect wheels for wear and tear. Ensure they rotate freely and are free from debris.
- **Fasteners:** Check all screws, bolts, and fasteners regularly to ensure they are tight. Do not overtighten.

6. TROUBLESHOOTING

If you encounter any issues with your rollator, refer to the following common solutions:

- **Brakes not engaging/disengaging properly:** Check the tension of the brake cables. They may need adjustment. Ensure no debris is obstructing the brake mechanism.
- **Wheels not rolling smoothly:** Inspect wheels for hair, lint, or other debris. Clean the wheel axles if necessary. Ensure wheels are securely attached.

- **Rollator feels unstable:** Verify that all adjustment pins (for handles and legs) are fully engaged. Ensure the frame is completely unfolded and locked into position.
- **Difficulty folding/unfolding:** Check for any obstructions in the folding mechanism. Ensure the seat is fully lifted before attempting to fold.

If problems persist, contact OasisSpace customer support for assistance.

7. SPECIFICATIONS

Feature	Specification
Brand	OasisSpace
Model	B07WR8Y6F4
Material	Aluminum
Color	Blue
Item Weight	23 Pounds
Weight Capacity	Up to 500 lbs
Product Dimensions (Unfolded)	9 x 20.6 x 31 inches
Folding Size	37" x 12" x 25.8"
Handle Height Adjustment	37" to 41" (11 levels)
Seat Height Adjustment	20.5" to 23.6"
Seat Area	15" x 18"
Wheel Size	8 inches (4 wheels)
Age Range	Senior
First Available	August 19, 2019

8. SAFETY WARNINGS

Legal Disclaimer: Statements regarding medical devices have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. Always consult with a healthcare professional before using any mobility aid.

- Always ensure all locking mechanisms are engaged before use, especially when sitting.
- Do not exceed the maximum weight capacity of 500 lbs.
- Use caution on uneven or sloped surfaces.
- Regularly inspect the rollator for any signs of wear, damage, or loose parts. Do not use if damaged.
- Keep hands and fingers clear of moving parts during folding and unfolding.

9. WARRANTY AND SUPPORT

For specific warranty information, product support, or to order replacement parts, please refer to the documentation included with your purchase or visit the official OasisSpace website. Contact information for customer service is typically provided in the warranty card or user manual.