

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [OrtoPrime](#) /

› [OrtoPrime Padded Portable Adjustable Commode Chair with Clip-on Bidet - Instruction Manual](#)

OrtoPrime 11251

OrtoPrime Padded Portable Adjustable Commode Chair - Instruction Manual

Model: 11251 | Brand: OrtoPrime

1. INTRODUCTION

This manual provides essential information for the safe and effective use of your OrtoPrime Padded Portable Adjustable Commode Chair. Designed for individuals with limited mobility, the elderly, or those requiring assistance with toileting, this chair offers comfort, stability, and hygiene. Please read all instructions carefully before assembly and use to ensure proper function and user safety.

2. SAFETY INFORMATION

- **Weight Capacity:** Do not exceed the maximum weight capacity of 100 kg (220 lbs).
- **Stability:** Ensure all four legs are firmly on a level surface before use. Verify that all locking mechanisms are engaged after height adjustment.
- **Assembly:** Assemble the chair according to the instructions. Do not use if any parts are missing or damaged.
- **Maintenance:** Regularly inspect the chair for wear, loose parts, or damage. Discontinue use if any issues are found and contact customer support.
- **Supervision:** Individuals with severe mobility limitations or cognitive impairments should be supervised during use.

3. PACKAGE CONTENTS

Verify that all the following components are present in your package:

- Commode Chair Frame with Padded Backrest and Armrests
- Padded Seat with Toilet Opening
- Removable Commode Bucket with Lid
- Clip-on Bidet Attachment
- Instruction Manual

4. SETUP AND ASSEMBLY

Follow these steps to assemble and prepare your commode chair for use:

1. **Unpack Components:** Carefully remove all parts from the packaging and place them on a clean, flat surface.
2. **Assemble Frame:** The main frame with the backrest and armrests is largely pre-assembled. Ensure all joints are secure.
3. **Adjust Height:** The chair features 7 height positions. Press the spring-loaded pins on each leg and slide the inner leg up or down to the desired height. Ensure the pins snap securely into place at the same height on all four legs for stability.
4. **Attach Padded Seat:** Place the padded seat onto the chair frame. Ensure it sits securely and is properly aligned.
5. **Insert Commode Bucket:** Slide the commode bucket into the designated rails beneath the seat opening. Ensure it is fully inserted and stable.
6. **Attach Bidet (Optional):** If using the bidet function, remove the padded seat. The clip-on bidet attachment can be secured to the toilet opening on the chair frame.



Figure 1: Diagram showing the dimensions and adjustable height range of the commode chair. The chair is adjustable from 43 cm to 59 cm in height, offering 7 height positions for user comfort and accessibility.

QUALITY MATERIALS



Figure 2: Close-up of the commode chair highlighting padded armrests and the removable seat, which can be lifted for cleaning or to access the bidet function.

EASY INSTALLATION

REMOVABLE BIDET INCLUDED



Figure 3: This image illustrates the removable commode bucket and the clip-on bidet feature, demonstrating easy installation and removal for hygiene.

5. OPERATING INSTRUCTIONS

The OrtoPrime Commode Chair can be used in two primary ways:

5.1. As a Standalone Commode Chair

- Ensure the commode bucket is correctly inserted beneath the seat.
- The padded seat can be used for comfort when not in use as a commode, allowing it to function as a regular chair.
- When needed, lift the padded seat to expose the toilet opening.
- After use, replace the padded seat and remove the commode bucket for emptying and cleaning.



Figure 4: The commode chair with its padded lid closed, appearing as a regular chair, suitable for bedroom use, offering discretion and elegance.

5.2. Over a Standard Toilet

- Adjust the chair's height so it comfortably fits over your existing toilet bowl. Ensure the chair's legs are stable and level.
- Remove the commode bucket. The chair will now function as a raised toilet seat with armrests for support.
- The clip-on bidet attachment can be used in this configuration for enhanced hygiene.



Figure 5: The commode chair positioned over a standard toilet, demonstrating its versatility and providing additional support.



Figure 6: A person comfortably seated on the commode chair, which is placed over a toilet, illustrating practical use and support.

6. MAINTENANCE AND CLEANING

Regular cleaning and maintenance will prolong the life of your commode chair and ensure hygiene.

- **Commode Bucket:** Empty and clean the commode bucket after each use. Use warm water and a mild disinfectant. Rinse thoroughly and dry before reinserting.
- **Chair Frame:** Wipe the metal frame with a damp cloth and mild soap. Dry immediately to prevent rust.
- **Padded Parts:** The padded seat, backrest, and armrests can be wiped clean with a damp cloth and mild, non-abrasive cleaner. Avoid harsh chemicals that may damage the material.
- **Bidet Attachment:** Clean the bidet attachment with warm water and a mild disinfectant. Ensure all nozzles are clear of debris.
- **Inspection:** Periodically check all screws, bolts, and pins to ensure they are tight and secure. Inspect the rubber tips on the legs for wear and replace if necessary.

7. TROUBLESHOOTING

- **Chair feels unstable:** Ensure all four legs are adjusted to the same height and the spring-loaded pins

are fully engaged. Check that the chair is on a flat, level surface.

- **Difficulty inserting/removing commode bucket:** Ensure the bucket is aligned correctly with the rails. Check for any obstructions.
- **Padded seat does not sit securely:** Verify that the seat is correctly aligned and placed firmly onto the chair frame.

8. SPECIFICATIONS

Feature	Specification
Model Number	11251
Brand	OrtoPrime
Product Dimensions	45D x 59H centimeters (45 x 50 x 59 cm)
Item Weight	8.81 Kilograms
Material	Foam (Padded parts)
Color	Blue
Minimum Height	43 Centimeters
Seat Depth	45 Centimeters
GTIN	08437019534042

9. WARRANTY AND SUPPORT

For warranty information, technical support, or any questions regarding your OrtoPrime Padded Portable Adjustable Commode Chair, please contact OrtoPrime customer service. Refer to the product packaging or the official OrtoPrime website for contact details.