

[Manuals.plus](#) /

› [Serta](#) /

› Serta Motion Essentials 1st Gen / Perfect 1.0 Adjustable Bed Replacement Power Cord Kit Instruction Manual

Serta Motion Essentials 1st Gen / Perfect 1.0

Serta Motion Essentials 1st Gen / Perfect 1.0 Adjustable Bed Replacement Power Cord Kit Instruction Manual

Model: Motion Essentials 1st Gen / Perfect 1.0

INTRODUCTION

This manual provides instructions for the installation and use of the Serta Replacement Power Cord Kit for Motion Essentials 1st Generation and Perfect 1.0 Adjustable Beds. This kit is designed to restore power to your adjustable bed base. Please read these instructions carefully before proceeding with installation.

Important Compatibility Note: This power cord kit is specifically designed for Serta Motion Essentials (1st Gen) and Perfect 1.0 adjustable bed models. It is crucial to verify that your bed base and remote control match the specifications and images provided in this manual to ensure compatibility. This kit does not include a remote control.

SAFETY INFORMATION

- Always unplug the bed from the wall outlet before performing any maintenance or installation.
- Do not operate the bed with a damaged power cord or plug.
- Keep all cords away from heated surfaces.
- Ensure the power cord is not pinched or routed in a way that could cause damage.
- This product is for indoor use only.

PACKAGE CONTENTS

Verify that your package contains the following items:

- 1x AC Adapter Box/Power Supply
- 1x Power Input Cord (connects AC adapter to wall outlet)
- 1x Output Cord (connects AC adapter to bed base control box)



Image: The complete Serta adjustable bed replacement power cord kit, including the AC adapter box, the power input cord for the wall outlet, and the output cord for the bed base.

SETUP AND INSTALLATION

Follow these steps to install the replacement power cord kit:

1. **Unplug Existing Power:** Ensure the adjustable bed is completely unplugged from the wall outlet before beginning.
2. **Locate Control Box:** Find the control box on your adjustable bed base. This is typically located underneath the bed frame.
3. **Disconnect Old Cord (if applicable):** If replacing an existing cord, carefully disconnect the old power input cord and the output cord from the control box and the AC adapter.
4. **Connect New Output Cord:** Connect the output cord from the new AC adapter box to the appropriate port on your bed's control box. Ensure a secure connection.
5. **Connect New Input Cord:** Connect the power input cord to the AC adapter box.
6. **Plug into Wall Outlet:** Once all connections to the bed base and AC adapter are secure, plug the power input cord into a grounded electrical wall outlet.
7. **Test Functionality:** Use your bed's remote control to test the adjustable functions (head, foot, etc.) to ensure power has been restored and the bed operates correctly.

Remote Control Compatibility Reference: It is essential that your existing remote control matches one of the designs shown below for this power kit to be compatible with your bed base. If your remote does not match, this kit may not work.



Image: A silver Serta adjustable bed remote control with multiple buttons for head, foot, massage, and preset positions. This remote is indicative of compatible 1st Gen models.



Image: A white, wired Serta Motion Essentials I remote control with basic head and foot adjustment buttons. This remote is also indicative of compatible 1st Gen models.

OPERATING INSTRUCTIONS

Once the power cord kit is correctly installed and connected, your adjustable bed should function normally using its original remote control. Refer to your bed's original user manual for specific operating instructions related to its features (e.g., head/foot articulation, massage functions, preset positions).

MAINTENANCE

- **Cleaning:** Wipe the power cord and adapter box with a dry, soft cloth. Do not use liquid cleaners.
- **Inspection:** Periodically inspect the power cord for any signs of damage, fraying, or kinks. If damage is found, discontinue use and replace the cord.
- **Storage:** If the bed is to be stored or moved, unplug the power cord from the wall outlet and carefully coil it to prevent damage.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Bed does not respond to remote.	No power to the bed base. Loose connection. Incorrect power cord kit for your model.	Ensure the power cord is securely plugged into both the wall outlet and the AC adapter box, and the adapter box is connected to the bed's control box. Check if the wall outlet is functional by plugging in another device. Verify that your bed model and remote match the compatibility information in this manual.
Intermittent power or flickering lights on control box.	Loose connection. Damaged cord or adapter.	Check all connections for tightness. Inspect the power cord and AC adapter for visible damage. If damaged, replacement may be necessary.

SPECIFICATIONS

- **Brand:** Serta
- **Model Compatibility:** Serta Motion Essentials 1st Gen, Serta Perfect 1.0
- **Connector Type:** Male-to-Female (for power input cord)
- **Voltage:** 110-120 Volts AC (Input)
- **Plug Format:** Type A (US Standard)

WARRANTY AND SUPPORT

For warranty information or technical support regarding this replacement power cord kit, please refer to the documentation provided with your original Serta adjustable bed or contact the seller directly. Keep your purchase receipt as proof of purchase.