

Dino-Lite AF4515ZT

Dino-Lite AF4515ZT USB Handheld Digital Microscope User Manual

Model: AF4515ZT

1. INTRODUCTION

This manual provides comprehensive instructions for the setup, operation, and maintenance of your Dino-Lite AF4515ZT USB Handheld Digital Microscope. The AF4515ZT is a versatile digital microscope designed for various applications, offering optical magnification from 20x to 220x, a 1.3-megapixel true image resolution, and advanced features like polarized light and Automatic Magnification Reading (AMR).

Please read this manual thoroughly before using the device to ensure proper functionality and to maximize your experience.





Figure 1: Dino-Lite AF4515ZT USB Handheld Digital Microscope

2. PRODUCT OVERVIEW

The Dino-Lite AF4515ZT features a robust design with intuitive controls for precise observation and analysis. Understanding its components is essential for effective use.




Figure 2: Labeled diagram of the Dino-Lite AF4515ZT components. This image illustrates the various parts of the microscope, including the USB 2.0 cable for connectivity, the MicroTouch sensor for image capture and LED control, the scroll lock to fix magnification, the magnification dial for adjustment, the magnification power indicator, the adjustable polarizer to reduce glare, LED lighting with flexible control, and interchangeable caps for different viewing needs.

Key Components:

- **Connectivity:** USB 2.0 Cable for connection to a computer.
- **MicroTouch:** A touch sensor trigger for various functions, including image capture and LED controls.
- **Scroll Lock:** Locks the current magnification power to prevent accidental changes.
- **Magnification Dial:** Used to adjust the optical magnification of the microscope.
- **Magnification Power Indicator:** Displays the current magnification level in use.
- **Adjustable Polarizer:** An integrated polarizer that can be adjusted to reduce or eliminate glare from reflective surfaces.
- **LED Lighting:** Features 8 LEDs with Flexible Light Control (FLC) for optimal illumination.
- **Interchangeable Caps:** Various cap designs are available (sold separately) to aid in providing the best possible image for specific applications.



Figure 3: Dino-Lite AF4515ZT dimensions and optional WiFi module. The image shows the microscope's dimensions (approximately 4.65 inches in length) and illustrates how it can be connected to an optional WiFi module (WF-20) for wireless operation, extending its length to approximately 5.5 inches. Accessories are sold separately.

3. KEY FEATURES

The Dino-Lite AF4515ZT is equipped with several advanced features to enhance your microscopy experience:

- **Optical Magnification:** Offers a variable optical magnification range from 20x to 220x.
- **Image Resolution:** Captures images at 1280 x 1024 (1.3MP) true resolution.
- **Video Frame Rate:** Provides smooth video at 30 frames per second (FPS).
- **Measurement Capabilities:** Perform various measurements including straight lines, circles, angles, arches, and grids in different units (inch/mil/mm/um).
- **Polarization Control:** The polarization feature can be switched on/off or adjusted to full or half polarization to effectively control glare on reflective surfaces.
- **Automatic Magnification Reading (AMR):** Displays the current magnification directly on the screen within the included software for Windows or Mac OS.
- **Comprehensive Software:** The included software offers rich features for saving, annotating, stitching, comparison, and analysis of captured images.
- **Operating System Compatibility:** Supports Windows (XP/Vista/7/8/10), Mac OS 10.5 and above, Android, and iOS (requires WF-20 Wi-Fi adapter).

Product Information

Key Features



AMR

Live magnification rate displayed in software



FLC

Control LED brightness and toggle by pairs



Adjustable Polarizer

Integrated polarizer eliminates glare



Wireless Solution

Optional WF-20 module

Main Features



20x - 220x

Variable magnification



USB 2.0

UVC Drivers for maximum compatibility



Windows and Mac

Windows 7, 8, 10 & Mac OS 10.9+.



Measurement & Annotation

Lines, arcs, angles, and much more



1280 x 1024

Resolution up to 1280 x 1024 pixels



Live viewing



Image and Video capture



Barcode reader



Language options



Drawing & annotation tools



Software and future updates included



Gridlines, crosshair, scales



Software Development Kit available

Figure 4: Summary of Dino-Lite AF4515ZT key features, including AMR, FLC, Adjustable Polarizer, Wireless Solution (optional), Variable Magnification (20x-220x), USB 2.0, Windows and Mac compatibility, Measurement & Annotation tools, 1280x1024 resolution, Live viewing, Image and Video capture, Barcode reader, Language options, Drawing & annotation tools, Software and future updates, Grids/crosshairs/scales, and Software Development Kit availability.

4. SETUP

Follow these steps to set up your Dino-Lite AF4515ZT microscope:

4.1. Package Contents Verification

Before proceeding, ensure all items are present in the package:

- Dino-Lite AF4515ZT Digital Microscope
- Software CD (DinoCapture 2.0)
- Quick Guide
- Calibration Target

- Storage Pouch
- USB Cable Tie
- Interchangeable Caps (e.g., open cap, closed cap, diffuser cap)



Figure 5: Contents of the Dino-Lite AF4515ZT package, including the microscope, software CD, quick guide, calibration target, storage pouch, USB cable tie, and various interchangeable caps.

4.2. Software Installation

1. Insert the included software CD into your computer's optical drive.
2. Follow the on-screen instructions to install DinoCapture 2.0 (for Windows) or DinoXscope (for Mac OS).
3. Alternatively, the latest software version can be downloaded from the official Dino-Lite website.
4. Ensure the software is fully installed before connecting the microscope.

4.3. Hardware Connection

1. Connect the USB cable of the Dino-Lite microscope to an available USB 2.0 port on your computer.
2. The operating system should automatically detect the device and install necessary drivers.
3. Once connected, launch the DinoCapture 2.0 or DinoXscope software.

5. OPERATING INSTRUCTIONS

This section details the basic operation of your Dino-Lite microscope.

5.1. Live Viewing

1. Place the microscope over the object you wish to observe.

2. Open the DinoCapture 2.0 (Windows) or DinoXcope (Mac OS) software. The live image from the microscope should appear on your screen.

5.2. Adjusting Magnification and Focus

1. Rotate the **Magnification Dial** on the microscope body to adjust the desired magnification level (20x to 220x).
2. Observe the **Magnification Power Indicator** to see the current magnification.
3. To focus, gently move the microscope closer to or further from the object until the image on the screen is clear. For fine adjustment, use a compatible stand (sold separately).
4. If you wish to maintain a specific magnification, engage the **Scroll Lock**.

5.3. Using LED Lighting

The microscope features 8 built-in LEDs. The brightness can be adjusted through the software interface or by using the MicroTouch sensor (if configured for LED control).

5.4. Adjusting Polarization

For reflective surfaces, rotate the **Adjustable Polarizer** ring located at the front of the microscope. This will reduce or eliminate glare, improving image clarity. You can adjust it for full or half polarization.


5.5. Capturing Images and Videos

To capture an image or record a video:

- Press the **MicroTouch** sensor on the microscope body once to capture a still image.
- Alternatively, use the capture buttons within the DinoCapture/DinoXcope software interface.
- For video recording, select the video recording option in the software and click to start/stop recording.

6. SOFTWARE FEATURES (DINO CAPTURE 2.0 / DINO XSCOPE)

The included software is a powerful tool for managing and analyzing your microscopic observations.



The included software.

DinoCapture 2.0 for Windows and DinoXcope for Mac OS provide the convenience to capture, record, annotate, measure, and more with the Dino-Lite.

Figure 6: The DinoCapture 2.0 software interface, showing live view, captured images in a sidebar, and various measurement and annotation tools applied to an image. The software allows for detailed analysis and documentation.

Key Software Functions:

- **Live Viewing:** Real-time display of the microscope's feed.
- **Image & Video Capture:** Save still images and record video footage.
- **Measurement Tools:** Perform precise measurements on captured images, including lines, circles, angles, and more.
- **Annotation:** Add text, arrows, and shapes to images for documentation.
- **Image Stitching:** Combine multiple images to create a larger field of view.
- **Comparison:** Compare multiple images side-by-side.
- **Automatic Magnification Reading (AMR):** Displays the current magnification on the screen.
- **Flexible Light Control (FLC):** Adjust LED brightness and toggle by pairs.
- **Compatibility:** Supports Windows (XP/Vista/7/8/10) and Mac OS 10.5 and above.

7. MAINTENANCE

Proper maintenance ensures the longevity and optimal performance of your Dino-Lite microscope.

7.1. Cleaning

- **Microscope Body:** Wipe the exterior with a soft, dry cloth. For stubborn dirt, use a slightly damp cloth with mild soap, then dry thoroughly. Avoid harsh chemicals or abrasive cleaners.

- **Lens:** Use a lens cleaning brush or a blower to remove dust. For smudges, use a specialized lens cleaning solution and a microfiber lens cloth. Do not touch the lens with your fingers.

7.2. Storage

- Store the microscope in its original packaging or the provided storage pouch when not in use.
- Keep the device in a cool, dry place, away from direct sunlight, extreme temperatures, and high humidity.
- Avoid storing in dusty environments.

8. TROUBLESHOOTING


If you encounter issues with your Dino-Lite AF4515ZT, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
No image displayed in software.	Microscope not connected or recognized. Software not launched or configured correctly.	Ensure USB cable is securely connected. Restart software and computer. Verify drivers are installed. Check if another application is using the camera.
Image is blurry or out of focus.	Incorrect focus or magnification setting. Lens is dirty.	Adjust the magnification dial and move the microscope closer/further from the object. Clean the lens as per maintenance instructions.
Image has glare or reflections.	Reflective surface.	Adjust the Adjustable Polarizer ring on the microscope.
LEDs are not lighting up.	LEDs turned off in software or via MicroTouch. Power issue.	Check LED settings in the software. Press MicroTouch to cycle LED states. Ensure USB connection provides sufficient power.
Software crashes or freezes.	Software conflict or corruption. Outdated drivers.	Update software to the latest version. Reinstall software. Ensure your computer meets minimum system requirements.

9. SPECIFICATIONS

Detailed technical specifications for the Dino-Lite AF4515ZT microscope:

Specification	Value
Model Name	AF4515ZT
Magnification Range	20x - 220x (Optical)
Image Resolution	1.3 Megapixel (1280 x 1024)
Video Frame Rate	30 FPS



[\[pdf\]](#) User Manual

Product Dino Lite Europe AF4515ZT Wired prieš 1 dieną — Flexible LED Control FLC The is equipped with premium features such as control target user manual N3C O Open cap C product 13 dino lite eu ru component eshop id |||

Dino-Lite Europe - **AF4515ZT** - Wired Phone: - Email: : Dino-Lite : **AF4515ZT** Short Description 1.3 Megapixel Edge sensor 20-220x magnification Adjustable polarizer Flexible LED Control FLC Automatic Magnification Reading AMR The Dino-Lite **AF4515ZT** is part of Universal...

lang:en score:47 filesize: 99.32 K page_count: 4 document date: 2022-06-10



[\[pdf\]](#)

Product Dino Lite Europe product 13 dino lite eu cs af4515zt eshop id |||

Dino-Lite Europe - **AF4515ZT** Zna ka: Kd produktu: **AF4515ZT** Phone: - Email: Short Description 1.3 Megapixel Edge sensor 20x - 220x magnification Adjustable polarizer Flexible LED Control FLC Automatic Magnification Reading AMR Popis The Dino-Lite **AF4515ZT** is part of Universal series and connect...

lang:en score:44 filesize: 170.98 K page_count: 4 document date: 2025-02-18



[\[pdf\]](#)

Product Dino Lite Europe product 13 dino lite eu es af4515zt eshop id |||

Dino-Lite Europe - **AF4515ZT** Phone: - Email: Marca: Cdigo del producto: **AF4515ZT** Short Description 1.3 Megapixel Edge sensor 20x - 220x magnification Adjustable polarizer Flexible LED Control FLC Automatic Magnification Reading AMR Descripcin The Dino-Lite **AF4515ZT** is part of Universal series...

lang:en score:44 filesize: 170.99 K page_count: 4 document date: 2025-03-11



[pdf]

Product Dino Lite Europe product 13 dino lite eu pl af4515zt eshop id |||

Dino-Lite Europe - **AF4515ZT** Brand: Product Code: **AF4515ZT** Phone: - Email: Short Description 1.3 Megapixel Edge sensor 20x - 220x magnification Adjustable polarizer Flexible LED Control FLC Automatic Magnification Reading AMR Description The Dino-Lite **AF4515ZT** is part of Universal series and c...

lang:en score:44 filesize: 170.96 K page_count: 4 document date: 2025-03-11




[pdf]

Product Dino Lite Europe AF4515ZT The main features of the are Wireless ready requires WF 20 adapter optional 20x ~ 220x magnification 1 3 Megapixel Edge sensor product 13 dino lite eu en af4515zt eshop id |||

Dino-Lite Europe - **AF4515ZT** Brand: Product Code: **AF4515ZT** Phone: - Email: Short Description 1.3 Megapixel Edge sensor 20x - 220x magnification Adjustable polarizer Flexible LED Control FLC Automatic Magnification Reading AMR Description The Dino-Lite **AF4515ZT** is part of Universal series and c...

lang:en score:42 filesize: 170.96 K page_count: 4 document date: 2025-07-01




[pdf]

Product Dino Lite Europe AF4515ZT Wired Phone Email Marke Artikelnummer Short Description 1 3 Megapixel Edge sensor 20 220x magnification product 13 dino lite eu de component eshop id |||

Dino-Lite Europe - **AF4515ZT** - Wired Phone: - Email: Marke: Dino-Lite Artikelnummer: **AF4515ZT** Short Description 1.3 Megapixel Edge sensor 20-220x magnification Adjustable polarizer Flexible LED Control FLC Automatic Magnification Reading AMR Beschreibung The Dino-Lite **AF4515ZT** is part of Univ...

lang:en score:41 filesize: 99.34 K page_count: 4 document date: 2022-06-26




[pdf] Datasheet

Product Dino Lite Europe Datasheet AF4515ZT USB 2.0 Digital Microscope 1.3M Pixels 20 → 220X Magnification RS A700000010613555 docs rs online b751 |||

Dino-Lite Europe - **AF4515ZT** - Wired Phone: - Email: Brand: Dino-Lite Product Code: **AF4515ZT** Short Description 1.3 Megapixel Edge sensor 20-220x magnification Adjustable polarizer Flexible LED Control FLC Automatic Magnification Reading AMR Description The Dino-Lite **AF4515ZT** is part of Univer...

lang:en score:39 filesize: 99.35 K page_count: 4 document date: 2023-10-12



[pdf]

Product Dino Lite Europe AF4515ZT Wired 1 day ago — fully wireless solution The is equipped with premium features such as Flexible LED control FLC integrated polarizer extensive measurement 76882 Datenblatt astroshop de Produkt |||

Dino-Lite Europe - **AF4515ZT** - Wired Phone: - Email: Marke: Dino-Lite Artikelnummer: **AF4515ZT** Short Description 1.3 Megapixel Edge sensor 20-220x magnification Adjustable polarizer Flexible LED Control FLC Automatic Magnification Reading AMR Beschreibung The Dino-Lite **AF4515ZT** is part of Univ...

lang:en score:38 filesize: 99.34 K page_count: 4 document date: 2023-08-11



[pdf]

AF4515ZT El Dino Lite Edge 1.3MP es un dispositivo todo en uno diseñado para cumplir con la mayoría de aplicaciones profesionales agelectronica lat s textos A |||

AF4515ZT Product Images Descripcin El Dino-Lite Edge 1.3MP**AF4515ZT** es un dispositivo todo-en-uno diseado para cumplir con la mayora de aplicaciones profesionales con aumentos de 20x - 220x. El **AF4515ZT** integra las caractersticas premium de la serie Edge AMR y FLC. El nuevo diseo modular listo para ...

lang:es score:36 filesize: 688.88 K page_count: 4 document date: 2021-02-10



[\[pdf\]](#)

Product Dino Lite Europe product 13 dino lite eu da component eshop |||

Dino-Lite Europe - **AF4515ZT** Brand: Product Code: **AF4515ZT** Phone: - Email: Short Description 1.3 Megapixel Edge sensor 20x - 220x magnification Adjustable polarizer Flexible LED Control FLC Automatic Magnification Reading AMR Description The Dino-Lite **AF4515ZT** is part of Universal series and c...

lang:en score:35 filesize: 170.96 K page_count: 4 document date: 2024-12-24



[\[pdf\]](#) Specifications Datasheet

Untitled 1 Datasheet 2204919 Dino Lite WF4515ZT Digital microscope Transmitted light 4712805477654 datasheet dino lite wf4515zt digital transmitted asset conrad media10 add 160267 c1 en 002204919DS00

DINO-LITE WIRELESS 100 WIRELESS 100 VERSATILE 100 QUALITY WIRELESS WITHOUT COMPROMISING ON ... LED control FLC AF4115ZTL Same as AF4115ZT, but with 10-140x Magnification Long working Distance **AF4515ZT** The Dino-Lite **AF4515ZT** is part of the Universal series and connects to your device via USB....

lang:en score:35 filesize: 1.06 M page_count: 6 document date: 2020-01-29



[\[pdf\]](#)

Untitled 1 datablad 2204918 dino lite wf4115ztl digitale microscoop doorvallend licht asset conrad media10 add 160267 c1 de 002204918DS00

DINO-LITE WIRELESS 100 DRAHTLOS 100 VIELSEITIG 100 QUALITT DRAHTLOS OHNE KOMPROMISSE IN SAC ... r mit 10-140-fache Vergrößerung Längerem Arbeitsabstand DINO-LITE WIRED -SERIE AF-SERIE AF4115ZT **AF4515ZT** AF4915ZT Das Dino-Lite AF4115ZT ist Teil der Universal-Serie und wird über USB an Ihr Gerät...

lang:de score:33 filesize: 1.06 M page_count: 6 document date: 2020-01-29



[\[pdf\]](#) Specifications Flyer Borchure

Dino Lite 3 fold flyer wf20 ENG Lambda indd Wireless Brochure WF4915ZT 1 3MP 20 200x Microscope wireless brochure LAMBDA lambdaphoto co uk s dinolite |||

DINO-LITE WIRELESS 100 WIRELESS 100 VERSATILE 100 QUALITY WIRELESS WITHOUT COMPROMISING ON ... LED control FLC AF4115ZTL Same as AF4115ZT, but with 10-140x Magnification Long working Distance **AF4515ZT** The Dino-Lite **AF4515ZT** is part of the Universal series and connects to your device via USB....

lang:en score:32 filesize: 6.32 M page_count: 6 document date: 2019-07-31




[\[pdf\]](#) Specifications Datasheet

Untitled 1 Data Sheet Microscope Digital element14 Australia 2864990 farnell datasheets |||

DINO-LITE WIRELESS 100 WIRELESS 100 VERSATILE 100 QUALITY WIRELESS WITHOUT COMPROMISING ON ... LED control FLC AF4115ZTL Same as AF4115ZT, but with 10-140x Magnification Long working Distance **AF4515ZT** The Dino-Lite **AF4515ZT** is part of the Universal series and connects to your device via USB....

lang:en score:25 filesize: 6.14 M page_count: 6 document date: 2019-04-24



[\[pdf\]](#) Datasheet

Product Dino Lite Europe Data sheet EN Microscope AF4515ZT WF 20 1 3MP 220x 8 LED 30 fps USB 2.0 WiFi 76862 Datenblatt micro pro Produkt |||

Dino-Lite Europe - WF4515ZT - Wireless Phone: - Email: Marke: Dino-Lite

Artikelnummer: WF4515ZT Short Description 1.3 Megapixel Edge sensor 20-220x

magnification Adjustable polarizer Flexible LED Control FLC Automatic Magnification

Reading AMR Beschreibung The Dino-Lite WF4515ZT is a Wireles...

lang:en **score:24** filesize: 98.63 K page_count: 4 document date: 2023-08-10