

## Stander B086VTHHVH

# Stander Let's Go Out Rollator User Manual

Model: B086VTHHVH

## 1. PRODUCT OVERVIEW

The Stander Let's Go Out Rollator is a lightweight, four-wheel Euro-style walker designed to provide mobility assistance and comfort. It features a durable frame, a comfortable seat, and a reliable braking system, making it an ideal aid for daily activities and travel.

### Key Features:

- **Large Wheels and Brakes:** Equipped with 10-inch front swivel wheels and 8-inch rear wheels for smooth navigation. Features a full locking-brake system for safety and stability, and built-in coil suspension for shock absorption on various terrains.
- **Comfortable Seat:** A 17-inch wide by 7-inch seat supports up to 300 pounds, complemented by a faux-leather backrest for additional support and comfort.
- **No Assembly Required:** Ready to use directly out of the box.
- **Adjustable Handles:** Handles can be adjusted from 30.5 to 36 inches high to suit individual user preferences.
- **Folding Design:** Weighs only 17 pounds and folds quickly into an 8-inch wide footprint for convenient storage and transport.
- **Stylish Design:** Features a two-tone design with faux-leather seat and backrest, molded steel accents, and a storage bag.



The Stander Let's Go Out Rollator, showcasing its sleek design and functional components.

## 2. SETUP AND INITIAL USE

The Stander Let's Go Out Rollator is designed for immediate use with no assembly required. Simply remove it from its packaging and unfold it according to the instructions below.

### Unfolding the Rollator:

1. Place the folded rollator on a flat, stable surface.
2. Gently pull the handles outwards until the frame fully extends and locks into place. Ensure all locking mechanisms are engaged.
3. Verify that the seat is securely positioned and flat.

### Adjusting Handle Height:

The handles can be adjusted to a comfortable height for the user, ranging from 30.5 to 36 inches.

1. Locate the height adjustment pins or levers on the handle shafts.
2. Press or release the pins/levers and slide the handles up or down to the desired height.
3. Ensure the pins snap securely into the adjustment holes on both sides to prevent slippage. The handles should be at wrist level when standing upright with arms relaxed.



A user adjusting the height of the rollator handles for personalized comfort and proper posture.

## 3. OPERATING INSTRUCTIONS

## Using the Brakes:

The rollator is equipped with a reliable braking system for safety and control.

- **To Slow Down or Stop:** Squeeze the brake levers located on the handles towards you. The harder you squeeze, the more braking force will be applied.
- **To Engage Parking Brakes:** Push the brake levers down until they click into the locked position. This will prevent the rollator from moving when you are sitting or standing still.
- **To Release Parking Brakes:** Pull the brake levers back up until they disengage from the locked position.

## Using the Seat:

The integrated seat provides a convenient resting spot.

- Ensure the rollator is on a flat, stable surface.
- Engage both parking brakes before attempting to sit down.
- Carefully turn around and sit on the seat, using the handles for support if needed.
- Always disengage the parking brakes before attempting to move the rollator again.



COMFORTABLE SLING SEAT WITH BACKREST



A user comfortably seated on the rollator's sling seat, demonstrating its use as a resting spot.

### Folding and Transporting:

The rollator's folding design makes it easy to store and transport.

- Ensure the storage bag is empty and any items are removed from the seat.
- Locate the folding mechanism, typically a strap or handle on the seat.
- Pull the strap/handle upwards to collapse the rollator frame.
- Once folded, the rollator can be easily carried or stored in compact spaces.



**EASY TO FOLD AND TAKE ANYWHERE**

A person demonstrating how easily the rollator folds for transport, making it convenient for travel or storage.

### Removable Storage Bag:

The rollator includes a removable storage bag for personal items.

- The bag attaches securely under the seat.
- It can be easily detached and carried using its shoulder strap.



## REMOVABLE BAG WITH SHOULDER STRAP

A user carrying the removable storage bag, highlighting its convenience and portability.

### 4. MAINTENANCE

Regular maintenance ensures the longevity and safe operation of your rollator.

- **Cleaning:** Wipe down the frame and seat regularly with a damp cloth and mild detergent. Avoid abrasive cleaners.
- **Wheel Inspection:** Periodically check the wheels for wear and tear. Ensure they spin freely and are free of debris.
- **Brake Check:** Test the brakes regularly to ensure they engage and disengage properly. If the brakes feel loose or too tight, consult a qualified technician.
- **Fastener Check:** Ensure all screws, bolts, and fasteners are tight. Do not overtighten.
- **Storage:** Store the rollator in a dry place away from extreme temperatures.

## 5. TROUBLESHOOTING

This section addresses common issues you might encounter with your rollator.

Problem	Possible Cause	Solution
Rollator does not fold/unfold smoothly.	Mechanism obstructed or stiff.	Check for debris. Apply a small amount of silicone lubricant to moving parts if necessary.
Brakes are not engaging or are too loose/tight.	Brake cable adjustment needed.	Consult the manufacturer's instructions or a qualified technician for brake adjustment. Do not attempt if unsure.
Wheels are squeaking or not rolling freely.	Dirt/debris in wheel bearings or loose axle.	Clean wheels and axles. Check for loose nuts or bolts.
Handles slip after adjustment.	Adjustment pins not fully engaged.	Ensure adjustment pins snap securely into place. Do not use if handles are not locked.

If you encounter issues not listed here or if solutions do not resolve the problem, please contact Stander customer support.

## 6. PRODUCT SPECIFICATIONS

Attribute	Value
Product Dimensions	24 x 25 x 30.5 inches (unfolded)
Item Weight	15.5 Pounds
Folded Footprint	8 inches wide
Handle Height Adjustment	30.5 to 36 inches
Seat Dimensions	17 inches (width) x 7 inches (depth)
Weight Capacity	Up to 300 pounds
Front Wheel Size	10 inches
Rear Wheel Size	8 inches
Material	Steel
Color	Black





## ADJUSTABLE HANDLES AND BRAKES

Visual representation of the rollator's key dimensions and features.

### 7. WARRANTY AND SUPPORT

For detailed warranty information and customer support, please refer to the official Stander website or the documentation included with your product.

A digital version of the user manual is available for download:

[Download User Manual  
\(PDF\)](#)

For further assistance, please contact Stander customer service through their official channels.






### 8. PRODUCT VIDEOS



No official product videos from the seller were found to embed in this manual.

© 2023 Stander. All rights reserved.  
For more information, visit the [Stander Store on Amazon](#).

## Related Documents - B086VTHHVH

	<p><a href="#">Wonder Rollator User Manual: Features, Safety, and Operation</a></p> <p>Comprehensive guide for the Stander Wonder Rollator, covering features, safety precautions, warranty, and operating instructions for this compact folding mobility aid.</p>
	<p><a href="#">Stander Wonder Tray 8800 Owners Manual and Assembly Instructions</a></p> <p>Official owners manual and assembly guide for the Stander Wonder Tray (Model 8800). Provides detailed safety warnings, product specifications, limited lifetime warranty information, and step-by-step instructions for installation.</p>
	<p><a href="#">Stander Bed Cane Assembly and Safety Instructions</a></p> <p>Comprehensive guide to the safe assembly and installation of the Stander Bed Cane, a mobility aid designed to assist individuals in standing from bed. Includes important safety warnings, warranty information, and step-by-step instructions.</p>
	<p><a href="#">Stander 7630 Wonder Tray Premium Owners Manual and Assembly Instructions</a></p> <p>Comprehensive owners manual and assembly instructions for the Stander 7630 Wonder Tray Premium, a mobility aid designed to assist users when standing from chairs or couches. Includes safety warnings, specifications, and step-by-step assembly guidance.</p>
	<p><a href="#">Stander Security Pole 1150: Safety Warnings, Specifications, and Assembly Instructions</a></p> <p>Comprehensive guide to the Stander Security Pole (Model 1150), including essential safety warnings, product specifications, limited warranty details, and step-by-step assembly instructions for safe and effective installation.</p>

Comprehensive owner's manual and assembly guide for the Stander 8800XL Wonder Tray XL. This document provides detailed assembly steps, safety warnings, product specifications, and limited lifetime warranty information for the Stander Wonder Tray XL, designed for use with chairs, couches, and sofas.