



Manuals.plus /

› ELENKER /

› ELENKER Knee Walker Instruction Manual

## ELENKER B0876QN2TF

# ELENKER Knee Walker Instruction Manual

Model: B0876QN2TF

Brand: ELENKER

## 1. OVERVIEW

---

The ELENKER Best Value Knee Walker is a steerable medical scooter designed as a crutch alternative for individuals recovering from foot, ankle, or leg injuries, or those undergoing foot or ankle surgery. It provides stable and comfortable mobility, suitable for both indoor and outdoor use.



Image: The ELENKER Knee Walker, showcasing its red frame, black knee pad, handlebars, and four wheels, including a front basket.

## 2. KEY FEATURES

- **Crutch Alternative:** Ideal for adults recovering from foot, ankle, or leg injuries, or post-surgery.
- **Durable Wheels:** Equipped with 10-inch PVC front wheels and 8-inch PVC rear wheels for smooth navigation.
- **Dual Braking System:** Features dual adjustable locking handbrakes and dual rear on-wheel brakes for enhanced safety and control.

- **Adjustable Design:** Fully adjustable handlebars (28.15" - 37.2") and knee pad (16.93" - 20.86") to accommodate various user heights.
- **Lightweight & Portable:** Weighs 22 lbs, designed for easy transport and storage with a compact folding front center frame.
- **High Weight Capacity:** Supports users up to 300 lbs.
- **Comfortable Knee Pad:** Memory foam knee pad evenly distributes pressure, reducing friction and skin irritation.
- **Detachable Cloth Bag:** Convenient for carrying personal belongings.

### Detailed Feature Views:



## Textured 10" Front Wheels

Featuring textured wheels, it is more non-slip, durable and provides greater stability for those with poor balance

Image: Close-up view of the knee walker's 10-inch textured front wheels, highlighting their size and durability for various terrains.

# Adjustable Design

Adjustable knee pad and handlebars to fit you well



Image: Diagram illustrating the adjustable design of the knee pad and handlebars, with height measurements.



# Dual Braking System

Supports immediate stop and slowdown

Image: Close-up of the dual braking system on the handlebars, showing the levers for immediate stop and slowdown.

## Locking Brake

This device is designed with ergonomic handgrips that make it simple to control. The handles include a brake to help safely slow down or stop. Easy-to-use locks, push down brakes to lock.



Image: Detailed view of the comfortable and soft memory foam knee pad, designed for pressure distribution.

## Sturdy Wheels

Knee walker feature 8" PVC wheels, it is suitable for indoor, outdoor.



Image: The detachable cloth bag, providing convenient storage for personal items.

### 3. SETUP

---

The ELENKER Knee Walker is designed for straightforward assembly and disassembly, making it easy to transport and store. While specific step-by-step instructions are typically included with the product packaging, here are general guidelines:

1. **Unpacking:** Carefully remove all components from the packaging. Ensure all parts listed in the product manual are present.
2. **Frame Assembly:** Connect the main frame components. The front center frame is designed to fold for compact storage. Ensure all locking mechanisms are securely engaged.
3. **Handlebar Installation:** Insert the handlebar post into the designated slot on the frame. Secure it using the provided locking mechanism.
4. **Knee Pad Attachment:** Attach the knee pad to its support structure. Ensure it is firmly seated.
5. **Wheel Check:** Verify that all four wheels are securely attached and rotate freely.
6. **Brake System Connection:** Ensure brake cables are properly routed and connected to the handbrakes and rear wheels. Test brake functionality before use.
7. **Adjustments:** Adjust the height of the handlebars and the knee pad to a comfortable and safe position for the user. Refer to the "Adjustable Design" feature for height ranges.
8. **Accessory Attachment:** If included, attach the detachable cloth bag to the front of the knee walker.

# Save Space

Compact design for easy transport and storage



Image: The knee walker demonstrating its compact design for easy transport and storage, with a close-up of the folding mechanism.

## 4. OPERATING INSTRUCTIONS

Using the ELENKER Knee Walker effectively and safely is crucial for recovery and mobility. Follow these guidelines for optimal use:

- **Positioning:** Place the injured leg with the knee bent onto the padded knee platform. Ensure your weight is evenly distributed.
- **Propulsion:** Use your uninjured foot to push off the ground, propelling the knee walker forward. Maintain a steady and controlled pace.
- **Steering:** The steerable front wheels allow for easy navigation. Use the handlebars to guide the direction of the walker.
- **Braking:**
  - Squeeze the handbrake levers on the handlebars to slow down or stop.
  - To engage the parking brake, squeeze the handbrake levers and push down the small locking tab on each lever. This

will lock the brakes for stationary use.

- To release the parking brake, pull the handbrake levers back towards the handlebars.

- **Terrain:** The knee walker is suitable for various terrains, including indoor surfaces, pavement, and even grass. Exercise caution on uneven or slippery surfaces.

- **Safety Precautions:**

- Always ensure the knee pad and handlebars are adjusted to your comfortable height before use.
- Avoid excessive speed, especially when turning or on inclines/declines.
- Be mindful of obstacles and transitions between different floor types (e.g., carpet to hard floor) to prevent sudden stops or falls.
- Do not exceed the maximum weight capacity of 300 lbs.

## Suitable for Various Terrains

Suitable for living room, office, outdoors and various roads.



Living room



Office



Outdoors



Image: A user demonstrating the knee walker's suitability for various terrains, including outdoor environments.

### 5. MAINTENANCE

Regular maintenance ensures the longevity and safe operation of your ELENKER Knee Walker:

- **Cleaning:** Wipe down the frame and knee pad regularly with a damp cloth and mild soap. Avoid abrasive cleaners.
- **Wheel Inspection:** Periodically check the wheels for wear and tear, and ensure they are free of debris that could impede smooth rolling.
- **Brake Check:** Test the brakes regularly to ensure they are responsive and effective. If they feel loose or ineffective, they may need adjustment.
- **Fastener Tightness:** Check all bolts, screws, and quick-release mechanisms to ensure they are securely tightened. Loose fasteners can affect stability.
- **Storage:** Store the knee walker in a dry, clean place away from extreme temperatures. Utilize its compact folding feature for convenient storage.

## 6. TROUBLESHOOTING

If you encounter issues with your ELENKER Knee Walker, refer to these common troubleshooting tips:

Problem	Possible Cause	Solution
Wheels not rolling smoothly	Debris in wheels, loose axles, worn bearings.	Clean wheels, check and tighten axle bolts, consider professional service if bearings are worn.
Brakes are ineffective or loose	Loose brake cables, worn brake pads, improper adjustment.	Adjust brake cable tension, inspect brake pads for wear and replace if necessary, ensure proper brake lever engagement.
Knee walker feels unstable	Loose frame connections, improper handlebar/knee pad height.	Check and tighten all frame fasteners. Adjust handlebar and knee pad height for proper balance and comfort.
Difficulty steering	Front wheel mechanism obstruction, loose steering column.	Clear any obstructions from the front wheel area. Ensure the steering column is securely fastened.

If the issue persists after attempting these solutions, please contact ELENKER customer support for further assistance.

## 7. SPECIFICATIONS

Attribute	Detail
Brand	ELENKER
Model (ASIN)	B0876QN2TF
Color	Red
Material	Polyvinyl Chloride (PVC wheels)
Package Dimensions	19.25 x 15 x 10.5 inches
Item Weight	21.6 Pounds
Weight Capacity	300 lbs
Handlebars Height Range	28.15" - 37.2"
Knee Pad Height Range	16.93" - 20.86"
Front Wheel Size	10 inches
Rear Wheel Size	8 inches
Age Range (Description)	Adult
Date First Available	May 4, 2018

## 8. WARRANTY AND SUPPORT

---

Specific warranty details for the ELENKER Knee Walker are typically provided with the product packaging or on the manufacturer's official website. Please refer to these resources for comprehensive warranty information.

For any questions, concerns, or technical support regarding your ELENKER Knee Walker, please contact ELENKER customer service. Contact information can usually be found in the product manual or on the ELENKER brand's official website.

