

Goplus B08BCNL97G

Goplus 3-in-1 Stand-Assist Folding Walker User Manual

1. INTRODUCTION

The Goplus 3-in-1 Stand-Assist Folding Walker is a versatile mobility aid designed to provide support and stability for individuals with limited mobility. It can be utilized as a standard walking aid, a stand-up assist, and a toilet safety rail. Constructed from durable aluminum alloy, this walker offers a stable structure and a high loading capacity of up to 400 lbs, ensuring reliable support. Its ergonomic design and adjustable features cater to various user needs, promoting independence and safety in daily activities.

2. KEY FEATURES

- **Stable Structure & High Loading Capacity:** Made of aluminum alloy, providing a firm structure and supporting up to 400 lbs.
- **Comfortable Handle & Ergonomic Design:** Durable PVC handles with two height options, ergonomically designed for comfort and wear resistance. FDA certified for safety.
- **Adjustable Height & Two Modes:** Features 8-level height-adjustable support feet and can be easily converted between fixed mode and interactive mode.
- **Tool-free Installation & Foldable Design:** Convenient one-click folding for compact storage and effortless transport.
- **Multifunctional Walker:** Suitable for use as a stand-up aid, walk aid, and toilet safety rail in various settings like hospitals, homes, and nursing facilities.

3. SETUP AND ASSEMBLY

The Goplus 3-in-1 Folding Walker is designed for tool-free installation. Follow these steps to set up your walker:

1. **Unfold the Walker:** Gently pull the sides of the walker apart until it locks into its open position. Ensure all locking mechanisms are securely engaged.
2. **Adjust the Height:** Press the snap buttons on each leg and slide the leg up or down to your

preferred height. Ensure both sides are adjusted to the same height for stability. There are 8 levels of height adjustment.

3. **Convert to Desired Mode:** The walker can be used in either fixed mode or interactive mode. Locate the mode adjustment mechanism (typically a button or lever near the top of the frame) and switch to your desired setting.
4. **Attach/Detach Wheels (Optional):** The walker comes with two front wheels. These can be attached or detached based on user preference for mobility or stability. To attach, insert the wheel legs into the designated slots until they click securely. To detach, press the release pin and pull the wheel legs out.

Once assembled, ensure all parts are securely fastened before use.



Image: An elderly man using the Goplus 3-in-1 Folding Walker for support while standing near a doorway.



Ergonomic
Design



Adjustable
Height



Toolless
Installation



High Loading
Capacity



FDA
Certification



Variable
Mode

Image: The Goplus 3-in-1 Folding Walker standing upright, highlighting its ergonomic design, adjustable height, toolless installation, high loading capacity, FDA certification, and variable mode features.

4. OPERATION

The Goplus 3-in-1 Folding Walker offers multiple functionalities to assist with daily mobility:

- **As a Walking Aid:** Use the walker to support your weight and maintain balance while walking. Ensure a steady gait and keep the walker close to your body.
- **As a Stand-Up Aid:** Position the walker in front of you when rising from a seated position. Use the lower handles for initial leverage and then the upper handles for full standing support.
- **As a Toilet Safety Rail:** Place the walker around a toilet to provide stable handholds for sitting down and standing up. Ensure the walker is stable and does not interfere with toilet use.

Always ensure the walker is on a stable, level surface before applying weight. Avoid using on uneven terrain or stairs without proper assistance.



FIXED WALKING MODE

Lift the left side of the walker, move the left foot forward, and then lift the right side and move the right foot forward

PRESS THE SNAP BUTTONS ON BOTH SIDES INTO THE UPPER HOLE TO USE IN FIXED MODE. LIFT THE WALKER TO MOVE FORWARD



Image: A composite image demonstrating the walker's three functions: as a walking assist, a getting-up assist from a chair, and a toilet safety rail.



Image: A close-up of the walker's leg, illustrating the 'One Click to Switch Two Modes' feature between Fixed Mode and Interactive Mode.

5. MAINTENANCE

Proper maintenance ensures the longevity and safe operation of your Goplus Folding Walker:

- **Cleaning:** Wipe down the frame and handles regularly with a damp cloth and mild soap. Avoid abrasive cleaners that may damage the finish.
- **Inspection:** Periodically check all screws, bolts, and locking mechanisms to ensure they are tight and secure. Inspect the rubber tips and wheels for wear and tear; replace if necessary.
- **Storage:** When not in use, fold the walker for compact storage. Store in a dry place away from extreme temperatures and direct sunlight.

6. TROUBLESHOOTING

If you encounter any issues with your Goplus Folding Walker, consider the following:

- **Instability:** Ensure all legs are adjusted to the same height and that the mode is correctly set (fixed or interactive). Check that all connections are secure.

- **Difficulty Folding/Unfolding:** Verify that no parts are obstructed and that the one-click folding mechanism is fully engaged or disengaged. Apply gentle, even pressure.
- **Worn Parts:** If rubber tips or wheels show significant wear, they should be replaced immediately to maintain safety and stability.

For persistent issues or concerns, please contact Goplus customer support.

7. SPECIFICATIONS

Product Dimensions	20.5"D x 18"W x 36"H (52.07 cm D x 45.72 cm W x 91.44 cm H)
Item Weight	5.5 Pounds (2.49 kg)
Material	Aluminum
Color	Silver
Maximum Weight Capacity	400 lbs



Image: A diagram illustrating the key dimensions of the Goplus 3-in-1 Folding Walker, including width, depth, and adjustable height range.

8. WARRANTY AND SUPPORT

For warranty information, product support, or to inquire about replacement parts, please refer to the contact information provided with your purchase or visit the official Goplus website. Keep your purchase receipt as proof of purchase for any warranty claims.