

[Manuals.plus](#) /

› [KAMER](#) /

› KAMER Extra Wire Saw 12 Metres User Manual

KAMER 103-311

KAMER Extra Wire Saw 12 Metres User Manual

Model: 103-311

INTRODUCTION

This manual provides essential instructions for the safe and effective use of the KAMER Extra Wire Saw 12 Metres. Please read this manual thoroughly before using the product and retain it for future reference.

The KAMER Extra Wire Saw is designed for specific applications, including adult bovine and embryotomy procedures. Its construction features an extremely wear-resistant wire, ensuring durability and performance.

PRODUCT OVERVIEW



Figure 1: KAMER Extra Wire Saw 12 Metres, showing the wire coiled within its dispenser box, labeled "KAMER -12m-" and "103-311".

The KAMER Extra Wire Saw comes in a convenient 12-meter dispenser box, designed for easy storage and deployment. The wire itself is highly durable and resistant to wear, making it suitable for demanding tasks.

SETUP

Before use, ensure you have the necessary wire saw handles (reference 103 312, sold separately) compatible with this wire saw.

1. **Unpack:** Carefully remove the wire saw dispenser box from its packaging.
2. **Prepare Wire:** Locate the end of the wire protruding from the dispenser box.
3. **Attach Handles:** Thread the ends of the wire through the designated attachment points on your wire saw handles (ref 103 312). Ensure the wire is securely fastened to prevent slippage during operation.
4. **Measure Length:** Pull out the required length of wire for your specific task. The dispenser box allows for easy measurement and cutting.

OPERATING INSTRUCTIONS

This wire saw is intended for specific professional applications. Always follow established safety protocols and guidelines relevant to your field of work.

1. **Safety First:** Always wear appropriate personal protective equipment (PPE), including gloves and eye protection, when handling and operating the wire saw.
2. **Positioning:** Position the wire saw correctly around the material to be cut. Ensure a clear path for the wire and that no obstructions are present.
3. **Cutting Motion:** Apply steady, controlled back-and-forth motions with the wire saw handles. Allow the wire's abrasive properties to do the work, avoiding excessive force.
4. **Maintain Tension:** Keep consistent tension on the wire throughout the cutting process to ensure efficient cutting and prevent tangling.
5. **Monitor Wear:** Periodically inspect the wire for signs of wear or damage during prolonged use. Replace the wire if significant wear is observed.

Important: This product is specifically designed for adult bovine and embryotomy procedures. Use outside of these intended applications may result in damage to the product or injury.

MAINTENANCE

- **Cleaning:** After each use, clean the wire and dispenser box to remove any debris or residue. Use a damp cloth and mild disinfectant if necessary. Ensure the wire is dry before storage.
- **Storage:** Store the KAMER Extra Wire Saw in its original dispenser box in a cool, dry place, away from direct sunlight and corrosive materials.
- **Inspection:** Before each use, inspect the wire for any kinks, fraying, or signs of excessive wear. A worn wire can break unexpectedly.
- **Disposal:** Dispose of worn or damaged wire responsibly according to local regulations for medical or veterinary waste, if applicable.

TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Wire not cutting effectively.	Wire is dull or worn. Incorrect cutting technique.	Replace the wire. Ensure proper tension and consistent back-and-forth motion.
Wire breaks during use.	Excessive force applied. Wire is old or damaged. Kink in the wire.	Reduce force. Inspect wire for damage before use and replace if necessary. Ensure wire is straight before cutting.
Wire tangles in dispenser.	Improper storage or handling.	Ensure wire is neatly coiled when returning to dispenser. Avoid pulling wire out at sharp angles.

SPECIFICATIONS

- Model:** 103-311
- ASIN:** B08NHLP2Q2
- Length:** 12 Metres
- Material:** Extremely wear-resistant wire
- Intended Use:** Adult bovine and embryotomy procedures
- Compatibility:** For use with wire saw handles, ref 103 312 (sold separately)
- Item Weight:** 3.52 ounces (approx. 100 grams)
- Package Dimensions:** 8.66 x 6.3 x 3.15 inches (approx. 22 x 16 x 8 cm)
- Manufacturer:** KAMER

WARRANTY AND SUPPORT

For information regarding warranty coverage or technical support for your KAMER Extra Wire Saw, please contact your point of purchase or the official KAMER distributor in your region.

Please have your product model (103-311) and ASIN (B08NHLP2Q2) available when contacting support.

