

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [Yuwell](#) /

› [Yuwell M103 Portable Mesh Nebulizer Instruction Manual](#)

Yuwell M103

Yuwell M103 Portable Mesh Nebulizer Instruction Manual

Model: M103 | Brand: Yuwell

1. INTRODUCTION

Thank you for choosing the Yuwell M103 Portable Mesh Nebulizer. This device is designed to efficiently deliver medication to the respiratory system. Its compact, lightweight, and quiet design makes it ideal for use at home or on the go. Please read this manual thoroughly before using the device to ensure safe and effective operation.

2. SAFETY INFORMATION

Important Safety Precautions:

- Always consult your physician or pharmacist to determine the appropriate medication, dosage, and nebulization regimen.
- Do not use the device for any purpose other than its intended use as a nebulizer.
- Keep the device out of reach of children and pets.
- Do not immerse the main unit in water or other liquids.
- Use only the accessories provided or recommended by Yuwell.
- Ensure the device is clean and dry before and after each use to prevent contamination.
- Do not operate the device if it has been dropped, damaged, or exposed to extreme temperatures.
- Disconnect the device from the power source before cleaning.
- If you experience any discomfort or adverse reactions during use, stop immediately and consult a medical professional.

3. PACKAGE CONTENTS

The Yuwell M103 Portable Mesh Nebulizer package typically includes:

- Yuwell M103 Main Unit
- Medication Cup
- Adult Mask
- Child Mask

- Mouthpiece
- USB Charging Cable
- Storage Pouch
- User Manual (this document)

4. PRODUCT OVERVIEW

Familiarize yourself with the components of your Yuwell M103 Nebulizer:



Image 1: Yuwell M103 Portable Mesh Nebulizer and its accessories. This image shows the compact main unit, along with the adult and child masks, and the mouthpiece, illustrating the complete set of components for medication delivery.

1. **Main Unit:** Contains the power button, indicator lights, and the mesh nebulization technology.
2. **Medication Cup:** Where liquid medication is placed for nebulization. Attaches to the main unit.
3. **Mesh:** The vibrating mesh inside the medication cup that converts liquid medication into a fine mist.

4. **Mouthpiece:** Used for direct inhalation of the mist.
5. **Masks (Adult/Child):** Used for inhalation, especially for those who find mouthpieces difficult to use, such as children.
6. **USB Charging Port:** For recharging the internal battery.

5. SETUP

Follow these steps to prepare your nebulizer for use:

1. **Charge the Device:** Before first use, fully charge the nebulizer using the provided USB cable. Connect the cable to the device's USB port and a suitable USB power adapter (not included). The indicator light will show charging status.
2. **Assemble the Medication Cup:** Ensure the medication cup is clean and dry.
3. **Add Medication:** Open the medication cup and carefully pour the prescribed liquid medication into it. Do not overfill.
4. **Attach Medication Cup:** Securely attach the medication cup to the main unit.
5. **Attach Accessory:** Choose either the mouthpiece, adult mask, or child mask and firmly attach it to the medication cup outlet.

6. OPERATING INSTRUCTIONS

Once assembled and filled, operate the nebulizer as follows:

1. **Power On:** Press and hold the power button on the main unit for a few seconds until the indicator light illuminates and nebulization begins.
2. **Inhale:**
 - **Using Mouthpiece:** Place the mouthpiece in your mouth, sealing your lips around it. Breathe in slowly and deeply through your mouth until all medication is nebulized.
 - **Using Mask:** Place the mask over your nose and mouth, ensuring a snug fit. Breathe normally until all medication is nebulized.
3. **Treatment Duration:** Continue inhalation until the medication cup is empty and no more mist is produced. The device may automatically shut off when the medication is depleted.
4. **Power Off:** If the device does not shut off automatically, press and hold the power button to turn it off.

Note: The device is designed for quiet operation, ensuring a comfortable experience for all users.

7. CLEANING AND MAINTENANCE

Proper cleaning and maintenance are crucial for the longevity and hygienic operation of your nebulizer.

After Each Use:

1. **Discard Remaining Medication:** Empty any leftover medication from the cup.
2. **Rinse Components:** Detach the medication cup, mouthpiece, and mask. Rinse these parts thoroughly under warm running water.
3. **Clean Mesh:** For the mesh, add a small amount of distilled water to the medication cup and run the device for 1-2 minutes to clear any residue.
4. **Dry:** Shake off excess water and allow all components to air dry completely on a clean cloth or paper towel. Do not wipe the mesh.

Weekly Cleaning (or as needed):

For more thorough disinfection, you may sterilize the medication cup, mouthpiece, and masks by boiling them in water for 5-10 minutes or using a medical disinfectant solution as recommended by your physician or the manufacturer.

Do NOT clean the main unit with water or immerse it. Wipe the main unit with a soft, damp cloth if needed, then dry immediately.

Storage: Store the clean and dry nebulizer and its accessories in the provided storage pouch in a cool, dry place, away from direct sunlight.

8. TROUBLESHOOTING

If you encounter issues with your Yuwell M103 Nebulizer, refer to the following table:

Problem	Possible Cause	Solution
Device does not turn on.	Low battery or no charge.	Charge the device using the USB cable.
No mist or weak mist production.	Insufficient medication. Medication cup or mesh is clogged. Medication cup not properly attached.	Add more medication (up to max fill line). Clean the medication cup and mesh thoroughly. Re-attach the medication cup securely.
Device stops during operation.	Medication depleted. Low battery.	This is normal; refill if more treatment is needed. Recharge the device.
Indicator light flashing.	Low battery warning.	Recharge the device.

If the problem persists after trying these solutions, please contact Yuwell customer support.

9. SPECIFICATIONS

Feature	Detail
Model Number	M103
Brand	Yuwell
Unit Count	1 Count
Item Weight	Approximately 97 grams
Package Dimensions	16.2 x 15.8 x 8 cm
Power Source	Rechargeable USB Battery
Nebulization Method	Mesh Technology

10. WARRANTY AND SUPPORT

For warranty information, please refer to the warranty card included with your product or visit the official Yuwell website. Keep your purchase receipt as proof of purchase.

If you require technical assistance, have questions about your device, or need to report a problem, please contact Yuwell

customer support through their official channels. Contact details can typically be found on the product packaging or the manufacturer's website.

