

MAIDeSiTe JC35M3

MAIDeSiTe Electric Bath Lift Handheld Controller JC35M3 Instruction Manual

Model: **JC35M3** | Brand: **MAIDeSiTe**

1. INTRODUCTION

The MAIDeSiTe Electric Bath Lift Handheld Controller JC35M3 is designed to provide convenient and safe operation of your MAIDeSiTe bath lift chair. This rechargeable, waterproof (IPX8) controller ensures ease of use and enhanced mobility for individuals requiring assistance with bathing. This manual provides detailed instructions for proper setup, operation, and maintenance to ensure optimal performance and user safety.

2. SAFETY INFORMATION

- **Legal Disclaimer:** Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.
- **Attention:** Our lift chair accessories come with a one-year warranty. If you need to purchase a controller, please contact us, and we can provide it for free within the warranty period.
- Ensure the controller is fully charged before use to prevent sudden power failure during operation. Refer to the indicator lights for battery status.
- Always ensure the bath lift is securely positioned with its suction cups before use.
- Do not immerse the charger in water. Only the handheld controller is waterproof.

3. PACKAGE CONTENTS

Please verify that all items are present in your package:

- MAIDeSiTe Electric Bath Lift Handheld Controller (JC35M3)

- Handheld Controller Charger Link (Note: Charger not included, but a dedicated link is provided for charging the remote.)
- Instruction Manual
- Bath Lift Base Unit (with integrated suction cups)
- Backrest and Seat Components (with blue padded cushions)

4. SETUP INSTRUCTIONS

4.1. Assembling the Bath Lift

1. Carefully remove all components from the packaging.
2. Place the main base unit on a flat, stable surface.
3. Attach the backrest to the base unit by aligning the hooks and pushing down until it clicks securely into place.
4. Extend the side flaps of the seat. These will automatically fold in when placed in a bathtub.

4.2. Connecting the Handheld Controller

1. Locate the connection port on the bath lift's main unit.
2. Plug the handheld controller's cable firmly into this port.
3. Ensure the connection is secure to prevent water ingress and ensure proper function.

4.3. Placing the Bath Lift in the Bathtub

1. Ensure the bathtub surface is clean and free of debris.
2. Carefully lower the bath lift into the bathtub. The suction cups on the base will adhere to the tub surface, providing stability.
3. Press down firmly on the seat to ensure all suction cups are securely attached.

Bath Lift



Image: MAIDeSITe Electric Bath Lift positioned securely in a bathtub. The blue padded seat and backrest are visible, with side flaps extended.

5. OPERATING INSTRUCTIONS

5.1. Handheld Controller Functions



Image: Close-up of the MAIDeSiTe Electric Bath Lift Handheld Controller, showing the red emergency stop button, green up/down arrows, and indicator lights.

- **Green Up Arrow:** Press and hold to raise the bath lift.
- **Green Down Arrow:** Press and hold to lower the bath lift.
- **Red Emergency Stop Button:** Press to immediately halt all movement of the bath lift.
- **Indicator Lights:**
 - **Green Light:** Indicates sufficient power for normal use.
 - **Blue Light:** Indicates insufficient power; charging is recommended.
 - **Red Light:** Indicates critically low power; immediate charging is required to prevent sudden power failure.

5.2. Using the Bath Lift

1. Use the handheld controller to raise the seat to a comfortable height for entry.
2. Carefully transfer onto the bath lift seat.
3. Use the green down arrow button to slowly lower yourself into the bathtub. The backrest will recline as the seat lowers.
4. To exit the tub, use the green up arrow button to raise the seat and backrest to a comfortable height for transfer.

- The side flaps of the seat will automatically adjust as the lift moves up and down, providing a wider surface for transfer when raised and folding in for a compact fit when lowered.

Your browser does not support the video tag.

Video: An official product video demonstrating the MAIDeSiTe Electric Bath Lift in operation, showing how the seat and backrest move up and down using the handheld controller. The video highlights the ease of use for entering and exiting a bathtub.

6. MAINTENANCE

- Charging:** When the blue or red indicator light appears on the handheld controller, it needs to be charged. Disconnect the controller from the bath lift and connect it to the dedicated charger link. Plug the charger into a wall outlet.
- Cleaning:** Clean the bath lift and controller with a mild, non-abrasive cleaner and a soft cloth. Do not use harsh chemicals.
- Storage:** The bath lift is foldable for convenient storage. After use, ensure it is dry before folding and storing in a cool, dry place.

7. TROUBLESHOOTING

Problem	Possible Cause	Solution
Bath lift does not operate.	Low battery or controller not connected.	Check indicator lights on the controller. If blue or red, charge the controller. Ensure the controller is securely plugged into the bath lift.
Bath lift is unstable in the tub.	Suction cups not properly adhered.	Ensure the bathtub surface is clean and dry. Press down firmly on the seat to activate the suction cups.
Controller lights are off.	Controller is not in use or completely discharged.	The lights only activate when a button is pressed. If no lights appear when pressing a button, the battery may be completely drained and requires charging.

8. SPECIFICATIONS

- Model:** JC35M3
- Brand:** MAIDeSiTe
- Waterproof Grade:** IPX8 (Handheld Controller)
- Battery:** 1 Lithium Ion battery (included)

- **Weight Capacity:** Up to 300 lbs (approx. 136 kg)
- **Minimum Height (when lowered):** Approximately 3 inches
- **Maximum Height (when raised):** Approximately 20.5 inches
- **Seat Cushion Width:** Approximately 12 inches
- **Backrest Cushion Length:** Approximately 21 inches
- **Package Dimensions:** 10.63 x 9.96 x 2.64 inches
- **Weight:** 1.41 Pounds

9. WARRANTY & SUPPORT

MAIDeSiTe lift chair accessories come with a **one-year warranty**. If you require a replacement controller within the warranty period, please contact MAIDeSiTe customer support for assistance. We are committed to providing high-quality products and support to our customers.

For further assistance or inquiries, please refer to the contact information provided with your purchase or visit the official MAIDeSiTe website.