

## Braun Braun Manual Nasal Aspirator

# Braun Manual Nasal Aspirator: User Manual

## INTRODUCTION

The Braun Manual Nasal Aspirator is designed to provide gentle and effective relief from nasal congestion for infants, toddlers, and children. This manual provides detailed instructions for the safe and proper use, cleaning, and maintenance of your device.



Image: The Braun Manual Nasal Aspirator, showing the main unit, flexible tube, and mouthpiece.

## PACKAGE CONTENTS

---

Before first use, please ensure all components are present in the package:

- 1 x Braun Manual Nasal Aspirator unit
- 2 x Soft nose tip sizes (one small for newborns, one larger for toddlers/kids)
- 5 x Disposable filters
- 1 x Carrying pouch

# Package contents



Image: Visual representation of the package contents, including the aspirator, two nose tips, five disposable filters, and a carrying pouch.

## SETUP

Follow these steps to prepare your nasal aspirator for use:

1. Remove all components from the packaging.
2. Select the appropriate nose tip size for your child. The smaller tip is for newborns, and the larger tip is for toddlers and older children.
3. Attach the selected nose tip firmly to the top of the mucus chamber.
4. Ensure a disposable filter is correctly placed inside the mucus chamber to prevent mucus from entering the tube.
5. Connect the flexible tube to the bottom of the mucus chamber and the mouthpiece.

For a visual guide on getting started, please watch the video below:



Video: Braun Manual Nasal Aspirator - Getting Started. This video demonstrates the initial assembly and preparation of the aspirator for use.

## OPERATING INSTRUCTIONS

To effectively clear your child's nasal passages:

1. Position your child comfortably.
2. Gently insert the nose tip into one of your child's nostrils, ensuring a snug seal.
3. Place the mouthpiece in your mouth.
4. Apply gentle, continuous suction through the mouthpiece. The mucus will be drawn from your child's nose into the device and captured in the mucus chamber, with the filter preventing it from reaching the tube.
5. Repeat the process for the other nostril as needed.
6. After each use, disassemble and clean the aspirator immediately.



Image: A parent demonstrating the use of the Braun Manual Nasal Aspirator on an infant, showing proper positioning and suction technique.

## MAINTENANCE AND CLEANING

Proper cleaning after each use is crucial for hygiene and device longevity:

1. **Disassembly:** Gently twist the mucus chamber to separate it from the main unit. Remove the nose

tip and mouthpiece from the tube. Unscrew the filter holder from the mucus chamber.

2. **Filter Disposal:** Discard the used disposable filter by gently shaking it out of the filter holder. Replace with a new filter for the next use.
3. **Washing:** The mucus chamber, mouthpiece, and nose tip are top-rack dishwasher safe. Alternatively, hand wash these pieces in hot, soapy water.
4. **Tube Cleaning:** The flexible tube should be thoroughly cleaned by boiling it in hot water to ensure proper sanitation.
5. **Drying and Reassembly:** Ensure all parts are completely dry before reassembling the aspirator. Store in the provided carrying pouch.



## Soft suction tips are easy to remove and clean

Image: A person cleaning the soft suction tips under running water, highlighting the ease of cleaning.

For a visual guide on cleaning, please watch the video below:



Video: Braun Manual Nasal Aspirator - How to Clean. This video demonstrates the disassembly and cleaning process for all components of the aspirator.

### SPECIFICATIONS

- **Model Name:** Braun Manual Nasal Aspirator

- **Material:** Plastic
- **Color:** Clear
- **Product Care Instructions:** Dishwasher Safe (top rack for chamber, mouthpiece, nose tip); Tube requires boiling.

## SAFETY INFORMATION

---

Please read and follow all safety instructions before using this product. Keep out of reach of children when not in use. Do not force the nose tip into the nostril. If your child experiences discomfort or bleeding, discontinue use and consult a healthcare professional. This device is intended for nasal aspiration only and should not be used for other purposes.