

## Easyinsmile IPR Strip System

# Easyinsmile Dental IPR Strip System Instruction Manual

Model: IPR Strip System (Black)

## 1. INTRODUCTION

This manual provides detailed instructions for the safe and effective use of the Easyinsmile Dental IPR Strip System. This system is designed for interproximal reduction (IPR), also known as enamel stripping, in orthodontic procedures. It is also suitable for descaling, polishing, grinding, and cleaning tooth surfaces.

The strips are crafted from wear-resistant corundum material, designed to be bendable to conform to natural tooth contours, thereby enhancing patient comfort and procedural efficiency. They are autoclavable and reusable.

## 2. SAFETY INFORMATION

**Important: For professional dental use only. Not for at-home use by individuals without proper training.**

- For external use only.
- Avoid contact with eyes. If irritation occurs, rinse immediately with water.
- Keep out of reach of children.
- Adult supervision is required for use in a professional setting.
- Do not use if the patient has sensitive teeth or gums, or if the patient is pregnant or nursing.
- If any discomfort or adverse effects are experienced by the patient, discontinue use immediately and consult a dental professional.
- Always inspect strips for damage before use. Discard any damaged strips.

## 3. PRODUCT OVERVIEW AND COMPONENTS

The Easyinsmile Dental IPR Strip System consists of 10 individual strips, each designed for specific interproximal reduction tasks. The strips are black in color, with a strip width of 3mm and a thickness of 0.1mm. They feature a serrated edge for initial contact opening.



Figure 1: Ten black EasyinSmile IPR strips. Each strip features a central abrasive section and handles for manipulation.



# 10 Pack



Figure 2: The 10-pack of IPR strips. An inset shows a strip being used on a dental model to demonstrate interproximal access.

Key features include:

- **Bendable Design:** Allows the strip to conform to the natural curvature of teeth, minimizing sharp corners and subgingival ledges.
- **Corundum Material:** Ensures durability and effective abrasion.
- **Autoclavable:** Can be sterilized at temperatures up to 270°F (132°C) for reuse.
- **Serrated Edge:** Facilitates safe and easy initial contact opening between teeth.

## 4. SETUP

1. **Sterilization:** Before first use and after each subsequent use, sterilize the IPR strips using an autoclave at temperatures up to 270°F (132°C). Ensure proper sterilization protocols are followed according to clinical guidelines.
2. **Inspection:** Visually inspect each strip for any signs of damage, wear, or debris before use. Discard any compromised strips to ensure patient safety and optimal performance.
3. **Patient Preparation:** Prepare the patient for the IPR procedure, ensuring adequate isolation and access to the interproximal areas.

## 5. OPERATING INSTRUCTIONS

The Easyinsmile Dental IPR Strip System is primarily used for interproximal reduction, descaling, polishing, grinding, and cleaning the surface of teeth. Follow these steps for effective application:

1. **Selection:** Choose the appropriate strip based on the required reduction or task. The black strips provided in this kit are 3mm wide and 0.1mm thick, featuring a serrated edge for initial contact opening.
2. **Initial Contact Opening (Serrated Edge):**




Figure 3: A close-up view of the serrated edge on an IPR strip, designed for safe initial contact opening.

Carefully insert the serrated edge of the strip into the interproximal contact point. Use gentle, controlled movements to initiate separation without causing trauma to the soft tissues.

3. **Interproximal Reduction/Shaping:**

Color: Black  
Strip width: 3MM  
Diamoand coating: No  
Diamond pattern: No  
Serrated: Yes  
Thickness: 0.1MM  
With holes: No




Figure 4: Demonstrates the use of an IPR strip on a dental model for interproximal reduction. The strip is positioned between two teeth.

Once initial contact is opened, use a back-and-forth motion to reduce the enamel as required. The bendable design allows the strip to adapt to the tooth's curved surface, ensuring a smooth reduction and preventing the formation of sharp edges or steps. Gradually expand the interproximal space according to the treatment plan.

4. **Descaling, Polishing, Grinding, and Cleaning:** For these applications, use the abrasive surface of the strip with appropriate pressure and motion to achieve the desired outcome. The fine abrasive surface is effective for removing minor calculus, polishing rough surfaces, or preparing tooth surfaces.
5. **Monitoring:** Regularly check the amount of reduction using appropriate gauges or visual assessment to avoid over-reduction.
6. **Rinsing:** Periodically rinse the area with water to remove debris and evaluate progress.

## 6. MAINTENANCE

Proper maintenance ensures the longevity and hygienic use of your Easyinsmile IPR strips.

- **Cleaning:** After each use, thoroughly clean the strips to remove all biological debris and contaminants. Use a brush and an appropriate cleaning solution.
- **Sterilization:** The strips are autoclavable. Sterilize them in an autoclave at temperatures up to 270°F

(132°C) for the recommended cycle duration. Ensure they are completely dry before storage.

- **Storage:** Store cleaned and sterilized strips in a dry, clean, and protected environment to prevent contamination and damage.
- **Inspection:** Before each use, inspect the strips for any signs of wear, damage, or dullness. Replace strips that show significant wear to maintain effectiveness and safety.




Figure 5: Key features of the IPR strips, including autoclavable and reusable properties.

## 7. TROUBLESHOOTING

While the Easyinsmile IPR Strip System is designed for reliability, here are some common issues and their resolutions:

- **Difficulty in Initial Contact Opening:** Ensure the serrated edge is correctly positioned. Apply gentle, consistent pressure. If resistance is high, consider using a thinner strip or a different technique as per clinical judgment.
- **Ineffective Reduction:** The strip may be worn out or clogged with debris. Clean the strip thoroughly or replace it with a new one. Ensure correct angulation and pressure are applied during use.
- **Strip Bending/Breaking:** This indicates excessive force or improper technique. Always use gentle, controlled movements. If a strip bends or breaks, immediately cease use and replace it.
- **Patient Discomfort:** Re-evaluate the technique and ensure the strip is conforming to the tooth contour. Use adequate local anesthesia if necessary. If discomfort persists, discontinue use and assess the patient's condition.

## 8. SPECIFICATIONS

Feature	Detail
Brand	Easyinsmile
Model	IPR Strip System (Black)
Quantity	10 pieces per kit
Strip Width	3mm
Thickness	0.1mm
Material	Corundum
Serrated	Yes
Autoclavable	Yes (up to 270°F/132°C)
Product Benefits	Brightening, Enamel Protection, Tartar Control (as per general product description, specific to IPR is space creation)
Package Dimensions	6.3 x 5.39 x 0.83 inches
Item Weight	1.76 ounces




Figure 6: Overview of different IPR strip types and their specifications. This kit specifically includes the black 3mm, 0.1mm serrated strips.

## 9. WARRANTY AND SUPPORT

Easyinsmile products are manufactured to high-quality standards. For specific warranty information or technical support, please refer to the contact details provided with your purchase or visit the official Easyinsmile website. Keep your purchase receipt for warranty claims.

# Exquisite packaging

## IPR Strip System+Instruction manual



Figure 7: The EasyinSmile IPR Strip System packaging, which includes the instruction manual.