

[Manuals.plus](#) /

> [Eschenbach](#) /

> Eschenbach Mobilux LED Handheld Magnifier with Mobase Stand (4X) (Model 1511-3) User Manual

Eschenbach 1511-3

Eschenbach Mobilux LED Handheld Magnifier with Mobase Stand (4X)

Model: 1511-3

INTRODUCTION

This manual provides instructions for the proper setup, operation, and maintenance of your Eschenbach Mobilux LED Handheld Magnifier with Mobase Stand. This device is designed to assist individuals with low vision by providing illuminated magnification for various tasks, offering both handheld and hands-free functionality.



The Eschenbach Mobilux LED Handheld Magnifier, featuring a black handle and a clear rectangular lens with LED illumination.

SETUP

1. Battery Installation

1. Locate the battery compartment on the handle of the magnifier.
2. Slide open the battery compartment cover.
3. Insert two (2) AA batteries, ensuring correct polarity (+/-) as indicated inside the compartment. The device comes with 2 AA batteries included.
4. Close the battery compartment cover securely.

2. Attaching the Mobase Stand

1. To convert your handheld magnifier into a hands-free device, align the base of the Mobilux magnifier with the slot in the Mobase Stand.
2. Gently slide the magnifier into the stand until it is securely seated.

3. Using Snap-on Filters

The Mobilux LED Handheld Magnifier includes three (3) snap-on filters. These filters can be attached to the lens area to adjust the light temperature or reduce glare, enhancing viewing comfort based on individual preference and lighting conditions.

OPERATING INSTRUCTIONS

1. Handheld Operation

1. Hold the magnifier by its ergonomic handle.
2. Locate the ON/OFF light switch on the handle. Slide the switch to the 'ON' position to activate the LED illumination.
3. Position the lens over the object or text you wish to magnify. Adjust the distance between the magnifier and the object until the image is clear and in focus.
4. Slide the switch to the 'OFF' position when finished.



A woman using the Eschenbach Mobilux LED Handheld Magnifier to read a book, demonstrating handheld operation.

2. Hands-Free Operation with Mobase Stand

Once the magnifier is inserted into the Mobase Stand, you have two positions for hands-free use:

- **Upright Position:** This position is ideal for tasks that require both hands, such as threading a needle, filing nails, or engaging in hobbies like crafting. The magnifier is angled to provide a clear view of the workspace beneath the lens.



A woman using the Eschenbach Mobilux LED Magnifier in its Mobase Stand, set in the upright position for detailed work like beading.

- **Flat Position:** For reading larger documents or books, the stand can be positioned flat. This allows you to easily slide the magnifier over the reading material, providing continuous magnification across the page.



A man using the Eschenbach Mobilux LED Magnifier in its Mobase Stand, set in the flat position for reading a newspaper.

MAINTENANCE

Cleaning the Lens

To ensure optimal clarity, keep the lens clean. Use the included Liberty Microfiber Cloth or another soft, lint-free cloth to gently wipe the lens surface. Avoid using abrasive cleaners, paper towels, or harsh chemicals, as these can scratch the lens.



The Liberty Microfiber Cloth, included for cleaning the magnifier lens.

Battery Replacement

When the LED illumination dims or stops working, it is time to replace the batteries. Follow the battery installation steps outlined in the Setup section. Always use fresh AA batteries.

Storage

When not in use, store the magnifier in its slip-on protective case to prevent scratches and damage to the lens and body.

TROUBLESHOOTING

- **Light does not turn on:**
 - Ensure the ON/OFF switch is fully in the 'ON' position.
 - Check if the batteries are inserted correctly with the correct polarity.
 - Replace the batteries with fresh AA batteries.
- **Image is blurry or unclear:**
 - Adjust the distance between the magnifier lens and the object being viewed.
 - Clean the lens surface using a microfiber cloth to remove any smudges or dust.

SPECIFICATIONS

Feature	Specification
Model Number	1511-3
Magnification	4x
Diopter	16D
Illumination	SMD LED
Lens Type	Aspheric
Lens Material	PXM Plastic
Lens Size	2.95 in. x 1.96 in. (75 mm x 50 mm)
Power Source	2 AA Batteries (included)
Item Weight	12.8 ounces
Handle Material	Plastic

WARRANTY INFORMATION

The Eschenbach Mobilux LED Handheld Magnifier with Mobase Stand is covered by a **Limited Lifetime Warranty** from Eschenbach. For detailed information regarding warranty terms and conditions, please refer to the official Eschenbach website or contact their customer service.

CUSTOMER SUPPORT

If you encounter any issues or have questions not addressed in this manual, please contact Eschenbach customer service for assistance. You may also reach out to the retailer from whom you purchased the product, LIBERTY Health Supply, for initial support.

For further information, visit the official Eschenbach website:www.eschenbach.com