

[Manuals.plus](#) /

› [ReSound](#) /

› ReSound Linx Quattro Charging Case User Manual

ReSound Linx Quattro

ReSound Linx Quattro Charging Case

USER MANUAL

Introduction

The ReSound Linx Quattro Charging Case is an expertly designed, portable charging solution for your ReSound Linx Quattro hearing aids. This intuitive and easy-to-use case ensures your hearing aids are always ready for use, providing convenience and protection. When hearing aids are placed inside, LED indicators clearly display their battery status.

What's in the Box

Upon opening your ReSound Linx Quattro Charging Case package, you should find the following items:

- 1 x GN ReSound Quattro Charger Case
- 1 x Power plug with cable
- User Guide and other documentation



This image shows the complete contents of the ReSound Linx Quattro Charging Case package. It includes the gray charging case, a power adapter, a USB cable, a small accessory bag, and the user guide booklet, all laid out on a white surface.

Setup

1. **Unpack the Charger:** Carefully remove the charging case, power plug, and cable from the packaging.
2. **Connect Power:** Insert the power cable into the charging case and plug the power adapter into a standard electrical outlet.
3. **Initial Charge:** It is recommended to fully charge the case before first use. The charger has an internal battery that can hold three full charge cycles for your hearing aids.



This image displays the ReSound Linx Quattro Charging Case in its gray color, positioned alongside its white retail box. The case is shown with its lid open, revealing the charging compartments.

Operating Instructions

1. **Placing Hearing Aids:** Open the charging case lid. Gently place your ReSound Linx Quattro hearing aids into their respective charging slots. Ensure they are seated correctly.
2. **Charging Process:** Close the lid. The LED lights on the hearing aids will pulsate, indicating that charging is in progress. Once fully charged, the LEDs will remain steadily on. A full charge typically takes 3 hours.
3. **Battery Status Indicators:**
 - **Hearing Aid Charge:** The five LED lights at the front of the charger indicate the charge level of the hearing aid with the lowest battery.
 - **Charger Battery Level:** The three LEDs at the back of the charger show the remaining battery level of the charging case itself.
4. **Removing Hearing Aids:** Once charging is complete, open the lid and remove your hearing aids. They will automatically turn on when removed from the charging dock.

5. **Portability:** The sleek design allows the case to be easily carried in a pocket or handbag for on-the-go charging. The hard case with a magnetic lid provides physical protection for the hearing aids while they charge.



The ReSound LiNX Quattro Charging Case is depicted in its closed state, showcasing its compact and sleek design. It is placed next to its white retail box, which features the ReSound LiNX Quattro branding.

Maintenance

- **Cleaning:** Use a soft, dry cloth to wipe the exterior of the charging case. Do not use liquids or abrasive cleaners.
- **Storage:** Store the charging case in a cool, dry place away from direct sunlight and extreme temperatures when not in use.
- **Protection:** Keep the lid closed when hearing aids are inside to protect them from dust and physical damage.

Troubleshooting

- **Hearing Aids Not Charging:**

- Ensure the charging case itself is charged or connected to a power source.
 - Verify that the hearing aids are correctly seated in their charging slots.
 - Check for any debris in the charging contacts of both the case and the hearing aids.
- **Charger LEDs Not Illuminating:**
 - Confirm the power cable is securely connected to both the charger and the power outlet.
 - Try a different power outlet or USB port.
 - **Short Battery Life:**
 - Ensure the charging case is fully charged before use.
 - If the issue persists, contact ReSound customer support.

Specifications

Feature	Detail
Brand	ReSound
Model Number	804608480
Color	Gray
Input Voltage	240 Volts
Output Voltage	240 Volts
Battery Charge Time	3 Hours (for hearing aids)
Connector Type	US power plug
Safety Standard	UL Certified

Warranty Information

This ReSound product is covered by a manufacturer's warranty. Please refer to the warranty card included in your product packaging or visit the official ReSound website for detailed information regarding warranty terms, conditions, and registration.

Support

For further assistance, technical support, or inquiries regarding your ReSound Linx Quattro Charging Case, please contact ReSound customer service. Contact details can typically be found on the official ReSound website or in the documentation provided with your product.