

iLiving ILG-622

iLIVING USA Foldable Upright Mobility Rollator (Model ILG-622) - Instruction Manual

Model: ILG-622

INTRODUCTION

This manual provides essential instructions for the safe and effective use, assembly, and maintenance of your iLIVING USA Foldable Upright Mobility Rollator, Model ILG-622. This lightweight walker is designed to enhance mobility and independence for seniors, featuring a foldable design for easy transport and storage, a comfortable padded seat and backrest, and dual brakes with a parking lock for safety.

The iLIVING USA Foldable Upright Mobility Rollator is constructed from durable aluminum, weighing approximately 15 pounds, and can support users up to 300 pounds. It includes 8-inch diameter wheels for smooth navigation across various surfaces.

SETUP AND ASSEMBLY

The iLIVING USA Foldable Upright Mobility Rollator is designed for quick and tool-free assembly. All components snap into place securely.

Assembly Steps:

- Unpack Components:** Carefully remove all parts from the packaging. You should have the main frame, two front legs with wheels, two rear legs, and a storage bag.
- Assemble the Main Frame:** Connect the two main frame sections. Ensure the pins click into place, securing the frame.
- Attach Front Legs:** Insert the front legs, which include the 8-inch swivel wheels, into the designated slots on the front of the main frame. Ensure the wheels are positioned on the outside for better balance. The pins should lock into place.
- Attach Rear Legs:** Insert the rear legs into the designated slots on the back of the main frame. Ensure the pins lock into place.
- Attach Storage Bag:** The included detachable fabric bag can be attached to the front of the rollator using the Velcro straps. Adjust the straps for a secure fit.

For a visual guide on assembly, please refer to the video below:

Video Description: This video demonstrates the step-by-step assembly process of the iLIVING USA Foldable Upright Mobility Rollator, including attaching the wheels and the storage bag. It highlights the tool-free setup and the final appearance of the assembled walker.

Product Images:



Image: All components of the iLIVING USA Foldable Upright Mobility Rollator laid out before assembly.



Image: The iLIVING USA Foldable Upright Mobility Rollator fully assembled, viewed from the front.

OPERATING INSTRUCTIONS

The iLIVING USA Foldable Upright Mobility Rollator is designed for ease of use and portability.

Using the Rollator:

- **Walking Support:** Use the handles for stability while walking. The 8-inch front swivel wheels and fixed rear wheels provide easy steering and glide smoothly over various outdoor surfaces.
- **Braking:** The rollator features dual brakes. Squeeze the brake levers on the handles to slow down or stop.
- **Parking Lock:** To engage the parking lock, push the brake levers down until they click into place. This locks the wheels, allowing the rollator to be used as a stationary chair. Pull the levers up to release the parking lock.
- **Seating:** When tired, engage the parking brake and use the padded seat and backrest for resting. The seat supports up to 300 pounds.
- **Storage Bag:** The detachable fabric bag can carry personal items, groceries, or other essentials.

Folding and Unfolding:

The rollator can be easily folded for storage and transport. To fold, lift the seat strap. The frame will collapse inward. To unfold, push the sides of the frame apart until the seat locks into place, ensuring it is fully open and stable before use.

Watch this video for a demonstration of the folding and unfolding mechanism and answers to common questions:

Video Description: This video addresses frequently asked questions about the iLIVING USA Foldable Upright Mobility Rollator, including how to easily fold and unfold the walker for transport and storage.

Product Images:



Image: The iLIVING USA Foldable Upright Mobility Rollator in a partially folded position, showing its compact design.



Image: The iLIVING USA Foldable Upright Mobility Rollator fully folded, ready for compact storage or transport.



Image: The iLIVING USA Foldable Upright Mobility Rollator with its front storage bag attached, providing convenient carrying space.

MAINTENANCE

Regular maintenance ensures the longevity and safe operation of your rollator.

- **Cleaning:** Wipe down the frame with a damp cloth and mild detergent as needed. The seat and backrest are designed for easy cleaning. Avoid abrasive cleaners.
- **Wheel Inspection:** Periodically check the wheels for wear and tear. Ensure they rotate freely and are free of debris.
- **Brake Check:** Regularly test the dual brakes and parking lock to ensure they are functioning correctly and provide adequate stopping power. If brakes feel loose or ineffective, consult a professional for adjustment.
- **Frame Integrity:** Inspect the aluminum frame for any signs of damage, cracks, or loose connections. Do not use the rollator if any structural damage is observed.

TROUBLESHOOTING

If you encounter any issues with your iLIVING USA Foldable Upright Mobility Rollator, please refer to the following common solutions:

- **Difficulty Folding/Unfolding:** Ensure no obstructions are preventing the frame from collapsing or expanding fully. Check that the seat mechanism is not jammed.
- **Wobbly Frame:** Verify that all pins are securely locked into place after assembly. Ensure all leg extensions are set to the same height if adjustable.
- **Brakes Not Engaging:** Check the brake cables for any kinks or damage. Ensure the brake levers are fully depressed or pulled to engage the parking lock. If the issue persists, professional adjustment may be required.
- **Wheels Not Rolling Smoothly:** Inspect wheels for trapped debris (hair, lint) and remove if present. Ensure wheel axles are not bent or damaged.

If problems persist after attempting these troubleshooting steps, please contact iLIVING USA customer support.

SPECIFICATIONS

Feature	Detail
Product Dimensions (Unfolded)	25"D x 24"W x 36"H
Item Weight	15 Pounds
Material	Aluminum
Color	Blue
Model Number	ILG-622
Weight Capacity	Up to 300 lbs
Wheel Diameter	8 inches (front and back)
Adjustable Height	31 - 36 inches (from floor)

WARRANTY AND SUPPORT

iLIVING USA is committed to customer satisfaction. While specific warranty details are not provided in this document, iLIVING USA offers local customer support to address all questions and concerns regarding your mobility rollator.

For assistance, please contact iLIVING USA customer support directly. Contact information can typically be found on the product packaging or the official iLIVING USA website.