



Manuals.plus /

› YOUSHARES /

› YOUSHARES Insulin Cooler Travel Case User Manual

YOUSHARES MPN_ISL_Case_Y6_PL_01

YOUSHARES Insulin Cooler Travel Case User Manual

Model: MPN_ISL_Case_Y6_PL_01

1. PRODUCT OVERVIEW

The YOUSHARES Insulin Cooler Travel Case is designed to provide a convenient and insulated solution for transporting insulin pens, vials, and other diabetic medication supplies. Its compact and portable design makes it ideal for short trips, ensuring your medication remains at a suitable temperature.

- **Convenient Organizer:** Ideal travel case for carrying medication, diabetes pens, vials, and insulin pen needles.
- **Compact Size:** Light and compact, making it convenient to carry.
- **Good Insulation:** Made of high-grade Oxford cloth material, lined with aluminum foil for effective temperature retention.
- **Durable Design:** Features a durable outer Oxford material and dual zipper pulls for secure closure.



Figure 1.1: YOUSHARES Insulin Cooler Travel Case (Purple) with included ice packs and internal view.



Figure 1.2: Key features of the travel case, including waterproof exterior and dual zippers.

2. SETUP INSTRUCTIONS

Before using your YOUSHARES Insulin Cooler Travel Case, it is essential to properly prepare the included ice packs to ensure optimal cooling performance.

2.1 Preparing the Ice Packs

1. Place the two reusable ice packs in a freezer for a minimum of 8 to 10 hours before you intend to use the insulin bag outside. Ensure they are completely frozen.
2. Once frozen, carefully insert the ice packs into the designated mesh pockets or elastic straps within the insulated compartments of the travel case.

FEATURES

The weight is only 104g.
Perfect for insulin dependent daily use.



Figure 2.1: Internal view of the case, ready for ice packs.

Video 2.1: Demonstration of preparing and inserting ice packs into the travel case.

3. OPERATING INSTRUCTIONS

After preparing the ice packs, you can now organize your medication and supplies within the travel case.

3.1 Organizing Your Supplies

1. Place your insulin pens, vials, or insulin needles into the designated compartments and elastic loops within the bag. The case is designed with various pockets and straps to keep items secure and organized.
2. Ensure all items are neatly arranged and the case is fully zipped to maintain the internal temperature.



Figure 3.1: Organized medical supplies inside the case.

Insulin Cooler Travel Case

- ▶ Great for traveling
- ▶ Perfect size for carry-on
- ▶ Keeping insulin cold
- ▶ TSA approved ice packs



 **Accessories
not included**

Figure 3.2: The travel case in a travel setting.

Video 3.1: Detailed demonstration of how to organize and pack the insulin cooler travel case.

4. MAINTENANCE

Proper maintenance will extend the life of your insulin cooler travel case.

4.1 Cleaning the Case

- Wipe the exterior and interior of the case with a damp cloth as needed.
- For stubborn spots, a mild soap solution can be used, followed by wiping with a clean, damp cloth.
- Allow the case to air dry completely before zipping it closed or storing it.
- Do not machine wash or tumble dry the case. Avoid using harsh chemicals or abrasive cleaners.

5. TROUBLESHOOTING

If you encounter issues with your insulin cooler travel case, consider the following common solutions:

5.1 Medication Not Staying Cool

- **Ice Packs Not Fully Frozen:** Ensure the ice packs are frozen for the recommended 8-10 hours before use. Partially frozen packs will not provide adequate cooling.
- **Case Not Sealed Properly:** Verify that all zippers are fully closed to maintain the insulated environment. Any gaps can lead to temperature loss.
- **Duration of Cooling:** The case is designed for short-distance travel. For extended periods, additional ice packs or re-freezing may be necessary.
- **External Temperature:** Extreme external temperatures can affect the cooling duration. Keep the case out of direct sunlight and in the coolest available environment.

6. PRODUCT SPECIFICATIONS

Feature	Detail
Model Number	MPN_ISL_Case_Y6_PL_01
Product Dimensions	8.2 x 3.7 x 1.2 inches
Product Weight	3.74 ounces
Materials	High-grade Oxford cloth, Aluminum foil lining
Included Accessories	2 x Ice Packs
Manufacturer	Heartorigin



 Waterproof



TSA approved ice packs



 Cute and stylish



 Dual zipper pulls

Figure 6.1: Product dimensions.

7. WARRANTY AND SUPPORT

For specific warranty information and customer support regarding your YOUSHARES Insulin Cooler Travel Case, please refer to the documentation provided with your purchase or contact the manufacturer, Heartorigin, directly. You may also reach out to Amazon customer service for assistance with your order.

Manufacturer: Heartorigin

For additional information, visit the YOUSHARES Store on Amazon.