

## YOUSHARES 1

# YOUSHARES Insulin Cooler Travel Case Instruction Manual

Model: 1 | Brand: YOUSHARES

## 1. INTRODUCTION

The YOUSHARES Insulin Cooler Travel Case is designed to provide a portable, refrigerated environment for temperature-sensitive medications such as insulin pens, vials, and diabetic supplies. This compact and lightweight case ensures your medical supplies remain cool during travel, camping, or daily activities. It is equipped with insulating materials and reusable ice packs to maintain the required temperature range for several hours.

Image: The YOUSHARES Insulin Cooler Travel Case in a typical usage scenario.

## 2. PRODUCT FEATURES

- **Insulated Design:** Constructed from Oxford cloth with a double-layer aluminum foil and insulating liner to create a thermal barrier.
- **Compact and Lightweight:** Designed for easy portability, fitting into luggage or carry-on bags.
- **Temperature Retention:** Capable of keeping medications cool for approximately 4 to 6 hours with the included ice packs.
- **TSA Approved Ice Packs:** Includes two reusable gel ice packs suitable for air travel.
- **Organized Storage:** Features designated compartments for insulin pens, vials, needles, and other small accessories.
- **Durable Construction:** Equipped with strong zippers and a robust exterior for daily use.

Image: Key features of the YOUSHARES Insulin Cooler Travel Case, including waterproof material, TSA approved ice packs, stylish

design, and dual zippers.

Image: Interior view highlighting the high-quality aluminum foil lining for insulation.

### 3. PACKAGE CONTENTS

Upon opening your YOUSHARES Insulin Cooler Travel Case package, you should find the following items:

- 1 x YOUSHARES Insulin Cooler Travel Case
- 2 x Reusable Gel Ice Packs

Image: Contents of the YOUSHARES Insulin Cooler Travel Case package, including the case and two ice packs.

### 4. SETUP INSTRUCTIONS

1. **Prepare Ice Packs:** Place the two reusable gel ice packs in a freezer for 6 to 8 hours until they are completely frozen. Ensure they are flat to freeze evenly.
2. **Insert Ice Packs:** Once frozen, carefully place the ice packs into the designated compartments within the cooler travel case. The case typically has mesh pockets or elastic straps designed to hold them securely.
3. **Load Medication:** Place your insulin pens, vials, needles, or other diabetic supplies into the main compartment of the case. Utilize the internal elastic loops and mesh pockets to keep items organized and secure.
4. **Close Case:** Zip the case closed securely using the dual zipper pulls to ensure optimal insulation.

Image: Interior view of the case demonstrating how to organize medical supplies.

Image: Example of medical supplies that can be stored in the case.

### 5. OPERATING INSTRUCTIONS

To ensure effective cooling and protection of your medication:

- **Pre-cooling:** Always ensure ice packs are fully frozen before use.
- **Temperature Duration:** The case is designed to keep medication cool for 4 to 6 hours. Actual duration may vary based on external ambient temperature and how frequently the case is opened.
- **Travel Considerations:** The included gel ice packs are TSA approved for air travel. It is recommended to keep the case in your carry-on luggage to maintain optimal temperature and for easy access.
- **Avoid Direct Sunlight:** For best results, keep the cooler case out of direct sunlight and in a cool environment whenever possible.

Image: The cooler case in various travel and outdoor settings.

Image: Demonstrating the convenient portability of the case.

Image: The cooler case being used during outdoor activities.

## 6. MAINTENANCE

- **Cleaning the Case:** Wipe the exterior of the case with a damp cloth and mild soap if necessary. Do not submerge the case in water. Allow it to air dry completely before storing.
- **Cleaning Ice Packs:** The gel ice packs can be wiped clean with a damp cloth. Do not puncture or open the ice packs.
- **Storage:** Store the case and ice packs in a cool, dry place when not in use. Ensure ice packs are completely dry before storing to prevent mold or mildew.

## 7. TROUBLESHOOTING

- **Medication not staying cool long enough:**
  - Ensure ice packs are fully frozen for 6-8 hours prior to use.
  - Minimize opening the case to reduce exposure to ambient air.
  - Keep the case out of direct sunlight and in the coolest available environment.
  - Consider using additional insulation or a larger cooler for extended travel periods beyond 4-6 hours.
- **Zipper sticking:**
  - Apply a small amount of zipper lubricant or a graphite pencil lead to the zipper teeth to improve smoothness.
  - Avoid forcing the zipper; ensure fabric is not caught in the teeth.

## 8. SPECIFICATIONS

Specification	Detail
Product Dimensions	8.27 x 3.94 x 1.58 inches
Item Weight	3.74 ounces
Model Number	1
Manufacturer	Heartorigin
Date First Available	November 8, 2021

Image: Product dimensions of the YOUSHARES Insulin Cooler Travel Case.

Image: Side profile of the cooler case.

## 9. IMPORTANT INFORMATION

### Indications

This insulin cooler travel case is intended for short-distance travel and daily use by individuals who require temperature-controlled storage for their diabetes medication.

### **Directions**

Follow the setup and operating instructions provided in this manual to ensure proper use and temperature maintenance for your medication.

### **Legal Disclaimer**

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

## **10. WARRANTY AND SUPPORT**

---

Information regarding product warranty and customer support is not provided within this instruction manual. Please refer to the product packaging or the manufacturer's official website for details on warranty coverage and how to contact customer support.

