

Chrome Cherry B09PVHCQVJ

Chrome Cherry Nurturally Double Electric Breast Pump User Manual

Model: B09PVHCQVJ

INTRODUCTION

Thank you for choosing the Nurturally Double Electric Breast Pump. This device is designed to support mothers in their breastfeeding journey by providing efficient and comfortable milk expression. Please read this manual thoroughly before first use to ensure proper operation, maintenance, and safety.

IMPORTANT SAFETY INFORMATION

- This product is made from BPA-free and non-toxic silicones, ensuring it is healthy and safe for use.
- Always inspect all pump components before each use for any signs of damage or wear. Replace damaged parts immediately.
- Do not immerse the pump motor unit in water or other liquids.
- Keep the pump and its accessories out of reach of children.
- Use only the provided power adapter and USB cable for charging.
- Discontinue use and consult a healthcare professional if you experience pain or discomfort during pumping.

PACKAGE CONTENTS

Your Nurturally Double Electric Breast Pump package includes the following components:



Image: The Nurturally Double Electric Breast Pump with its control unit, two breast shield assemblies, bottles, tubing, USB cable, and belt clip, all displayed with the product box.



Image: A detailed view of all components: 1x pump unit, 2x breast shield assemblies, 2x milk collection bottles, 1x USB cable, 1x belt clip, 1x tubing, and 4x duckbill valves.

- 1x Pump Motor Unit with Touch Screen Display
- 2x Breast Shield Assemblies (including flanges and massage cushions)
- 2x Milk Collection Bottles with Lids
- 1x USB Charging Cable
- 1x Adhesive Belt Clip
- 1x Tubing Set
- 4x Duckbill Valves

SETUP AND ASSEMBLY

Follow these steps to assemble your breast pump:

1. Ensure all components are clean and dry before assembly.
2. Screw the breast shield assembly onto the milk collection bottle.
3. Insert a duckbill valve into the bottom of the breast shield assembly.

4. Connect one end of the tubing to the breast shield assembly and the other end to the pump motor unit. If using both breast shields, connect both tubes to the pump unit.
5. To use the pump hands-free, attach the adhesive belt clip to the back of the pump motor unit.

Video: An overview demonstrating the assembly of the Nurturally Portable Electric Breast Pump components.



Image: A user demonstrating the hands-free capability of the Nurturally pump, attached to a pumping bra while working.

Nurturally

Portable Electric Breast Pump



Two Pumping
Modes -
Dos Extracción
Modos



Touch Screen
Display -
Pantalla Táctil

Mode Level Indicator

Indicador de nivel de modo

Mode Icon

Icono de modo

Mode Level Changing

Cambio de nivel de modo

Battery Level Indicator

Indicador de nivel
de batería

Mode Changing

Cambio de modo



Image: A detailed view of the pump's touch screen interface showing mode and intensity indicators, alongside the detachable belt clip for portability.

CHARGING THE PUMP

The Nurturally breast pump features a built-in USB rechargeable battery. To charge the device:

1. Connect the USB charging cable to the pump motor unit.
2. Plug the other end of the USB cable into a compatible USB power source (e.g., wall adapter, computer, power bank).
3. The battery indicator on the touch screen display will show the charging status. A full charge can last up to 10 pumping sessions.
4. The pump can be used while it is charging.

Nurturally

Portable Electric Breast Pump

Rechargeable
Battery -
Batería recargable



The Battery
operates up to 10 sessions
on a single full charge

La batería funciona hasta 10 sesiones en un
sola carga completa

Image: The Nurturally pump unit being charged via a USB cable connected to a laptop, highlighting its rechargeable feature.

OPERATING THE PUMP

The pump offers two distinct modes and multiple intensity levels for personalized comfort and efficiency.

Power On/Off

- Press and hold the power button (circle icon with a line) on the touch screen display to turn the pump on or off.

Pumping Modes

The pump features two modes:

1. **Stimulation Mode:** This mode mimics the fast, fluttering sucking pattern of a baby to help stimulate milk flow.
2. **Expression Mode:** This mode mimics the slower, steady suck-and-swallow pattern of a baby

during active milk flow.

Switching Modes and Adjusting Intensity

- Once powered on, the pump will typically start in Stimulation mode.
- Press the mode change button (circular arrow icon) to switch between Stimulation and Expression modes.
- Use the '+' and '-' buttons to adjust the intensity level within each mode. The touch screen display will show the current mode and intensity level.
- Experiment with different levels to find what is most comfortable and effective for you.

MAINTENANCE AND CARE

Proper cleaning and maintenance are crucial for the longevity and hygiene of your breast pump.

Cleaning Instructions

- **Before First Use and After Each Use:** Disassemble all parts that come into contact with breast milk (breast shields, bottles, valves, tubing). Wash these parts in warm, soapy water. Rinse thoroughly with clean water.
- **Sterilization:** Sterilize milk-contacting parts regularly according to manufacturer guidelines for baby feeding equipment (e.g., boiling, steam sterilizer).
- **Motor Unit:** Wipe the pump motor unit with a clean, damp cloth. Do not immerse it in water.
- Ensure all parts are completely dry before reassembling or storing to prevent mold or bacterial growth.

Naturally

Portable Electric Breast Pump

Massage Cushions

Made of soft silicone, it reduces pain and stimulates when expressing

Cojines de masaje de suave silicona, reduce el dolor y estimula la extracción



Food Grade Materials

Made of BPA free and non-toxic silicones, safe and healthy for babies

El producto está hecho de siliconas libres de BPA y no tóxicas. Son materiales de calidad alimentaria seguros y saludables para los bebé

Image: The breast shield and bottle, emphasizing the soft silicone massage cushions and the use of BPA-free, non-toxic, food-grade materials.

TROUBLESHOOTING

If you encounter issues with your Naturally Double Electric Breast Pump, refer to the following common problems and solutions:

- **No Suction or Weak Suction:**

- Ensure all pump components are correctly assembled and securely connected.
- Check that the duckbill valves are properly seated and not damaged or torn.
- Verify that the tubing is securely attached to both the breast shield assembly and the motor unit, with no kinks or blockages.
- Confirm that the breast shield forms a tight seal around your breast.
- Increase the suction intensity using the '+' button on the control panel.

- **Pump Not Turning On:**

- Ensure the pump is adequately charged. Connect it to a power source using the USB cable.
- Press and hold the power button for a few seconds to activate the device.

- **Unusual Noise:**

- While the pump is designed to be quiet, any unusual or loud noises may indicate a problem. Check for proper assembly and ensure no parts are loose.
- If the noise persists, contact customer support.

SPECIFICATIONS

Feature	Detail
Brand	Chrome Cherry
Model	B09PVHCQVJ
Operation Mode	Electric
Material	Silicone
Power Source	1 Lithium Ion battery (included), USB rechargeable
Battery Life	Up to 10 sessions on a single full charge
Dimensions	9.1 x 7.7 x 3.5 inches (Package)
Weight	1.5 Pounds (Package)
Manufacturer	Guangdong Bruce Technology Co., Ltd.

WARRANTY AND SUPPORT

For warranty information or technical support, please refer to the product packaging or contact Chrome Cherry customer service directly. Keep your purchase receipt as proof of purchase.