

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Golden Technologies](#) /
- › [Golden Technologies Buzzaround Carry On Mobility Scooter User Manual](#)

Golden Technologies GB120

Golden Technologies Buzzaround Carry On Mobility Scooter User Manual

Model: GB120

INTRODUCTION

This manual provides essential information for the safe and effective operation, maintenance, and care of your Golden Technologies Buzzaround Carry On Fold-Flat Travel Lithium Battery Lightweight Mobility Scooter. Please read this manual thoroughly before using your scooter to ensure proper understanding of its features and functions.

KEY FEATURES

- **Fold-Flat Design:** Easily folds for compact storage and transport, with a folded height of only 12 inches.
- **Lightweight Construction:** Heaviest piece weighs 47.5 lbs (without seat and batteries), facilitating easy handling.
- **Lithium-Ion Battery:** Provides an impressive drive range of approximately 18 miles on a single charge.
- **High Weight Capacity:** Supports users up to 300 lbs.
- **Full Lighting Package:** Equipped with front and rear lighting for enhanced visibility and safety.
- **Adjustable Tiller:** Features an infinitely adjustable tiller for personalized comfort.
- **Flip-Back Armrests:** Allows for easy entry and exit from the seat.
- **Convenience Features:** Includes a horn, front basket for personal items, and a drink holder.

SETUP

Unpacking and Initial Assembly

Your Buzzaround Carry On scooter is designed for easy assembly. It disassembles into two main pieces for transport. Upon receiving, carefully remove all components from the packaging.

1. Place the main base unit on a flat, stable surface.
2. Unfold the tiller assembly until it locks into the upright position.
3. Attach the seat by aligning it with the seat post and securing it firmly. Ensure the seat is fully engaged and stable.
4. Attach the front basket to the designated mounting point on the tiller.



Image: The Golden Technologies Buzzaround Carry On Mobility Scooter in its fully assembled state, showcasing its compact design and front basket.

Battery Charging

The scooter is equipped with a 15AH Lithium-Ion battery pack. Before first use, fully charge the battery.

- Locate the charging port on the tiller.
- Connect the off-board battery charger to the charging port and then to a standard electrical outlet.
- Allow the battery to charge until the indicator light on the charger signals a full charge. This may take several hours.
- Always use the charger provided by Golden Technologies.

OPERATING INSTRUCTIONS

Starting and Driving

1. Ensure the battery is fully charged and properly installed.
2. Insert the key into the ignition and turn it to the 'ON' position.
3. Adjust the tiller to a comfortable position using the infinite adjustable tiller mechanism.
4. To move forward, gently push the throttle lever forward. To reverse, gently pull the throttle lever backward.
5. The maximum speed is 3.60 mph. Adjust speed using the speed control dial on the tiller.
6. To stop, release the throttle lever. The electro-mechanical braking system will engage automatically.

Maneuvering and Safety

- The scooter has a turning radius of 48.75 inches. Practice in an open area to familiarize yourself with its maneuverability.
- Utilize the horn to alert others when necessary.

- Ensure the front and rear lights are on when operating in low-light conditions.
- The flip-back armrests allow for easier transfer onto and off the scooter. Simply lift the armrests upwards.



Image: Side view of the Golden Technologies Buzzaround Carry On Mobility Scooter, highlighting its compact frame and four-wheel design.

MAINTENANCE

Regular Care

- **Cleaning:** Wipe down the scooter regularly with a damp cloth. Avoid using abrasive cleaners or excessive water.
- **Battery Care:** Always store the battery in a charged state. If storing for an extended period, charge it at least once a month.
- **Tire Inspection:** Periodically check the tires for wear and proper inflation (if applicable, though these are solid tires).
- **General Inspection:** Regularly check all nuts, bolts, and fasteners to ensure they are secure.

Disassembly for Transport/Storage

The scooter is designed for easy disassembly into two pieces for transport or compact storage, such as under a bed (folded height 12 inches).

1. Remove the seat by lifting it off its post.
2. Fold down the tiller assembly.
3. Separate the main base unit into its two primary components as per the quick-release mechanism.



Image: The Golden Technologies Buzzaround Carry On Mobility Scooter in its folded configuration, ready for compact storage or transport.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Scooter does not turn on.	Battery not charged, key not fully inserted, power switch off.	Charge battery, ensure key is fully inserted and turned to 'ON', check power switch.
Reduced range or speed.	Low battery charge, worn battery, rough terrain, excessive weight.	Fully charge battery, consider battery replacement if old, operate on smoother surfaces, ensure weight is within capacity.
Unusual noises during operation.	Loose components, motor issue.	Check for loose parts and tighten. If noise persists, contact customer support.

For issues not listed or if solutions do not resolve the problem, please contact Golden Technologies customer support or refer to the detailed troubleshooting guide in your product's full manual.

SPECIFICATIONS

Feature	Detail
Model Number	GB120
Brand	Golden Technologies
Top Speed	3.60 mph
Turning Radius	48.75 inches
Drive Range	18.00 miles
Weight Capacity	300 lbs
Heaviest Piece Weight (without seat and batteries)	47.5 lbs
Folded Dimensions (D x W x H)	12"D x 19.5"W x 12"H
Overall Length	38 inches
Overall Width	19.5 inches
Ground Clearance	4.00 inches

Feature	Detail
Battery Type	15AH Lithium-Ion (included)
Wheel Type	Four-Wheeled
Front Wheel Size	6 inches
Rear Wheel Size	7 inches
Seat Width	15 inches
Seat Depth	13.75 inches
Max Seat to Ground Height	22.75 inches
Max Seat to Deck Height	19.5 inches
Drive Type	Rear-Wheel Drive
Braking System	Electro-Mechanical
Adjustable Tiller	Yes
Flip-Back Armrests	Yes
Headlight	Yes
Basket Included	Yes



Image: Diagram showing key dimensions of the Golden Technologies Buzzaround Carry On Mobility Scooter.

WARRANTY & SUPPORT

Your Golden Technologies Buzzaround Carry On Mobility Scooter comes with a **1-Year In-Home Service contract** at no additional cost. This service provides professional assistance for any covered issues directly at your home.

For warranty claims, technical support, or service inquiries, please contact Golden Technologies customer service directly. Refer to the contact information provided in your product packaging or on the official Golden Technologies website.

Always ensure to register your product to activate your warranty and receive important updates.



