



Manuals.plus /

› AudioNova /

› AudioNova Clean & Dry Box 3.2 Hearing Aid Dryer and UV Sterilizer User Manual

AudioNova TC-CDB-20

AudioNova Clean & Dry Box 3.2 User Manual

Model: TC-CDB-20

Brand: AudioNova

1. INTRODUCTION

The AudioNova Clean & Dry Box 3.2 is an advanced device designed to provide optimal care for your hearing aids. It effectively dries and sterilizes your devices, protecting sensitive electronic components from moisture damage, earwax buildup, and harmful microorganisms. Regular use helps to significantly extend the lifespan and maintain the performance of your hearing aids.



A view of the AudioNova Clean & Dry Box 3.2, a compact device designed for drying and sterilizing hearing aids. The lid is open, revealing the internal compartment where hearing aids are placed for treatment.

2. SAFETY INFORMATION

Please read all safety instructions carefully before using the AudioNova Clean & Dry Box 3.2 to ensure safe and effective operation.

- **General Hygiene:** Always ensure your hands are clean and dry before handling your hearing aids or the Clean & Dry Box.
- **Protection Against Humidity:** Remove your hearing aid before showering or swimming. If it gets wet, dry it only in the designated drying box. Never use external heat sources like a hairdryer or radiator.
- **Avoid Shocks:** Be careful not to drop your hearing aid. Even a fall from a short height can damage the device.
- **Protection Against Heat:** High temperatures can damage the hearing aid. Do not expose it to strong heat, such as in a sauna or direct sunlight.
- **Power Supply:** Use only the provided micro-USB cable and charger adapter. Do not use damaged cables or adapters.

- **Children and Pets:** Keep the device out of reach of children and pets.



Hygiène

L'hygiène de l'oreille est importante pour une audition parfaite. Avant de toucher votre appareil auditif, assurez-vous que vos doigts soient propres et secs.

Protection contre l'humidité

Retirez votre aide auditive avant de prendre une douche ou de nager. Si elle est mouillée, séchez-la uniquement dans la boîte de séchage prévue à cet effet. N'utilisez jamais de sources de chaleur telles qu'un sèche-cheveux ou un radiateur.

Éviter les chocs

Veillez à ne pas faire tomber votre aide auditive. Même un impact à faible hauteur peut endommager l'appareil.

Protection contre la chaleur

Les températures élevées endommagent l'appareil auditif. Ne l'exposez jamais à une forte chaleur, par exemple dans un sauna ou en plein air, aux rayons du soleil.

An informational graphic providing essential care tips for hearing aids, covering hygiene practices, protection from humidity, advice on avoiding physical shocks, and precautions against high temperatures.

3. PACKAGE CONTENTS

Upon receiving your AudioNova Clean & Dry Box 3.2, please ensure all the following items are included:

- 1x AudioNova Clean & Dry Box 3.2 (Drying box for hearing aids)
- 1x Micro-USB cable
- 1x Charger adapter
- 1x User Manual (this document)



The AudioNova Clean & Dry Box 3.2 displayed alongside its retail packaging, highlighting the product's design and branding.

4. PRODUCT FEATURES

- **Effective Drying:** The Drybox dries your hearing aids at 30 degrees Celsius, protecting sensitive electronic components from damage caused by moisture and earwax, thereby extending their lifespan.
- **UV-C Sterilization:** Features a UV-C program for hygienic cleaning, eliminating 99.9% of harmful microorganisms such as bacteria, fungi, and viruses. This prevents ear infections and protects both your hearing aids and your health.
- **Easy & Automatic Use:** Simple one-button operation initiates a complete 5-hour disinfection and drying cycle.
- **Wide Compatibility:** Compatible with all in-ear and behind-the-ear hearing aids, including rechargeable models. It can also be used for reusable earplugs and small jewelry items that require cleaning.
- **Dual Device Capacity:** This box can clean and dry two devices simultaneously.



Affichage de l'état



USB-C + Prise EU



Séchante



UV-C



1 cycle de 5 heures



An image detailing the key features of the AudioNova Clean & Dry Box 3.2, including a status display, USB-C connectivity, drying function, UV-C sterilization, and a 5-hour operating cycle. A person wearing a hearing aid is shown in the background.



The AudioNova Clean & Dry Box 3.2 with its removable drying tray lifted, illustrating how hearing aids are placed inside for drying and UV-C sterilization.

5. SETUP

1. Unpack the AudioNova Clean & Dry Box 3.2 and all accessories.
2. Place the device on a stable, flat surface away from direct sunlight and heat sources.
3. Connect the micro-USB cable to the port on the back of the Clean & Dry Box.
4. Plug the charger adapter into a standard electrical outlet.
5. The device is now ready for use.

6. OPERATING INSTRUCTIONS

Follow these steps to dry and sterilize your hearing aids:

1. Open the lid of the Clean & Dry Box.
2. If your hearing aids use disposable batteries, open the battery compartment to allow for better drying. For rechargeable hearing aids, ensure they are clean and ready for the drying process.

3. Place your hearing aids into the drying compartment. The device can accommodate up to two hearing aids simultaneously.
4. Close the lid securely.
5. Press the Start button located on the device. A light indicator will illuminate, signifying that the drying and disinfection process has begun.
6. The device will automatically complete a 5-hour cycle. Once the process is finished, the light indicator will turn off, and the device will power down.
7. Open the lid and remove your clean, dry, and sterilized hearing aids.



SÉCHER ET DÉSINFECTER AVEC LA DRY BOX

Placez l'aide auditive dans l'appareil de séchage avec le compartiment de la pile ouvert. Appuyez sur le bouton de démarrage. Le témoin lumineux s'allume. Une fois le processus de séchage terminé le boîtier s'éteint automatiquement.

The AudioNova Clean & Dry Box 3.2 with its lid open, illustrating the process of placing a hearing aid inside for drying and disinfection. The accompanying text explains to place the hearing aid with the battery compartment open and press the start button.

7. MAINTENANCE

Proper maintenance ensures the longevity and effectiveness of your hearing aids and the Clean & Dry Box.

- **Daily Cleaning of Hearing Aids:** Daily cleaning is crucial. Use specially designed cleaning cloths or sprays. Never use alcohol, solvents, or other harsh cleaning products on your hearing aids.
- **Battery Handling:** Whenever you are not wearing your hearing aids, remove the batteries to protect your device from potential damage.
- **Safe Storage:** Keep your hearing aid out of reach of children and animals. Always transport it in its case. When not using your hearing aid, the battery compartment should be open (for non-rechargeable models).
- **Sterilization:** Microorganisms are part of our daily lives. Special formulas in cleaning and maintenance products protect against impurities and effectively prevent external ear infections.
- **Cleaning the Dry Box:** Wipe the interior and exterior of the Clean & Dry Box with a soft, dry cloth as needed. Do not use abrasive cleaners or immerse the device in water.



Un nettoyage quotidien

Un nettoyage en quotidien est particulièrement important. Utilisez des chiffons spécialement prévus à cet effet ou le spray de nettoyage, jamais d'alcool, de solvants ou d'autres produits de nettoyage.



Les piles auditives!

Chaque fois que vous ne portez pas vos appareils auditifs, retirez les piles afin de protéger votre appareil auditif d'un éventuel endommagement



Rangement en toute sécurité

Gardez l'aide auditive hors de portée des enfants et des animaux. Transportez-la toujours dans son étui. Lorsque vous n'utilisez pas votre aide auditive, le compartiment de la pile doit être ouvert



Stérilisation

Les micro-organismes nous accompagnent au quotidien, certains sont indésirables. C'est aussi le cas de l'appareil auditif et de son étui. Des formules spéciales dans les produits de nettoyage et d'entretien protègent des impuretés et préviennent efficacement les infections de votre oreille externe.

An informational graphic outlining key maintenance practices for hearing aids, including daily cleaning routines, proper handling of hearing aid batteries, secure storage methods, and the importance of sterilization.

8. TROUBLESHOOTING

If you encounter any issues with your AudioNova Clean & Dry Box 3.2, please refer to the following common solutions:

• Device Not Turning On:

- Ensure the micro-USB cable is securely connected to both the device and the power adapter.
- Verify that the power adapter is properly plugged into a working electrical outlet.
- Try using a different USB port or power outlet.

• Drying/UV Cycle Not Starting:

- Make sure the lid is fully closed before pressing the Start button. The device has a safety mechanism that prevents operation if the lid is open.
- Press the Start button firmly.

• Hearing Aids Still Damp After Cycle:

- Ensure the battery compartments of your hearing aids (if applicable) were open during the drying cycle.
- Avoid placing excessively wet hearing aids into the box. Gently wipe off any visible moisture before placing them inside.
- Ensure the drying compartment is not overcrowded, allowing for proper air circulation.

• UV Light Not Visible:

- The UV-C light is designed to operate safely within the closed box. It may not be visible from the outside. If the drying cycle starts, the UV-C function is likely active.

If the problem persists after attempting these solutions, please contact customer support.

9. SPECIFICATIONS

| | |
|-----------------------------|--|
| Model Number | TC-CDB-20 |
| Brand | AudioNova |
| Manufacturer | Wake up International B.V. |
| Drying Temperature | Approximately 30°C |
| Sterilization Method | UV-C Light |
| Cycle Duration | 5 hours (automatic) |
| Power Input | Micro-USB |
| Compatibility | All in-ear and behind-the-ear hearing aids (including rechargeable), reusable earplugs, small jewelry. |

10. WARRANTY INFORMATION

Specific warranty terms and conditions for the AudioNova Clean & Dry Box 3.2 may vary by region and retailer. Please refer to the warranty card included with your product or contact your point of purchase for detailed information regarding warranty coverage and claims.

11. SUPPORT

For further assistance, technical support, or inquiries regarding your AudioNova Clean & Dry Box 3.2, please visit the official AudioNova website or contact their customer service department. Contact details can typically be found on the product packaging or the manufacturer's website.