

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Joylux](#) /
- › [Joylux vFit Device, Photonic Gel & Wipes Set Instruction Manual](#)

Joylux BL-3010-05

Joylux vFit Device, Photonic Gel & Wipes Set Instruction Manual

Model: BL-3010-05

INTRODUCTION

This manual provides instructions for the safe and effective use of the Joylux vFit device, Photonic Gel, and Photonic Wipes. This set is designed to support feminine wellness through a combination of red light technology, gentle warmth, and sonic pulsation, complemented by specialized gel and cleaning wipes.

SAFETY INFORMATION

- Not for use on the body.
- Read all terms and conditions before use.
- For adult use only.

WHAT'S INCLUDED

Your Joylux vFit 3-Piece Set includes:

- One (1) vFit Device
- One (1) Joylux Photonic Gel (2.5 fl oz / 75ml)
- One (1) Pack of Joylux Photonic Wipes (30 count)



Image: The complete Joylux vFit device, Photonic Gel, and Photonic Wipes set.

SETUP

Charging the vFit Device

Before first use, fully charge your vFit device. Refer to the device's specific charging instructions for details.

App Installation and Pairing

Download the Joylux app from your device's app store. Follow the in-app instructions to pair your vFit device. The app provides personalized plans and tracks usage.

Video: An overview of the Joylux vFit Gold device, highlighting its features and benefits for intimate health, including its non-invasive nature, red-light technology, and app integration.

OPERATING THE vFIT DEVICE

Application of Photonic Gel

Apply a small amount of Joylux Photonic Gel to the surface of the vFit device before each use. This gel is pH-balanced and designed to enhance comfort and device performance.

PHOTONIC GEL



MOISTURIZING

Promotes lasting hydration



USE WITH vFIT™

Designed for use with vFit™ Intimate Wellness device



PERFORMANCE

Enhances performance of vFit by up to 45%



COMFORT

Enhances comfort of the vFit device



Image: A tube of Joylux Photonic Gel, emphasizing its moisturizing properties and use with the vFit device.

Device Usage

Insert the vFit device as instructed by the Joylux app. The device utilizes red light technology, gentle warmth, and sonic pulsation.

Follow the personalized session duration and frequency recommended by the Joylux app. The app will guide you through each session.



Image: The Joylux vFit device, illustrating its natural, easy, quick, and effective features.

vFIT™

Easy-To-Use App Tracker



Calendar with alerts
& reminders



Helpful educational
content



Daily health check-
ins & reports



Image: The Joylux app tracker interface on a smartphone, showing calendar, educational content, and daily health check-ins.

Video: Testimonials from users of the Joylux vFit device, discussing their experiences with the red light menopause device.

MAINTENANCE

Cleaning the vFit Device

After each use, clean the vFit device thoroughly with Joylux Photonic Wipes. These wipes are specifically formulated to be safe for the device's silicone and electronic components, preventing damage from harsh household cleaners.

PHOTONIC WIPES

MILD FORMULA

- ✓ 99.5% effective at removing the most common germs
- ✓ Protects the silicone, lights and electronics from chemical & water damage
- ✓ Safe and easy to use



Image: Joylux Photonic Wipes, highlighting their mild formula and effectiveness for device cleaning.

Video: A user demonstrating the proper use of vFit Photonic Wipes for cleaning the device, emphasizing their importance for device care.

SPECIFICATIONS

Feature	Detail
Product Dimensions	10 x 4 x 4 inches; 10.58 ounces
Item Model Number	BL-3010-05
Batteries	1 Lithium Ion batteries required (included)
UPC	850012528197
Manufacturer	Joylux
Skin Type	Sensitive (for associated products)

Feature	Detail
Material Feature	Natural (for associated products)
Unit Count	1.0 Count (for the device)
Number of Items	3 (for the set)

TROUBLESHOOTING

If you experience any issues with your vFit device, Photonic Gel, or Photonic Wipes, please refer to the Joylux app's support section or contact Joylux customer service.

- **Device Not Turning On:** Ensure the device is fully charged.
- **App Connectivity Issues:** Verify Bluetooth is enabled on your smartphone and the device is within range. Restart both the app and the device if issues persist.
- **Discomfort During Use:** Ensure adequate application of Photonic Gel. If discomfort continues, discontinue use and consult a healthcare professional.

WARRANTY AND SUPPORT

For warranty information and customer support, please visit the official Joylux website or contact their customer service directly.

The product has a 30-day return/replacement policy from the date of purchase.