

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [MySmile](#) /

› [MySmile Teeth Whitening Kit with LED Light User Manual](#)

MySmile B0B46GGJ9B

MySmile Teeth Whitening Kit with LED Light User Manual

Model: B0B46GGJ9B

1. PRODUCT OVERVIEW

The MySmile Teeth Whitening Kit with LED Light is designed for at-home teeth whitening. It utilizes a non-sensitive formula with Carbamide Peroxide gel and an LED light accelerator to help remove stains and whiten teeth effectively.

The advertisement features a woman on the left with a bright smile. In the upper right, two circular icons show a tooth labeled 'before' (yellowish) and a tooth labeled 'after' (bright white). Below these icons, a white arrow points to the text 'YES, IT WORKS!'. Another white arrow points down to a white LED whitening device with a clear U-shaped tray, which is being held in a hand. The device has 'MySmile' written on it. In the bottom left, a circular graphic contains the text 'UP TO 10 SHADES WHITER AFTER 9 APPLICATIONS'. The brand name 'MySmile' is written in the bottom right corner.

Image 1.1: The MySmile Teeth Whitening Kit, illustrating potential before and after whitening results with the LED device.

Powerful Whitening Formula



Image 1.2: Close-up of the MySmile whitening gel syringe, highlighting its non-sensitive, enamel-safe, and effective whitening formula.

2. PACKAGE CONTENTS

- LED Light Device
- Mouth Tray
- Teeth Whitening Gel Syringes (3 units)

3. SETUP AND PREPARATION

1. **Charge the LED Device:** Ensure your LED light device is fully charged before first use.

2. **Attach Mouth Tray:** Securely attach the clear mouth tray to the LED light device.
3. **Prepare Gel Syringe:** Remove the cap from a whitening gel syringe.



Image 3.1: Demonstrates the process of applying the whitening gel from the syringe into the mouth tray.

4. OPERATING INSTRUCTIONS

Follow these steps for an effective teeth whitening session:

1. **Apply Gel:** Carefully apply a thin, even layer of whitening gel along the top and bottom sections of the mouth tray where it will contact your teeth.
2. **Insert Device:** Place the mouth tray into your mouth, ensuring your teeth are fully covered by the gel.
3. **Activate LED Light:** Press the power button on the LED light device to turn it on. The blue light will

illuminate.

4. **Whitening Session:** Keep the device in your mouth for **10 minutes**. The device features an automatic timer and will beep and turn off once the session is complete.
5. **Remove and Rinse:** After 10 minutes, remove the device from your mouth. Rinse your mouth and the device thoroughly with water.

Your browser does not support the video tag.

Video 4.1: This video demonstrates the application of the whitening gel into the mouth tray and the subsequent use of the LED light device for a whitening session.



Image 4.1: A user demonstrating the application of the MySmile LED teeth whitening device for a 10-minute fast whitening result.

5. MAINTENANCE AND CARE

Proper maintenance ensures the longevity and hygiene of your MySmile Teeth Whitening Kit:

- **Cleaning the Mouth Tray:** After each use, detach the mouth tray from the LED light device. Rinse the mouth tray thoroughly under cool running water to remove any residual gel.
- **Cleaning the LED Device:** Gently wipe the LED light device with a damp cloth. **Do not submerge**

the LED device in water to prevent damage to the electronic components.

- **Storage:** Store all components of the kit in a clean, dry, and cool place away from direct sunlight.

6. SAFETY INFORMATION AND WARNINGS

Please read and adhere to the following safety guidelines:

- **Avoid Ingestion:** Do not ingest the whitening gel. If swallowed, seek medical attention immediately.
- **Children:** This product is not intended for use by children. Keep out of reach of children.
- **Sensitivity:** The product may cause temporary tooth or gum sensitivity. If sensitivity persists or worsens, discontinue use and consult a dental professional.
- **Proper Use:** Always follow the directions for safe and proper use as outlined in this manual.
- **Medical Disclaimer:** Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

7. SPECIFICATIONS

Feature	Detail
Brand Name	MySmile
Product Benefits	Whitening
Item Form	Pen (for gel applicator)
Included Components	Teeth Whitening Gel
Age Range Description	Adult
Set Name	7 Piece Set
Number of Items	2 (Kit components)
Active Ingredients	Carbamide peroxide
Material Type Free	Gluten Free
Material Features	Gluten Free
Flavor	Unflavored
ASIN	B0B46GGJ9B
Date First Available	June 15, 2022

8. TROUBLESHOOTING

If you encounter any issues with your MySmile Teeth Whitening Kit, please consider the following general advice:

- **Device Not Turning On:** Ensure the LED light device is fully charged.
- **No Whitening Effect:** Verify that the gel is applied correctly and evenly to your teeth. Ensure the LED light is active during the entire 10-minute session. Consistent use as directed is key for visible

results.

- **Increased Sensitivity:** If you experience discomfort, reduce the frequency of use or discontinue temporarily. Consult a dental professional if symptoms persist.

For further assistance, please refer to the contact information provided in the Warranty & Support section.

9. WARRANTY AND SUPPORT

Specific warranty details for the MySmile Teeth Whitening Kit are not provided in the product information. For warranty claims, technical support, or customer service inquiries, please contact MySmile directly through their official website or the retailer where the product was purchased.

You can often find support information on the product packaging or by visiting the [MySmile Store on Amazon](#).