

Dynarex 10602-7

Dynarex Ever-Soft Air Cushion User Manual

Model: 10602-7

Brand: Dynarex

1. INTRODUCTION

The Dynarex Ever-Soft Air Cushion is engineered to provide superior comfort and pressure redistribution for individuals requiring extended periods of sitting. Its innovative air cell design helps prevent and alleviate pressure wounds, enhancing overall seated comfort and posture.



Figure 1.1: Dynarex Ever-Soft Air Cushion (Pack of 7) and packaging.

2. SETUP INSTRUCTIONS

2.1. Cushion Inflation

The Ever-Soft Air Cushion comes with a convenient hand pump for inflation. Locate the inflation valve on the cushion. Insert the pump nozzle into the valve and inflate the cushion to your desired firmness. The independent, interconnected air cells will redistribute air to provide a gentler and more supportive seat. Do not overinflate.

HAND PUMP INCLUDED

Included for quick inflation of the air cushion for transport and standard wheelchairs.



WE CARE LIKE FAMILY™

Figure 2.1: Hand pump for quick inflation of the air cushion.



Figure 2.2: The cushion includes a hand pump and a repair kit.

2.2. Cover Attachment

The cushion cover is made of heavy-duty fabric with a non-skid bottom. Place the cushion inside the cover. The cover features adjustable locking straps designed to keep the cushion securely in place on a wheelchair or

other seating surface. Fasten these straps around the frame of the chair to prevent shifting.

CUSHION COVER

The cushion cover is a heavy-duty fabric with a non-skid bottom, can be removed, and has adjustable locking straps to keep the cushion in place.



WE CARE LIKE FAMILY™

Figure 2.3: Cushion cover with adjustable straps for secure placement on a wheelchair.

3. OPERATING INSTRUCTIONS

The Dynarex Ever-Soft Air Cushion is versatile and can be used in various settings to increase seated comfort and prevent pressure wounds. Place the inflated cushion on the desired seating surface, such as a wheelchair, car seat, or other chairs within the home. Ensure the non-skid bottom of the cover is facing down to maintain stability. Adjust the air fill amount as needed for customized comfort and support.

AIR CUSHION DESIGN

The Dynarex Ever-Soft Air Cushion is designed with independent, interconnected air cells.



Figure 3.1: The air cushion features independent, interconnected air cells for pressure redistribution.

4. MAINTENANCE

4.1. Cushion Cover Care

The Ever-Soft Air Cushion cover is machine washable. Remove the cover from the cushion before washing. Follow the care instructions on the cover's label for best results. Allow the cover to air dry completely before reattaching it to the cushion.

5. TROUBLESHOOTING

No specific troubleshooting steps are provided for this product. If you experience issues with your Dynarex Ever-Soft Air Cushion, please refer to the contact information in the Support section.

6. SPECIFICATIONS

Product Dimensions	18 x 18 x 2 inches
Item Model Number	10602-7
Weight	2 Pounds
Weight Capacity	300 lb
Color	Black
Shape	Rectangular
Manufacturer	Dynarex Corporation
Country of Origin	China



Figure 6.1: Dynarex Ever-Soft Air Cushion, Model 10602-7.

7. WARRANTY AND SUPPORT

Specific warranty information for the Dynarex Ever-Soft Air Cushion is not provided in this manual. For details regarding warranty coverage or to contact customer support, please visit the official Dynarex website or refer to the product packaging.

You can also visit the [Dynarex Store on Amazon](#) for more information and product inquiries.