Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [SENSSUN](#) /

› [SENSSUN Smart Body Fat Scale User Manual](#)


SENSSUN IF1090D

SENSSUN Smart Body Fat Scale

Model: IF1090D

INTRODUCTION

The SENSSUN Smart Body Fat Scale (Model: IF1090D) is a high-precision digital personal scale designed to measure body weight and analyze various body composition metrics. Utilizing Bioelectric Impedance Analysis (BIA) technology, it provides comprehensive data to help you track your health and fitness progress. This manual provides essential information for the proper setup, operation, and maintenance of your scale.

The SENSSUN Smart Body Fat Scale, displaying a weight of 68.8 kg on its LCD screen.



IMPORTANT SAFETY INSTRUCTIONS

Please read all instructions carefully before use and save them for future reference.

- The scale is designed for body weight measurement and body fat analysis. Do not use it for any other purpose.
- Do not use this scale if you are pregnant. Do not use this scale if you have a pacemaker or other internal medical device.
- Always place the scale on a hard, flat surface. Do not use on carpet or soft surfaces, as this may affect accuracy.
- The glass platform is slippery when wet. Ensure the platform is dry before stepping on.
- The scale is a high-precision measuring device. Never jump or stomp on the scale.
- Do not use the scale if it is damaged. Do not attempt to disassemble or repair it.
- Keep the scale out of reach of children.
- This scale is for household use only. Do not use for commercial or medical purposes.
- Do not expose the scale to extreme temperatures or humidity.

Battery Safety

- **User-Friendly Design:** Automatic calibration, auto on/off function, clear LCD display, low battery indicator, and overload indicator.
- **Durable Construction:** Built with high-performance explosion-proof glass.
- **Capacity:** Supports up to 150 kg / 330 lb with 0.1 kg / 0.2 lb increments.

The scale provides 15 essential body composition measurements for comprehensive health tracking.

Multifunctional Device

Key features of the SENSUN Smart Body Fat Scale, highlighting its robust design and smart functionalities.

SETUP

1. Battery Installation

1. Open the battery compartment on the back of the scale.
2. Insert 2 x AAA batteries, ensuring correct polarity (+ and -).
3. Close the battery compartment securely.




2. App Download and Setup

The SENSUN Smart Body Fat Scale works with the free MovingLife app. Follow these steps to download and set up the app:

1. Scan the QR code on the product packaging or visit http://www.camry.com.cn/app_download_movinglife.html to download the MovingLife app from the App Store (iOS) or Google Play Store (Android).

2. Install the app on your smartphone.

3. Open the MovingLife app and create a user profile. Enter your personal information such as gender, age, height, and activity level for accurate body composition analysis.

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Instructions for downloading and setting up the MovingLife app.

3. Connecting to the App (Pairing)

1. Ensure Bluetooth is enabled on your smartphone.
2. Open the MovingLife app.
3. Step onto the scale barefoot. The scale will automatically turn on and display your weight.
4. The app will automatically detect and connect to the scale. Once connected, your measurements will be synchronized to the app.
5. If the connection fails, ensure the app is open and try stepping on the scale again.

Taking Measurements

1. Place the scale on a hard, flat, and stable surface.
2. Step onto the scale barefoot, ensuring your feet are placed on the four metal electrodes.
3. Stand still until your weight is displayed and stabilized on the LCD screen.
4. If the app is open and connected, all body composition data will be automatically transferred and recorded.
5. You can weigh yourself at any time without your phone. The data will synchronize with the app as soon as they are connected again.


Using the MovingLife App

The MovingLife app allows you to track your progress, view detailed body composition reports, and set goals.

- **Data Tracking:** View daily, weekly, monthly, and yearly trends for all 15 body metrics.
- **Goal Setting:** Set personal health and fitness goals within the app.
- **User Management:** The app supports multiple user profiles, making it suitable for families.
- **Data Synchronization:** Syncs data with Apple Health, Google Fit, and FitBit for integrated health management.

Track Your Progress

You can know all the data changes from the app on a daily or monthly, and set goals to further manage your fitness plan.

★ Our Scale
(Data are measured by contacting the bioelectrode with your foot.)

Others
(Only the weight can be measured, other data are obtained by the algorithm.)

Track your progress over time with detailed graphs and reports in the MovingLife app.


Baby Mode

The scale features a convenient Baby Mode, allowing you to weigh infants accurately.

1. Weigh yourself first and let the measurement stabilize.
2. While still on the scale, pick up your baby and hold them.
3. The scale will automatically calculate and display your baby's weight by deducting your weight from the combined total.

Baby Mode

Let you hold your baby while weighing and deduct your weight from you and your baby's total weight.

The Baby Mode feature allows for easy and accurate weighing of infants.

CARE & MAINTENANCE

- Always use the scale on a solid and flat surface to ensure accuracy.
- The glass surface can be slippery when wet. Keep it dry before and after use.
- To clean the scale, wipe it with a damp cloth and mild detergent. Do not use abrasive cleaners or immerse the scale in water.
- DO NOT strike, shake, or drop the scale, as this may damage the sensors or glass platform.
- Store the scale in a cool, dry place away from direct sunlight and extreme temperatures.

Care and maintenance guidelines for the SENSUN Smart Body Fat Scale.


TROUBLESHOOTING

Error Code/Issue	Description	Solution
Err1	Overload. The weight exceeds the maximum capacity (150 kg / 330 lb).	Step off the scale immediately to avoid damage.
Err2	Low battery.	Replace the 2 x AAA batteries.
Err3	Calibration error.	Step off the scale, wait for it to turn off, then step on again to re-calibrate automatically. Ensure the scale is on a hard, flat surface.
Err4	Data transmission error.	Ensure Bluetooth is enabled on your phone and the MovingLife app is open. Try restarting the app or your phone.

Error Code/Issue	Description	Solution
No display/Scale not turning on	Batteries are dead or incorrectly installed.	Check battery polarity or replace batteries.
Inaccurate readings	Scale is on an uneven surface, or not properly calibrated.	Place the scale on a hard, flat surface. Step on and off to allow auto-calibration. Ensure feet are clean and dry.

SPECIFICATIONS

Model Number	IF1090D
Brand	SENSSUN
Maximum Capacity	150 kg / 330 lb
Measurement Increments	0.1 kg / 0.2 lb
Units	Kilograms (kg), Pounds (lb)
Power Source	2 x AAA batteries
Product Dimensions	28 x 24 x 2 cm (11 x 9.45 x 0.79 inches)
Item Weight	1.15 kg (1150 Grams)
Display Type	LCD with Backlit Display
Connectivity	Bluetooth 4.0
App Compatibility	MovingLife App (iOS 8.0+ / Android 4.3+)

Dimensions of the SENSSUN Smart Body Fat Scale.

WARRANTY & SUPPORT

SENSSUN provides professional support for its products. If you encounter any issues or have questions, please contact the support team.

- **Warranty:** The product comes with a 1-year warranty.
- **Support Response:** Expect a response from the SENSSUN support team within 24 hours.

Contact Information:

Manufacturer: ZHONGSHAN CAMRY ELECTRONIC CO., LTD.

Manufacturer Email: am-contact@camry.com.cn

EU Representative: E-CrossStu GmbH, Mainzer Landstr 89, Frankfurt am Main, Germany 80329

EU Representative Email: e-crosstu@outlook.com

UK Representative: YH Consulting Limited, C/O YH Consulting Limited Office 147, Centurion House, London Road, Staines-upon-Thames, Surrey, London. TW18 4AX

UK Representative Email: H2Yhuk@gmail.com

Model: IF1090D	Smart Body Fat Scale
Max.: 150kg/330lb	Battery: 2 x 1.5V AAA
d= 100g/0.2lb	S/N: E- IF 202406

Manufacturer: ZHONGSHAN CAMRY ELECTRONIC CO., LTD.

Address: Baishawan Industrial Park, Qiwan Road, East District,
ZHONGSHAN, GUANGDONG, CHINA, 528400
Email: am-contact@camry.com.cn

E-CrossStu GmbH
Mainzer Landstr.69, Frankfurt am Main, Germany 60329
Tel: +49 6933296767
Email: e-crossstu@outlook.com

YH Consulting Limited
C/O YH Consulting Limited Office 147, Centurion House,
London Road, Staines-upon-Thames, Surrey, London,
TW18 4AX
Tel: +44 07514-677868 Email: H2YHUK@gmail.com

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation

	WARNING KEEP OUT OF REACH OF CHILDREN Swallowing can lead to chemical burns, perforation of soft tissue, and death. Severe burns can occur within 2 hours of ingestion. Seek medical attention immediately.	
--	---	--

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Contact details for manufacturer and representatives, including FCC compliance information.