

JUMMICO MP

JUMMICO Power Lift Recliner Chair with Massage (Model MP) - User Manual

Brand: JUMMICO | Model: MP

1. INTRODUCTION

This manual provides essential information for the safe and efficient use of your JUMMICO Power Lift Recliner Chair with Massage, Model MP. Please read these instructions thoroughly before assembly and operation to ensure proper function and longevity of your product.

2. SAFETY INFORMATION

Important Safety Warnings:

- Ensure the user's weight does not exceed 300 pounds (136 kg). Exceeding this limit may cause damage to the chair and potential injury.
- Always follow instructions for the safe use of reclining and tilting features. Keep hands and feet clear of moving parts.
- Do not allow children or pets to play on or around the chair, especially when it is in motion.
- Ensure the power cord is not pinched or damaged. Do not operate the chair if the cord is compromised.
- Unplug the chair from the power outlet before cleaning or performing any maintenance.
- Avoid spilling liquids on the chair, especially near electrical components.

3. PACKAGE CONTENTS

Your JUMMICO Power Lift Recliner Chair is typically shipped in 1-2 packages. Please ensure all components are present before beginning assembly. If you receive only one package, please wait for the second package to arrive before proceeding.

- Chair Backrest
- Chair Seat Base with Armrests
- Power Adapter and Cables

- Remote Control for Lift and Recline
- Remote Control for Massage Function
- User Manual (this document)
- Any necessary assembly hardware (usually pre-attached or minimal)



Figure 3.1: The recliner chair is typically shipped in two separate packages.

4. SETUP AND ASSEMBLY

Assembly of the JUMMICO Power Lift Recliner Chair requires minimal steps and no additional tools. Follow these instructions carefully:

1. **Unpack Components:** Carefully remove all parts from the packaging. Place them on a soft, clean surface to prevent scratches.
2. **Attach Backrest:** Slide the backrest onto the metal brackets located on the seat base. Ensure it clicks securely into place.
3. **Connect Cables:** Locate the power cables extending from the backrest and seat base. Connect them according to the color-coded or labeled connectors. Connect the power adapter to the chair's main power input.
4. **Connect Remote Controls:** Plug in the remote controls for the lift/recline and massage functions into their respective ports, usually located on the side pocket or near the power input.
5. **Position Chair:** Place the chair in its desired location, ensuring there is adequate space behind and in front for reclining and lifting movements.
6. **Plug In:** Connect the power adapter to a wall outlet.



Figure 4.1: The fully assembled JUMMICO Power Lift Recliner Chair.

5. OPERATING INSTRUCTIONS

5.1 Power Lift and Recline Function

The chair features an electrically controlled lift and recline mechanism, operated by a dedicated remote control.

- **To Recline:** Press and hold the "DOWN" button on the lift remote control. The backrest will recline, and the footrest will extend. Release the button when the desired position is reached.
- **To Lift/Stand:** Press and hold the "UP" button on the lift remote control. The chair will gently tilt forward and lift, assisting you to a standing position. Release the button when the desired height is reached.
- **To Return to Upright Position:** From a reclined position, press and hold the "UP" button. The footrest will retract, and the backrest will return to an upright position.



Figure 5.1: Motorized control remote and chair positions.

POWER LIFT RECLINER CHAIR

Our lift chair helps the senior stand up easily

I	II	III	IV
Lifting sit's comfortable to lie flatwitch	Very comfortable to sit	It's easy to stand up	Give you different enjoyment

Figure 5.2: The power lift function assists users in standing up easily.

5.2 Massage Function

The chair includes a four-point massage system (back, waist, thighs, legs) with adjustable modes and intensity, controlled by a separate remote.

- **Power On/Off:** Press the "Power Button" on the massage remote to activate or deactivate the

massage function.

- **Select Massage Area:** Use the "Back Massage," "Waist Massage," "Thigh Massage," and "Leg Massage" buttons to select specific areas for massage. Multiple areas can be active simultaneously.
- **Adjust Intensity:** The "Three Vibration Intensity" button cycles through low, medium, and high intensity levels.
- **Change Mode:** The "Three Vibration Mode" button cycles through different massage patterns.
- **Set Timer:** Use the "Timer" button to set the massage duration to 15, 30, or 60 minutes. The massage will automatically stop after the selected time.



Figure 5.3: Massage function remote control and massage points.

6. MAINTENANCE

Proper care and maintenance will extend the life of your JUMMICO Power Lift Recliner Chair.

- **Cleaning:** For general cleaning, wipe the PU leather surface with a clean, dry cloth. For spills, immediately wipe with a clean, damp cloth and then dry thoroughly. Avoid harsh chemicals or abrasive cleaners.
- **Inspection:** Periodically inspect all electrical cords and connections for damage. Ensure all mechanical parts move smoothly and are free from obstructions.
- **Weight Limit:** Adhere strictly to the maximum weight recommendation of 300 pounds to prevent structural damage.



Figure 6.1: Side pocket for convenient storage.



Figure 6.2: Integrated cup holder in the armrest.

7. TROUBLESHOOTING

If you encounter issues with your JUMMICO Power Lift Recliner Chair, refer to the following common problems and solutions:

- **Chair Does Not Power On:**

- Ensure the power cord is securely plugged into both the chair and a functional wall outlet.
- Check if the power adapter is properly connected to the chair.
- Verify that the wall outlet is receiving power (test with another appliance).

- **Lift/Recline Function Not Working:**

- Confirm the power is connected and the chair is receiving electricity.
- Check the connection of the lift/recline remote control to the chair.
- Inspect for any obstructions preventing the chair's movement.
- If the lift function fails, and the chair is out of warranty, some users have successfully replaced the lift motor actuator. This typically involves searching for the part number (e.g., "385436504600" as mentioned in a customer review) and replacing it. However, DIY repair is not officially recommended and should only be attempted by individuals with appropriate technical skills and safety precautions.

- **Massage Function Not Working:**

- Ensure the massage remote control is properly connected to the chair.
- Verify the chair is powered on.
- Check if the massage unit's power cord is securely connected.

- **Unusual Noises During Operation:**

- Ensure no objects are caught in the chair's mechanism.
- If noises persist, contact customer support. Do not attempt to lubricate or disassemble internal components unless instructed by support.

If these troubleshooting steps do not resolve the issue, please contact JUMMICO customer service for further assistance.

8. SPECIFICATIONS

Feature	Detail
---------	--------

Brand	JUMMICO
Model Name	Power Lift Recliner Chair
Model Number	MP
Color	Black
Product Dimensions (D x W x H)	30"D x 31.8"W x 43"H
Item Weight	76.3 pounds (34.6 kg)
Maximum Weight Recommendation	300 Pounds (136 kg)
Material	PU Leather, High-Density Foam (Sponge), Metal Frame
Special Features	Power Lift, Reclining, 4-Point Massage (Back, Waist, Thighs, Legs), Cup Holders, Side Pockets
Product Care Instructions	Wipe Clean



Figure 8.1: Product dimensions in various positions.

9. WARRANTY AND SUPPORT

JUMMICO is committed to providing quality products and customer satisfaction.

- **Warranty Information:** Specific warranty details may vary. Please refer to the documentation included with your purchase or contact JUMMICO customer service for current warranty terms.
- **Customer Support:** For any questions, concerns, or technical assistance, please contact JUMMICO customer service. They provide professional support for all customers.
- **Contact Information:** Refer to your purchase receipt or the JUMMICO brand store on Amazon for the most up-to-date contact details.

© 2023 JUMMICO. All rights reserved.
This manual is subject to change without notice.