

Generic ADI 1025

Applikon ADI 1025 Bio Console Instruction Manual

Model: ADI 1025 | Brand: Generic

1. INTRODUCTION

This instruction manual provides essential information for the safe and efficient operation, setup, and maintenance of the Applikon ADI 1025 Bio Console. This device is designed for precise control and monitoring in various biological applications. Please read this manual thoroughly before operating the unit to ensure proper function and to prevent damage or injury.

2. SAFETY INFORMATION

Adherence to safety guidelines is paramount for the protection of personnel and equipment. Always observe the following precautions:

- Electrical Safety:** Ensure the unit is connected to a properly grounded power outlet with the correct voltage as specified on the device label. Do not operate with damaged power cords.
- Ventilation:** Position the console in an area with adequate ventilation to prevent overheating. Do not block ventilation openings.
- Liquid Spills:** Avoid liquid spills on or into the console. In case of a spill, immediately disconnect power and allow the unit to dry completely before re-connecting.
- Gas Connections:** Ensure all gas connections (Air, CO2) are secure and leak-free. Use appropriate pressure regulators for gas supply.
- Biological Hazards:** When working with biological samples, follow all standard laboratory safety procedures, including proper handling of biohazardous materials and waste disposal.
- Qualified Personnel:** Operation and maintenance should only be performed by trained and qualified personnel.

3. PRODUCT OVERVIEW AND COMPONENTS

The Applikon ADI 1025 Bio Console integrates various control and monitoring functions for bioreactor systems. Key components are detailed below:

Bio Console

29

ADI 1025

OVERLAY AIR


OVERLAY AIR

OVERLAY CO₂

POWER

STIRRER MOTOR

SENSOR CABLES

The control panel features three flowmeters with yellow liquid and black dials. The first two are labeled 'AIR' and the third 'CO₂'. To the right are several control buttons: 'OVERLAY AIR', 'AIR', 'CO₂', and 'OVERLAY'. Below these are two ports labeled 'STIRRER MOTOR' and 'SENSOR CABLES'. A green power button is at the top right, and a white power inlet with a yellow warning symbol is at the bottom right.

12875-ays-3
applied
RELIABLE INSTRUMENTS
made in Holland

Bio Console ADI 1025

PUMP 2

PUMP 3

PUMP 1

CONDENSER FROM

REACTOR FROM

TO

This panel contains three pump sections labeled 'PUMP 1', 'PUMP 2', and 'PUMP 3'. Each section has a central blue port and several smaller ports. At the bottom right, there are two sets of ports: 'CONDENSER FROM' and 'REACTOR FROM', with a 'TO' port below them.

Figure 3.1: Front view of the Applikon ADI 1025 Bio Console, showing the main control panel, gas flow meters, and pump connection points.




Figure 3.2: Detailed view of the front panel, highlighting the Air and CO₂ flow meters, Overlay Air controls, Power button, Stirrer Motor connector, and Sensor Cables input.




Figure 3.3: Close-up of the auxiliary power outlet located on the front panel, typically used for external devices.



Figure 3.4: Lower front panel connections, including ports labeled "CONDENSER FROM" and "REACTOR FROM/TO", indicating fluid or gas pathways.



GEORCIOT
ADJ: ZS10100115 Ver: 2003 CE
U: 220-240 V AC 50/60 Hz
P: 1600 VA T: SAVY 10A
I/O SIGNALS

GEORCIOT
ADJ: ZS10100115 Ver: 2003 CE
U: 220-240 V AC 50/60 Hz
P: 1600 VA T: SAVY 10A

Figure 3.5: Rear view of the console, displaying the main power input, outputs to bio controller and power unit, gas inlets (CO2, Air, Overlay Air), I/O signals port, and water inlet/outlet for flow condenser.



Figure 3.6: Detailed view of one of the rear product labels, showing "Applikon Dependable Instruments", model ADI Z510100110, year 2003, voltage 220-240 / 110-120 Vac 50/60 Hz, power 1600 VA, and fuse rating T 5A/T 10A.



Figure 3.7: Detailed view of a second rear product label, showing model ADI Z310250010, year 2003, voltage 220-240 / 110-120 Vac 50/60 Hz, power 1600 VA, and fuse rating T 8A / T 16A. Note the variation in model number and fuse rating compared to Figure 3.6, indicating potential product variations.



Figure 3.8: Close-up of the main green Power button and the multi-pin connector for the Stirrer Motor on the front panel.

Key Components:

- **Front Panel:** Features gas flow meters (Air, CO₂), overlay air controls, main power button, stirrer motor connector, sensor cable input, and auxiliary power outlet.
- **Pump Connections:** Designated areas for connecting external pumps (PUMP 1, PUMP 2, PUMP 3).
- **Rear Panel:** Includes main power input, outputs for bio controller and power unit, gas inlets, I/O signals port, and connections for flow condenser water.
- **Internal Components:** (Not visible externally) Control circuitry for precise regulation of parameters.

4. SETUP

Proper setup is crucial for the optimal performance and longevity of the ADI 1025 Bio Console.

1. Unpacking and Inspection:

Carefully remove the console from its packaging. Inspect the unit for any signs of shipping damage. Report any damage to the carrier and supplier immediately.

2. Placement:

Place the console on a stable, level surface in a laboratory environment. Ensure there is sufficient clearance around the unit for proper ventilation, especially at the rear panel where air vents are located. Avoid direct sunlight, excessive heat, or high humidity.

3. Power Connection:

Before connecting to power, verify that the voltage selector (if present, or check the rear label as shown in Figure 3.6 and 3.7) matches your local power supply (220-240 Vac or 110-120 Vac, 50/60 Hz). Connect the main power cord to the appropriate inlet on the rear panel (see Figure 3.5) and then to a grounded electrical outlet.

4. **Gas Connections:**

Connect regulated gas lines for Air and CO₂ to the corresponding inlets on the rear panel (see Figure 3.5). Ensure all connections are tight and leak-free. The front panel features flow meters for precise control of these gases (see Figure 3.2).

5. **External Device Connections:**

- **Stirrer Motor:** Connect the stirrer motor cable to the designated multi-pin connector on the front panel (see Figure 3.8).
- **Sensor Cables:** Connect relevant sensor cables (e.g., pH, DO, temperature) to the "SENSOR CABLES" input on the front panel.
- **Pumps:** Connect external pumps to the designated pump control ports (PUMP 1, PUMP 2, PUMP 3) on the front panel.
- **Condenser/Reactor:** Connect tubing for the condenser and reactor to the ports on the lower front panel (see Figure 3.4) and the water inlet/outlet for the flow condenser on the rear panel (see Figure 3.5).
- **I/O Signals:** If integrating with external control systems, connect the I/O signals cable to the D-sub connector on the rear panel (see Figure 3.5).

5. OPERATING INSTRUCTIONS

This section outlines the basic steps for operating the Applikon ADI 1025 Bio Console.

1. **Powering On:**

Once all connections are secure, press the green "POWER" button on the front panel (see Figure 3.8) to turn on the console. The indicator light will illuminate.

2. **Gas Flow Control:**

Adjust the flow rates for Air and CO₂ using the respective knobs above the flow meters on the front panel (see Figure 3.2). Monitor the flow rates directly on the rotameters. The "OVERLAY AIR" control allows for separate air flow to the headspace of the bioreactor.

3. **Stirrer Motor Operation:**

Ensure the stirrer motor is properly connected. The console provides control for the stirrer motor, typically through an external bio controller connected via the "TO BIO CONTROLLER" port on the rear (see Figure 3.5).

4. **Sensor Monitoring:**

If sensors are connected, their signals will be processed by the console and typically relayed to an external bio controller for display and data logging.

5. **Pump Control:**

The console provides interfaces for controlling external pumps (PUMP 1, PUMP 2, PUMP 3). Refer to the specific pump manuals for their operation, as the console acts as a control interface.

6. **Powering Off:**

To turn off the console, press the green "POWER" button again. Disconnect the main power cord from the wall outlet if the unit will not be used for an extended period.

6. MAINTENANCE

Regular maintenance ensures the longevity and reliable performance of your Applikon ADI 1025 Bio Console.

- **Cleaning:**

Wipe the exterior of the console with a soft, damp cloth. For stubborn stains, a mild detergent can be used. **Do not use abrasive cleaners or solvents.** Ensure no liquids enter the ventilation openings or connectors.

- **Connection Checks:**

Periodically inspect all electrical, gas, and fluid connections for wear, damage, or leaks. Replace any worn or damaged tubing and cables immediately.

- **Ventilation:**

Ensure that the ventilation grilles on the rear panel remain clear of dust and debris to prevent overheating. Use a soft brush or compressed air to clean them if necessary.

- **Calibration:**

While the console itself primarily acts as a control unit, ensure that any connected sensors (pH, DO, etc.) and pumps are calibrated according to their respective manufacturer's instructions to maintain accuracy.

- **Storage:**

If storing the unit for an extended period, disconnect all cables, clean the unit, and store it in a dry, dust-free environment within the specified operating temperature range.

7. TROUBLESHOOTING

This section provides solutions to common issues you might encounter with the Applikon ADI 1025 Bio Console. For problems not listed here, contact technical support.

Problem	Possible Cause	Solution
Unit does not power on.	<ul style="list-style-type: none">• Power cord not connected.• No power from outlet.• Internal fuse blown.	<ul style="list-style-type: none">• Ensure power cord is securely connected to the console and a live outlet.• Check the circuit breaker or try a different outlet.• Contact qualified service personnel for fuse replacement. Do not attempt to replace fuses yourself unless explicitly instructed and trained.
Gas flow meters show no flow.	<ul style="list-style-type: none">• Gas supply off or empty.• Gas lines kinked or disconnected.• Flow control knob closed.	<ul style="list-style-type: none">• Verify gas cylinder is open and has sufficient pressure.• Inspect gas lines for obstructions or disconnections.• Rotate the flow control knob counter-clockwise to open the valve.
Stirrer motor not operating.	<ul style="list-style-type: none">• Motor not connected.• Issue with external bio controller.	<ul style="list-style-type: none">• Ensure stirrer motor cable is securely connected to the front panel port.• Check the settings and status of the external bio controller.

8. SPECIFICATIONS

The following specifications apply to the Applikon ADI 1025 Bio Console:

Parameter	Value
-----------	-------

Parameter	Value
Model Number	ADI 1025
Manufacturer Part Number (MPN)	1025
Brand	Generic (Applikon Dependable Instruments)
ASIN	B0BZXLP362
Color	Black (console housing is light grey/beige)
Date First Available	February 15, 2019
Input Voltage	220-240 Vac / 110-120 Vac (selectable/configurable)
Frequency	50/60 Hz
Power Consumption	1600 VA
Fuse Rating (Example 1)	T 5A / T 10A (refer to specific unit label)
Fuse Rating (Example 2)	T 8A / T 16A (refer to specific unit label)
Year of Manufacture (Example)	2003 (refer to specific unit label)

Note: Specific electrical ratings and model variations may differ. Always refer to the product label on your specific unit for precise specifications.

9. WARRANTY AND SUPPORT

Information regarding the warranty period and specific support contacts for the Applikon ADI 1025 Bio Console is not available in this document. For warranty claims, technical support, or service inquiries, please contact the original vendor or manufacturer directly. Ensure you have your product's serial number and purchase details ready when contacting support.

No official product videos from the seller were provided in the available data for inclusion in this manual.




© 2024 Generic. All rights reserved. This manual is for informational purposes only.



Bio Bundles 1 - 15 Liter (US Version)

for Microbial Applications



HARDWARE & INSTALLATION MANUAL

Apr 9 2004



[\[pdf\]](#) User Manual Instructions

G MANUALSB Manuals CE Bundles USA Microbial Hardware Inst Hardware1 7 Installation 110 wpd JVe for 39 Applications manual contains information and warnings which have to be followed by the user ensure safe operation retain equipment in condition Before switching on make sure that it is set line voltage This designed bi o process control must not used other purposes! Caution BioBundle 1 15L InstallationManual archive resources coleparmer Manual s

Bio Bundles 1 -15 Liter US Version for Microbial Applications rnD- --- - HARDWARE INSTAIAT ... , The Netherlands, hereby certifies that the following set of instruments: ADI 1010 BIO CONTROLLER, **ADI 1025 BIO CONSOLE** and ADI 1032 STIRRER CONTROLLER meets the requirements of the EU Directives 89/...

lang:en score:29 filesize: 6.34 M page_count: 69 document date: 2004-06-07



[\[pdf\]](#) User Manual

15LiterBioBundleSystem Cole Parmer Laboratory Pilot Plant Projects Applikon Biotechnology Ordering Information 15 Liter BioBundle System Z4BUND0015C cell culture Z4BUND0015M microbial leaflet15L archive resources coleparmer Manual s

15 Liter BioBundle System Complete fermentation system for microbial and cell culture applications ... y years of successful usage. Applikon Biotechnology ADI 1010 Bio Controller Applikon Biotechnology **ADI 1025 Bio Console** High precision motor pH, dissolved O2 and temperature sensors 15 Liter BioBu...

lang:en score:26 filesize: 2.63 M page_count: 6 document date: 2004-11-16



[\[pdf\]](#) User Manual

7LiterBioBundleSystem Cole Parmer Each BioBundle system component is thoroughly tested and evaluated based on many years of successful usage in a wide range applications Below partial listing the components leaflet7L archive resources coleparmer Manual s

7 Liter BioBundle System Complete fermentation system for microbial and cell culture applications ... y years of successful usage. Applikon Biotechnology ADI 1010 Bio Controller Applikon Biotechnology **ADI 1025 Bio Console** High precision motor pH, dissolved O2 and temperature sensors Easy access hea...

lang:en score:26 filesize: 2.57 M page_count: 6 document date: 2004-11-16




[pdf]

COL01181 SUB Verco CC Micocarriers App Note 03 24 17 indd sub vero cell culture on microcarriers app note assets thermofisher TFS Assets BPD Application Notes

APPLICATION NOTE HyPerforma Single-Use Bioreactor HyPerforma Single-Use Bioreactor for Vero cell c ... It User Guide UG003 and equipped with a pH/DO control system Applikon ADI 1010 Bio Controller and **ADI 1025** Bio Console . Direct air flow to the sparger was set to 0.1 L/min and the overlay gas flow ...

lang:en score:13 filesize: 1.26 M page_count: 4 document date: 2017-03-24




[pdf]

Applicatione Note SUB Verco CC Micocarriers COL01181 Feb 2016 null cq textimage jsp b1839a52 643c 491f a6e5 8ae833dc5aa5MainParsystextimage 8a5a max width 100% img 230px Form Success Thermo Fisher Scientific US thermofisher content dam LifeTech Documents |||

APPLICATION NOTE HyPerforma Single-Use Bioreactor HyPerforma Single-Use Bioreactor for Vero cell c ... It User Guide UG003 and equipped with a pH/DO control system Applikon ADI 1010 Bio Controller and **ADI 1025** Bio Console . Direct air flow to the sparger was set to 0.1 L/min and the overlay gas flow ...

lang:en score:13 filesize: 1.27 M page_count: 4 document date: 2017-03-24



[pdf] Catalog

Semiconductors Integrated Circuits Temperature Controllers Sensors and Transducers con 07TL 898 common leocom jp cat catalogpage

Farnell Page 898 Date: 06-09-06 time:19:47 898 Semiconductors - Integrated Circuits Temperature Co ... or Transmitter Cerdip AD693AQ 694 ADI 694 ADI 694 ADI 1021 **ADI 1025** LT 16 DIL 16 DIL 16 SOIC 16 QSOP 8 DIL 4mA to 20mA Transmitter

lang:en score:12 filesize: 54.04 K page_count: 1 document date: 0000-00-00