



Manuals.plus /

› Aidisy /

› Aidisy Portable Nebulizer (Green) User Manual

## Aidisy Aidisy Portable Nebulizer Green

# Aidisy Portable Nebulizer (Green) User Manual

Model: Aidisy Portable Nebulizer Green

## 1. INTRODUCTION

Thank you for choosing the Aidisy Portable Nebulizer. This device is designed to provide efficient and convenient respiratory care for both adults and children. It utilizes advanced ultrasonic nebulization technology to convert liquid medication into fine mist particles for effective inhalation. Please read this manual thoroughly before use to ensure proper operation and maintenance.

## 2. SAFETY INFORMATION

- Do not use oil-based medications or solutions not intended for nebulization.
- Use only purified water for cleaning the device.
- Keep the device out of reach of children without adult supervision.
- Consult a healthcare professional for appropriate medication dosage and suitability.
- Ensure the device is fully charged before each use for optimal performance.
- Do not immerse the main unit of the nebulizer in water or any other liquid.
- If the device malfunctions or shows signs of damage, discontinue use immediately and contact customer support.

## 3. PACKAGE CONTENTS

The Aidisy Portable Nebulizer package includes:

- 1 x Aidisy Portable Nebulizer Main Unit
- 1 x Adult Mask
- 1 x Child Mask
- 1 x Mouthpiece

- 1 x USB-C Charging Cable
- 1 x User Manual (Multi-language)

## 4. PRODUCT OVERVIEW

Familiarize yourself with the components of your Aidisy Portable Nebulizer:



Image: The Aidisy Portable Nebulizer showing its main unit, medication cup, power button, and various attachments including an adult mask and mouthpiece. A person is shown using the nebulizer with a mask.

- **Main Unit:** Contains the power button, indicator lights, and USB-C charging port.
- **Medication Cup:** Removable cup for holding liquid medication or purified water.
- **Nebulizer Head:** Where the mist is generated.
- **Power Button:** Used to turn the device on/off and select nebulization modes.
- **Indicator Lights:** Display charging status and operating mode.
- **Masks (Adult & Child):** For inhalation through the nose and mouth.

- **Mouthpiece:** For direct oral inhalation.

## 5. SETUP

### 5.1 Charging the Device

1. Connect the provided USB-C charging cable to the nebulizer's charging port.
2. Plug the other end of the cable into a compatible USB power source (e.g., laptop, USB wall adapter, power bank).
3. The indicator light will turn **yellow** while charging.
4. Once fully charged, the indicator light will turn **blue**. A full charge provides approximately 120 minutes of continuous use.



Image: Diagram illustrating the charging process and indicator light colors: yellow for charging, blue for fully charged and ready for use, and alternating yellow/blue for automatic cleaning mode.

### 5.2 Filling the Medication Cup

1. Open the lid of the medication cup.
2. Carefully pour the prescribed liquid medication into the cup. Do not overfill beyond the maximum fill line.
3. Close the lid securely.

# Easy To Operate



①

Open the lid of the medicine cup



②

Inject liquid medicine or purified water



③

Install the mouthpiece



④

Short press the power button

Image: A four-step visual guide demonstrating how to operate the nebulizer: 1. Open the lid of the medicine cup. 2. Inject liquid medicine or purified water. 3. Install the mouthpiece. 4. Short press the power button.

## 5.3 Attaching Accessories

1. Select the appropriate accessory: adult mask, child mask, or mouthpiece.
2. Securely attach the chosen accessory to the nebulizer head. Ensure it is firmly connected to prevent leaks.

# Multiple Masks

Adapt to different atomization needs



Image: The Aidisy Portable Nebulizer shown with three different attachments: an adult mask, a mouthpiece, and a child's mask, indicating adaptability for various users.

## 6. OPERATING THE NEBULIZER

### 6.1 Powering On and Selecting Mode

1. With the medication cup filled and accessory attached, short press the power button once to turn on the device.
2. The device will start in the **Strong Mode** (nebulization speed: 0.25 ml/min), suitable for adults.
3. Press the power button a second time to switch to **Mild Mode** (nebulization speed: 0.15 ml/min), suitable for children.
4. Press the power button a third time to switch to **Weak Mode** (nebulization speed not specified, but for gentler use).

### 6.2 Inhalation

- Hold the device with the mask firmly against your face, covering your nose and mouth, or place the

mouthpiece in your mouth.

- Breathe normally and deeply to inhale the mist.
- The device operates quietly, with a running noise level of less than 25dB, allowing for discreet use.

# Keep Quiet

Running Noise  $\leq 25\text{dB}$



Image: The Aidisy Portable Nebulizer operating quietly (less than 25dB) on a bedside table next to a sleeping child, highlighting its low noise level.

## 6.3 Automatic Shutdown

For safety and to prevent dry burning, the nebulizer is equipped with an intelligent chip that automatically stops operation after 20 minutes of continuous use. If further nebulization is required, simply restart the device.

## 7. MAINTENANCE

### 7.1 Self-Cleaning Mode

To prevent clogging and ensure optimal performance, use the self-cleaning mode regularly:

1. After each use, empty any remaining medication from the cup.

2. Add a small amount of purified water to the medication cup.
3. Long press the power button for 5 seconds. All three indicator lights will illuminate, and then alternate **yellow** and **blue**, indicating the self-cleaning mode is active.
4. Allow the device to run through the cleaning cycle.

## Cleaning Mode

Long press for 5 seconds , Eliminate drug residues



Image: The Aidisy Portable Nebulizer demonstrating its cleaning mode, activated by a 5-second long press of the power button, designed to eliminate drug residues. The image also highlights the 20-minute automatic shutdown feature.

## 7.2 Cleaning the Medication Cup and Accessories

1. Remove the medication cup from the main unit.
2. Rinse the medication cup thoroughly with purified water to remove any medication residue.
3. Wash the masks and mouthpiece with warm water and a mild, non-abrasive soap.
4. Rinse all accessories thoroughly to remove soap residue.
5. Allow all components to air dry completely before reassembling or storing.

## 7.3 Storage

Store the clean and dry nebulizer and its accessories in a cool, dry place, away from direct sunlight and extreme temperatures.

## 8. TROUBLESHOOTING

Problem	Possible Cause	Solution
No mist or weak mist produced	Low battery Medication cup improperly filled or attached Nebulizer head or medication cup clogged Medication not suitable for nebulization	Charge the device fully. Ensure medication cup is correctly filled and securely attached. Perform self-cleaning or manually clean components. Consult your healthcare professional regarding medication compatibility.
Device not turning on	Low battery Power button not pressed correctly	Charge the device fully. Press the power button firmly and briefly.
Liquid leaking from the device	Medication cup not securely closed or attached Medication cup overfilled	Ensure the medication cup is properly closed and firmly attached to the main unit. Do not exceed the maximum fill line in the medication cup.

## 9. SPECIFICATIONS

Feature	Detail
Brand	Aidisy
Model	Portable Nebulizer Green
ASIN	B0C8146L46
Dimensions (L x W x H)	6 x 4 x 11 cm
Weight	250 grams
Battery Type	2 Lithium-ion (included)
Charging Interface	USB-C
Nebulization Rate (Strong Mode)	0.25 ml/min
Nebulization Rate (Mild Mode)	0.15 ml/min
Particle Size	< 5 $\mu$ m
Operating Noise Level	$\leq$ 25 dB
Continuous Use Time (Fully Charged)	Up to 120 minutes
Automatic Shut-off	20 minutes

## 10. WARRANTY AND SUPPORT

The Aidisy Portable Nebulizer comes with a standard manufacturer's warranty. For specific warranty details, please refer to the documentation included with your purchase or contact Aidisy customer support.

If you encounter any issues or have questions regarding your device, please contact our customer support team for assistance. Contact information can typically be found on the product packaging or the official Aidisy website.