

Rayovac 13AE-60

Rayovac Extra Advanced Size 13 Hearing Aid Battery User Manual

Model: 13AE-60

INTRODUCTION

This manual provides essential information for the safe and effective use of Rayovac Extra Advanced Size 13 Hearing Aid Batteries. These batteries are designed to provide reliable power for your hearing devices. Rayovac Extra Advanced Size 13 batteries are Zinc Air, mercury-free, and feature an orange tab for easy identification. Each pack contains 60 batteries (10 cards of 6 batteries each).

PRODUCT FEATURES

- Specifically designed for hearing aids.
- Long-lasting power with a 2-3 year shelf life.
- Mercury-free Zinc Air chemistry.
- Identifiable by the orange tab.
- Manufactured in the USA.

SETUP AND INSTALLATION

Follow these steps to properly install your Rayovac Extra Advanced Size 13 hearing aid batteries:

1. **Prepare the Battery:** Remove a single battery from its packaging. Locate the orange protective tab on the battery.
2. **Activate the Battery:** Peel off the orange protective tab. This exposes the battery to air, activating it. For optimal performance, allow the battery to sit for approximately 1-2 minutes after removing the tab before inserting it into your hearing aid. This allows the battery to fully activate.
3. **Open Hearing Aid Battery Compartment:** Gently open the battery compartment of your hearing aid. Refer to your hearing aid's specific manual for instructions on opening its battery door.
4. **Insert Battery:** Insert the activated battery into the compartment with the correct polarity. Most hearing aids have a diagram indicating the correct orientation (+/-). Ensure the flat side of the battery

faces the correct direction as indicated by your hearing aid manufacturer.

5. **Close Compartment:** Carefully close the battery compartment door. Do not force it. If the door does not close easily, recheck the battery's orientation.
6. **Test Hearing Aid:** Turn on your hearing aid and test its functionality.

NEW!
NOUVEAU!

RAYOVAC[®]

EXTRA[®] **ADVANCED**

LONG LASTING GUARANTEED!
LA PLUS DURABLE GARANTIE!



6 PK

MADE IN THE
USA

Image: Rayovac Extra Advanced Size 13 Hearing Aid Battery packaging, showing six batteries in a circular dispenser. The orange tab is visible on each battery.

OPERATING INSTRUCTIONS

Once installed, Rayovac Extra Advanced Size 13 batteries operate continuously until their power is depleted. The lifespan of the battery will vary depending on the hearing aid model, usage patterns, and environmental conditions.

Battery Life: The typical battery life for Size 13 batteries can range from 3 to 10 days, but this is an estimate. Factors such as streaming audio, using wireless features, and the power requirements of your specific hearing aid will affect battery longevity.

Low Battery Indication: Your hearing aid will typically provide an audible or visual alert when the battery is low. Replace the battery promptly when this indication occurs to ensure uninterrupted performance.

MAINTENANCE AND STORAGE

Proper maintenance and storage can extend the shelf life of your batteries and ensure optimal performance:

- **Storage:** Store unused batteries in their original packaging at room temperature (approximately 20-25°C or 68-77°F) in a dry place. Avoid extreme temperatures, direct sunlight, and high humidity. Do not store batteries in a refrigerator.
- **Handling:** Always handle batteries with clean, dry hands. Oils and dirt can affect battery performance. Avoid touching the activated side of the battery.
- **Disposal:** Dispose of used batteries responsibly according to local regulations. Do not dispose of batteries in household waste. Many communities have recycling programs for batteries.
- **Keep Out of Reach of Children:** Batteries can be a choking hazard. Store them securely away from children and pets. If swallowed, seek immediate medical attention.



Image: Two packs of Rayovac Extra Advanced Size 13 Hearing Aid Batteries with the outer box, displaying product information and barcode.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Hearing aid not working after battery replacement.	Battery not fully activated, incorrect polarity, or faulty battery.	Ensure the orange tab was removed and the battery sat for 1-2 minutes. Check battery orientation. Try a new battery.
Short battery life.	High hearing aid usage, environmental factors, or battery nearing end of shelf life.	Reduce streaming or wireless features if possible. Ensure proper storage. Check the expiration date on the battery packaging.
Difficulty opening battery packaging.	Secure packaging design.	Use scissors to carefully cut open the packaging if needed.

SPECIFICATIONS

Attribute	Detail
Model Number	13AE-60
Battery Size	13 (Orange Tab)
Chemistry	Zinc Air (Mercury Free)
Voltage	1.45V
Quantity	60 batteries (10 x 6-count cards)
Equivalent Models	13A, A13, 13AE, DA13, 13HPX, AC13 P13, PR13H, ZA13, PR48, B0134, B26PA, DA13H, DA13N, 13HP, AC13E, AC13EZ, ME8Z, L13ZA, W13ZA, S13A, 13SA
Shelf Life	2-3 years
Manufacturing Origin	Made in USA
Item Weight	4 ounces
Package Dimensions	4.37 x 2.52 x 2.05 inches

WARRANTY AND SUPPORT

Specific warranty information for Rayovac Extra Advanced Size 13 Hearing Aid Batteries is typically provided on the product packaging or through the point of purchase. For detailed warranty terms or customer support, please refer to the official Rayovac website or contact their customer service directly.

Rayovac Official Website: www.rayovac.com