



[Manuals.plus](#) /

> [FEELLIFE](#) /

> FEELLIFE AEROGO-2025 Portable Mesh Nebulizer User Manual

FEELLIFE AEROGO-2025

FEELLIFE AEROGO-2025 Portable Mesh Nebulizer User Manual

Model: AEROGO-2025 | Brand: FEELLIFE

1. PRODUCT OVERVIEW

The FEELLIFE AEROGO-2025 is a state-of-the-art portable mesh nebulizer designed for efficient and quiet medication delivery. Utilizing advanced AIMESH technology and a German AI chip, this device atomizes liquid medication into a fine mist, ensuring targeted treatment for respiratory conditions in both children and adults. Its compact and lightweight design, coupled with convenient USB charging, makes it an ideal solution for home use and travel, providing relief anytime, anywhere.

2. PRODUCT COMPONENTS

The AEROGO-2025 nebulizer comes with all necessary components for immediate use and versatile application.



Image: The FEELLIFE AEROGO-2025 nebulizer shown with its included accessories: a mouthpiece, an adult mask, and a child mask, along with the main device unit.

- **Main Unit:** The core device containing the power button, battery, and nebulization mechanism.
- **Medication Cup:** Detachable cup for holding liquid medication, featuring the mesh atomizer.
- **Mouthpiece:** For direct inhalation of medication.
- **Adult Mask:** A larger mask designed for adult users.
- **Child Mask:** A smaller mask suitable for children and infants.
- **USB Charging Cable:** For recharging the device's internal battery.
- **Storage Bag:** A convenient bag for carrying and protecting the nebulizer and its accessories.

3. KEY FEATURES

3.1 AIMESH Technology

The AEROGO-2025 utilizes advanced AIMESH technology, featuring a microporous membrane made from polymer structure

materials and integrated with German AI chip technology. This innovative design ensures efficient atomization of liquid medication through polymer membrane vibrations. It achieves a high proportion of respirable particles (MMAD: 3 μm) and an aerosolization rate of ≥ 0.4 ml/min, guaranteeing targeted treatment with minimal residual medication (less than 0.1 ml).



AIMESH-TECHNOLOGIE
MMAD $\leq 3 \mu\text{m}$

Obere Atemwege **Luftröhre** **Bronchien & Alveolen**

Direkt zum betroffenen Bereich

Image: Visual representation of how AIMESH technology delivers medication directly to the upper airways, trachea, bronchi, and alveoli.



 $\leq 25\text{dB}$

FLÜSTERLEISER BETRIEB

Verabschieden Sie sich vom Lärm.
Geräuschpegel <25 dB

Image: The nebulizer in operation, demonstrating its quiet performance, allowing for use without disturbing sleep.

3.2 Quiet Nebulization

The FEELLIFE nebulizer operates with a noise level below 25 dB, making it nearly silent. This quiet operation, facilitated by the German AI chip, ensures a comfortable and convenient inhalation experience. It can be used at home, in the office, or in a classroom without disturbing others, making it suitable for use even during sleep.

3.3 Lightweight and Portable Design

The AEROGO-2025 is designed to be exceptionally small and elegant, allowing it to be easily stored in a school bag or clothing pocket. It is USB-rechargeable, providing up to 3.5 hours of continuous use on a full charge. Its portable design, recognized by a German Red Dot Award, ensures you can perform inhalations anytime, anywhere.



Image: The compact nebulizer being placed into a bag, highlighting its portability for on-the-go use.

UNIVERSELL ZUGÄNGLICHES DESIGN

Für unterschiedliche Zerstäubungsbedürfnisse.



Image: Demonstrates the versatility of the nebulizer with a mouthpiece, adult mask, and child mask, suitable for all family members.

3.4 Versatile Home Use

The FEELLIFE nebulizer includes masks and mouthpieces for both children and adults, making it suitable for the entire family while ensuring hygienic use. The medication cup is easily removable for cleaning, making it perfect for daily travel and home care.

Note: This device is not suitable for essential oils or viscous liquids.

4. SETUP INSTRUCTIONS

Follow these steps to set up your FEELLIFE AEROGO-2025 nebulizer for first use or after cleaning.

- 1. Charge the Device:** Before first use, fully charge the nebulizer using the provided USB cable. Connect the cable to the device and a compatible USB power source (e.g., laptop, power bank, AC adapter). The green light will blink during charging and turn solid green when fully charged. A full charge takes approximately 1.5 hours and allows for about 6 uses.
- 2. Open the Lid:** Gently open the lid of the medication cup.
- 3. Fill with Liquid:** Carefully pour the prescribed liquid medication into the medication cup. Ensure the liquid level does not exceed the maximum fill line. *Do not use essential oils or viscous liquids.*

4. **Select and Attach Accessory:** Choose the appropriate accessory (mouthpiece, adult mask, or child mask) for the user. Attach it securely to the nebulizer's outlet.
5. **Prepare for Use:** Ensure all connections are secure and the device is held upright.



Image: A visual guide illustrating the five easy steps for setting up and operating the nebulizer: charging, opening lid, filling with liquid, selecting mask, and pressing the button.

5. OPERATING INSTRUCTIONS

Operating the AEROGO-2025 nebulizer is simple and intuitive.

1. **Power On:** Press the power button located on the front of the device. The blue light will illuminate, indicating the device is working and nebulization has started.
2. **Inhale Medication:** Place the mouthpiece in your mouth or position the mask over your nose and mouth. Breathe normally and deeply to inhale the fine mist.
3. **Monitor Operation:** The device will continue to nebulize until the medication cup is empty or the battery is low.
4. **Power Off:** Once nebulization is complete or if you need to stop, press the power button again to turn off the device. The blue light will turn off.
5. **Battery Indicator:**
 - **Blue Light Blinks 5 Times:** Indicates low battery, device will shut down soon.
 - **Green Light Blinks:** Charging in progress.
 - **Blue Light Solid:** Device is operating.
 - **Green Light Solid:** Fully charged.

EIN-KNOPF-START

**Vollständige Aufladung 1,5 Stunden,
6-mal verwendbar.**

 **Blaues Licht blinkt 5-mal**
Batterie schwach, Abschaltung

 **Grünes Licht blinkt**
Ladevorgang läuft

 **Blaues Licht leuchtet**
Bei der Arbeit

 **Grünes Licht leuchtet**
Vollständig geladen



Image: Illustration of the single power button and the different LED light indicators for charging status, low battery, and active operation.

6. MAINTENANCE AND CLEANING

Proper maintenance and cleaning are essential for the longevity and hygienic operation of your nebulizer.

1. After Each Use:

- Remove the medication cup from the main unit.
- Discard any remaining medication.
- Rinse the medication cup, mouthpiece, and mask with clean, warm water.
- Shake off excess water and allow all parts to air dry completely on a clean surface.

2. Daily Disinfection (if used frequently):

- After rinsing, soak the medication cup, mouthpiece, and mask in a solution of 70% isopropyl alcohol or a mild disinfectant solution recommended for medical devices for 10-15 minutes.
- Alternatively, a solution of white vinegar and distilled water (1 part vinegar to 3 parts distilled water) can be used for

soaking.

- Rinse thoroughly with sterile water or distilled water to remove any disinfectant residue.
- Allow all parts to air dry completely.

3. Cleaning the Main Unit:

- Wipe the exterior of the main unit with a soft, damp cloth.
- *Do not submerge the main unit in water or any liquid.*

4. Storage:

- Store the clean and dry nebulizer and its accessories in the provided storage bag in a cool, dry place.

Important Notes:

- Do not use abrasive cleaners or harsh chemicals.
- Ensure all parts are completely dry before reassembling or storing to prevent bacterial growth.
- The device is not suitable for essential oils or viscous liquids, as these can damage the mesh and affect performance.

7. TROUBLESHOOTING

If you encounter issues with your FEELLIFE AEROGO-2025 nebulizer, refer to the table below for common problems and their solutions.

Problem	Possible Cause	Solution
Device does not turn on.	Battery is depleted.	Charge the device fully using the USB cable. Ensure the green light blinks during charging and turns solid green when complete.
No mist or very weak mist.	Insufficient medication in the cup. Medication cup or mesh is clogged. Incorrect liquid used (e.g., essential oils, viscous liquids). Medication cup not properly attached.	Add more medication, ensuring it is above the minimum level. Clean the medication cup and mesh thoroughly as per maintenance instructions. Only use prescribed liquid medications suitable for nebulization. Ensure the medication cup is securely attached to the main unit.
Blue light blinks 5 times and device shuts off.	Low battery.	Recharge the device.
Medication cup is difficult to remove.	Secure magnetic connection.	Apply firm, steady pressure when pulling the medication cup away from the main unit. It is designed to have a secure fit.

If the problem persists after trying these solutions, please contact FEELLIFE customer support.

8. SPECIFICATIONS

Feature	Detail
Model	AEROGO-2025

Feature	Detail
Brand	FEELLIFE
Technology	AIMESH (Microporous Membrane with German AI Chip)
Particle Size (MMAD)	3 µm
Aerosolization Rate	≥0.4 ml/min
Residual Medication	<0.1 ml
Noise Level	<25 dB
Power Source	USB Rechargeable
Full Charge Time	Approx. 1.5 hours
Usage per Full Charge	Up to 3.5 hours (or 6 uses)
Product Dimensions	5.5 x 4 x 9.2 cm
Product Weight	380 grams
Country of Origin	China

9. WARRANTY AND SUPPORT

FEELLIFE products, including the AEROGO-2025 nebulizer, are manufactured to high-quality standards and come with a standard warranty. For specific warranty details, please refer to the warranty card included in your product packaging or contact FEELLIFE customer service.

For any questions, technical support, or warranty claims, please contact FeellifeDirect, the official seller, or visit the FEELLIFE brand store on Amazon:

Seller: FeellifeDirect

FEELLIFE Amazon Store: [Visit Store](#)