

Manuals.plus /

- › Phanpy /
- › Phanpy Wearable Breast Pump Milk Collector Cup S2 Instruction Manual

Phanpy AM9951

Phanpy Wearable Breast Pump Milk Collector Cup S2 Instruction Manual

Model: AM9951

INTRODUCTION

Thank you for choosing the Phanpy Wearable Breast Pump Milk Collector Cup S2. This manual provides essential information for the safe and effective use, assembly, cleaning, and maintenance of your milk collector cups. Please read these instructions thoroughly before first use and retain them for future reference.



Image: Two Phanpy S2 Milk Collector Cups, demonstrating their capacity and design.

PRODUCT OVERVIEW AND COMPONENTS

The Phanpy Wearable Breast Pump Milk Collector Cup S2 is designed for comfortable and efficient milk collection. Each cup weighs approximately 0.25 lb and features a BPA-free, food-grade silicone flange for skin-friendly contact. The design incorporates a 97° 3D simulated breastfeeding angle to aid milk production. The cups are self-standing for convenience and include a lanyard and carrying rope for portability.

The S2 cups are compatible with Phanpy E-Rosy breast pumps, Phanpy E-Shine breast pumps, and may be suitable for other breast pump models. Please verify compatibility before use.

Included Components:



Image: Detailed breakdown of the Phanpy S2 Milk Collector Cup's individual parts.

1. Breast Shield
2. Newcup (Milk Collection Cup)
3. Duckbill Valve
4. 20mm Breast Pump Flange Insert (24mm flange is integrated)
5. Top Cover
6. Tee Joint Connector
7. Silicone Airbag
8. Tubing (not always included with cups, but shown for connection)

360° Wrap your Breast

Comfortable and Free



Two sizes



20mm



24mm



Food-grade silicone

Image: The 24mm flange and optional 20mm insert, highlighting the food-grade silicone material.

SETUP AND INSTALLATION

Proper assembly ensures optimal performance and prevents leaks.

1. **Prepare Components:** Ensure all components are clean and dry before assembly.
2. **Insert Duckbill Valve:** Firmly push the duckbill valve into the bottom opening of the Newcup. Ensure it is securely seated.
3. **Attach Breast Shield:** Align the breast shield with the Newcup and press firmly to create a tight seal.
4. **Insert Flange Insert (Optional):** If using the 20mm insert, place it inside the 24mm breast shield, ensuring it sits flush.
5. **Connect Top Cover:** Place the silicone airbag onto the Tee joint connector, then attach the Tee joint connector to the top opening of the Newcup. Secure the top cover over the Tee joint connector.
6. **Connect Tubing:** Attach one end of the breast pump tubing to the Tee joint connector on the top cover and the other end to your compatible breast pump unit.

Easier Installation



Image: Visual guide for assembling the milk collector cup components.

OPERATING INSTRUCTIONS

Follow these steps for effective milk collection:

1. **Position the Cup:** Place the assembled milk collector cup against your breast, ensuring your nipple is centered within the flange opening. The 97° 3D simulated breastfeeding design aims for a comfortable fit.
2. **Secure in Bra:** Position the cup securely within your nursing bra. The compact and lightweight design allows for discreet use.
3. **Start Pumping:** Turn on your compatible breast pump. Adjust the suction and cycle settings according to your comfort and milk flow.
4. **Monitor Milk Collection:** The cup has clear measurement markings (up to 6 oz / 180 ml) to monitor your milk output.
5. **Finish Pumping:** Once pumping is complete, carefully turn off your breast pump.
6. **Remove and Pour:** Gently remove the cups from your bra. The self-standing design allows you to place them on a flat

surface. Carefully open the top cover and pour the collected milk into a sterile storage container.

The included lanyard and carrying rope can be used to secure the cups during activities, enhancing portability.

CLEANING AND MAINTENANCE

Regular cleaning and sterilization are crucial for hygiene and product longevity.

1. **Disassemble:** After each use, fully disassemble all parts of the milk collector cup.
2. **Rinse:** Rinse all components immediately with cold water to remove milk residue.
3. **Wash:** Wash all parts in warm, soapy water using a dedicated brush for breast pump parts. Ensure all milk film is removed.
4. **Rinse Thoroughly:** Rinse all parts thoroughly under clean running water.
5. **Sterilize:** The components are high temperature-resistant and boiling water-resistant. You can sterilize them using:
 - **Boiling:** Place parts in a pot of boiling water for 5 minutes.
 - **Steam Sterilizer:** Use a steam sterilizer according to its instructions.
 - **Dishwasher:** The cups are dishwasher safe. Place parts on the top rack of your dishwasher.
6. **Dry:** Allow all parts to air dry completely on a clean paper towel or drying rack. Avoid using cloth towels as they may leave lint.
7. **Store:** Store dry components in a clean, covered container until next use.

Easy to Clean



Image: The milk collector cup being rinsed, illustrating the ease of cleaning.

TROUBLESHOOTING

If you encounter issues with your Phanpy Milk Collector Cups, consider the following:

- **Reduced Suction:**

- Ensure all parts are correctly assembled and securely sealed. A loose connection can cause suction loss.
- Check the duckbill valve for tears or damage. Replace if necessary.
- Verify that the tubing is properly connected to both the cup and the pump unit, and that there are no kinks or blockages.

- **Leaking Milk:**

- Confirm that the breast shield is properly sealed against the Newcup and that the top cover is securely in place.
- Inspect the silicone flange and duckbill valve for any damage or wear that could compromise the seal.
- Ensure the cup is not overfilled beyond its maximum capacity (6 oz / 180 ml).

- When pouring, ensure the top cover is fully sealed to prevent spills.
- **Discomfort During Pumping:**
 - Ensure your nipple is centered in the flange.
 - If the 24mm flange feels too large, try using the included 20mm insert. Correct flange sizing is crucial for comfort and efficiency.
 - Adjust your breast pump's suction settings to a comfortable level.

If problems persist, please contact Phanpy customer support for assistance.

SPECIFICATIONS

Feature	Detail
Brand	Phanpy
Model Number	AM9951
Material	Silicone (BPA-free, food-grade)
Included Components	Milk Collector Cups (2 pieces), 24mm Flange (integrated), 20mm Insert
Single Cup Weight	0.25 Pounds (approx. 113 grams)
Single Cup Capacity	6 oz / 180 ml
Flange Angle	97° 3D simulated breastfeeding
Dishwasher Safe	Yes
Sterilization	High temperature-resistant, boiling water-resistant, steam sterilizer compatible

WARRANTY AND SUPPORT

For warranty information or customer support regarding your Phanpy Wearable Breast Pump Milk Collector Cup S2, please refer to the official Phanpy website or contact their customer service department directly.

Manufacturer: Phanpy

ASIN: B0CBRRBNMS

Date First Available: July 13, 2023

You can often find contact details on the product packaging or by visiting thePhanpy Brand Store on Amazon.

