



[Manuals.plus](#) /

› [LandTale](#) /

› LandTale Adjustable Shower Stool with Back User Manual

LandTale LT-B02

LandTale Adjustable Shower Stool with Back User Manual

Model: LT-B02

INTRODUCTION

The LandTale Adjustable Shower Stool with Back is designed to provide a safe, stable, and comfortable seating solution for individuals who require assistance while showering. This includes seniors, the elderly, disabled individuals, those recovering from injury, or pregnant women. Its robust design and adjustable features ensure a secure and relaxing shower experience.



Image: The LandTale Adjustable Shower Stool with Back, featuring a white seat and backrest, blue padded insert, and silver aluminum legs with blue anti-slip feet.

SAFETY INFORMATION

- **Weight Capacity:** Do not exceed the maximum weight capacity of 500 lbs (227 kg).
- **Stable Surface:** Always place the shower stool on a flat, level, and non-slip surface within the shower or bathtub.
- **Leg Adjustment:** Ensure all four legs are adjusted to the same height and the locking pins are securely engaged before use.
- **Regular Inspection:** Periodically check all components for wear, damage, or loose connections. Discontinue use if any damage is observed.
- **Cleaning:** Clean the stool regularly to prevent soap scum buildup, which can reduce friction.
- **Children:** Keep children away from the shower stool unless supervised by an adult.

PACKAGE CONTENTS

Verify that all components are present before beginning assembly:

- Shower Stool Seat with Padded Insert
- Backrest
- 4 x Height-Adjustable Legs with Anti-Slip Feet
- Assembly Hardware (pre-attached or minimal, tool-free design)

Tool-Free Assemble In 5 Minutes



Image: Close-up view of the underside of the shower stool seat, showing the tool-free knobs for attaching legs and crossbars.

SETUP AND ASSEMBLY

The LandTale Shower Stool is designed for quick and tool-free assembly, typically taking only a few minutes.

1. **Attach Legs:** Insert each of the four height-adjustable legs into the designated slots on the underside of the shower stool seat. Secure them by tightening the blue hand-turn knobs until firm.
2. **Attach Crossbars:** Connect the crossbars between the legs as indicated by the design, securing them with the remaining blue hand-turn knobs. These crossbars enhance stability.
3. **Install Backrest (Optional):** If desired, slide the backrest into the two vertical receivers at the rear of the seat. Ensure it clicks securely into place. The backrest is removable for user preference or space constraints.
4. **Adjust Height:** Before first use, adjust the height of all four legs to your desired level. Press the spring-loaded pins on each leg and slide the inner tube up or down until the pin locks into the desired hole. Ensure all legs are set to the same height for stability.

Removable BackRest

Choose What you need before shower



Image: Two shower stools, one with the backrest attached and one without, demonstrating the removable backrest feature.

Height Adjustment in 4 Level

Height Adjustable from 15" to 18"



Image: Two shower stools illustrating the height adjustment from 15 inches to 18 inches.

OPERATING INSTRUCTIONS

Using the Shower Stool

- Place the shower stool securely in your shower or bathtub. Ensure all four anti-slip feet are firmly on the surface.
- Carefully sit on the stool, ensuring your weight is centered.
- The soft, durable padded seat provides comfort during use.
- The built-in side handles on the seat facilitate easier lifting and repositioning of the stool.

Built-in Shower Head Holder

This shower stool features a convenient built-in holder for your shower head, allowing for hands-free showering.

Built-in Shower Head Holder



Image: The shower stool in a shower setting, with a shower head placed in the integrated holder on the side of the seat.

MAINTENANCE AND CARE

- **Cleaning:** Clean the shower stool regularly with mild soap and water. Rinse thoroughly and wipe dry to prevent water spots and mineral buildup. Avoid abrasive cleaners or harsh chemicals.
- **Drying:** Allow the stool to air dry completely after each use to maintain hygiene and prevent mildew.
- **Inspection:** Periodically check all connections, screws, and the condition of the anti-slip feet. Ensure the height adjustment pins are functioning correctly.
- **Storage:** When not in use, store the shower stool in a dry place. Its compact size allows for easy storage in smaller bathrooms or tubs.

TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Stool feels wobbly or unstable.	Legs not adjusted to the same height; locking pins not fully engaged; uneven surface.	Ensure all four legs are set to the exact same height and the spring-loaded pins are fully extended and locked into the adjustment holes. Place the stool on a flat, level surface.
Stool does not fit in narrow tub.	Tub dimensions are smaller than the stool's footprint.	The LandTale stool is designed to fit smaller tubs. Verify your tub's internal dimensions against the stool's footprint (14" x 12.7" at base). The backrest can be removed to save space if needed.
Anti-slip feet are slipping.	Soap residue or excessive water on the shower floor; worn feet.	Ensure the shower floor is clean and free of soap residue. Clean the anti-slip feet regularly. If feet are worn, contact customer support for replacements.

Heavy Duty Shower Stool

Creative Structure with Stability and Sturdiness



Image: Diagram highlighting the 500LB weight capacity, 1.4mm thick aluminum tubes, and anti-slipped feet for enhanced stability.

Perfect Size Stool

For Small Shower Tub and Bathrooms



Image: Detailed dimensions of the shower stool, showing seat width (17.4"), seat depth (10.6"), and base footprint (14" x 12.7").

SPECIFICATIONS

Feature	Detail
Model Number	LT-B02
Brand	LandTale
Material	Aluminum frame, Plastic seat and backrest, Rubber anti-slip feet
Color	White seat/backrest with Blue padded insert and blue accents
Maximum Weight Capacity	500 lbs (227 kg)
Adjustable Height Range	15 inches to 18 inches (from floor to seat)
Seat Dimensions	Approximately 17.4 inches (width) x 10.6 inches (depth)

Feature	Detail
Base Footprint	Approximately 14 inches (width) x 12.7 inches (depth)
Item Weight	6.78 pounds (3.08 kg)
Assembly	Tool-Free Assembly

WARRANTY AND SUPPORT

LandTale is committed to providing high-quality products. For any questions, concerns, or assistance with your shower stool, please contact LandTale customer support through the retailer where the product was purchased. Please have your model number (LT-B02) and purchase information ready when contacting support.

For more information about LandTale products, you may visit the official LandTale store on Amazon: [LandTale Amazon Store](#).