

## Mobiclinic Fortuna Basic

# Mobiclinic Fortuna Basic Electric Patient Lift Instruction Manual

Model: QA-00561/04

## 1. INTRODUCTION

---

The Mobiclinic Fortuna Basic electric patient lift is designed to safely transfer individuals with reduced mobility between various positions, such as from a wheelchair to a bed or sofa. This device significantly reduces the physical strain on caregivers while ensuring the patient's comfort and safety. Its robust design and quality construction make it a reliable solution for facilitating patient care.



*Mobiclinic Fortuna Basic electric patient lift with harness.*

## 2. SAFETY INFORMATION

Before operating the Mobiclinic Fortuna Basic patient lift, please read and understand all safety instructions. Failure to do so may result in injury to the patient or caregiver.

- **Maximum Weight Capacity:** Do not exceed the maximum weight capacity of 135 kg (297 lbs).
- **Emergency Stop:** The lift is equipped with an emergency stop button. Familiarize yourself with its location and function. In case of an emergency, press this button immediately to halt all lift operations.
- **Secure Harness:** Always ensure the patient is securely positioned in the harness and all attachment points are correctly fastened before any transfer.
- **Locking Wheels:** Engage the safety brakes on the rear wheels when the lift is stationary or during patient transfer to prevent unintended movement.
- **Clear Path:** Ensure the area around the lift is clear of obstacles before moving the patient.
- **Supervision:** Never leave a patient unattended in the lift.
- **Professional Use:** This device should be operated by trained caregivers or individuals who have thoroughly read and understood this manual.



*The emergency stop button and non-slip features enhance safety during operation.*

### 3. PACKAGE CONTENTS

Upon unpacking, please verify that all components are present and undamaged:

- Main lift frame (pre-assembled base)
- Lift arm
- Electric actuator (motor)
- Harness support bar
- Control unit with cables
- Patient transfer harness
- Assembly tools (wrenches)
- User Manual

### 4. ASSEMBLY INSTRUCTIONS

---

The Mobiclinic Fortuna Basic patient lift comes with a pre-assembled base for simplified installation. Follow these steps for complete assembly. A visual guide is also provided in the video below.

1. **Attach the Base of the Structure:** Secure the main vertical support to the pre-assembled base using the provided bolts and wrenches.
2. **Install the Lift Arm:** Attach the main lift arm to the top of the vertical support, securing it with the appropriate bolt.
3. **Install the Electric Actuator:** Connect the electric actuator (motor) to the lift arm and the vertical support. Ensure it is firmly in place.
4. **Attach the Harness Support Bar:** Fix the horizontal bar, which holds the harness, to the end of the lift arm.
5. **Attach the Harness:** Hook the patient transfer harness onto the four designated hooks on the harness support bar.
6. **Integrate the Control Unit:** Mount the control unit onto the main structure.
7. **Connect Control and Power Cables:** Plug the actuator cable, power cable, and control cable into their respective ports on the control unit. Ensure all connections are secure.

Your browser does not support the video tag.

*Official Mobiclinic assembly guide for the Fortuna Basic electric patient lift, demonstrating step-by-step installation.*



*The lift features lightweight materials, adjustable height, and a secure harness for patient safety.*

## 5. OPERATING INSTRUCTIONS

Operating the Mobiclinic Fortuna Basic patient lift is straightforward with its intuitive controls.

- **Power On:** Ensure the lift is connected to a power source and the battery is charged.
- **Positioning the Patient:** Carefully position the patient in the transfer harness while they are seated or lying down. Ensure the harness straps are evenly distributed and securely attached to the lift's hooks. The lift's 4 handles provide maximum safety and a greater sense of space for the user.
- **Raising the Patient:** Use the control unit to slowly raise the patient. Monitor the patient's comfort and position throughout the process. The lift arm can rotate 180 degrees for flexible maneuvering.
- **Moving the Lift:** Once the patient is clear of the original surface, disengage the wheel brakes and carefully push the lift to the desired location. The wheels are designed for easy movement, and the base can adapt to different door widths.
- **Lowering the Patient:** Position the patient over the target surface (bed, chair). Engage the wheel brakes. Use the control unit to slowly lower the patient onto the surface. Once the patient is safely seated, unhook the harness.

# 4 POIGNÉES POUR UNE SÉCURITÉ MAXIMALE

Donne à l'utilisateur une plus grande sensation d'espace



*The four handles and safety lock on the wheels ensure maximum safety and comfort during transfers.*



*The lift's design allows for easy movement and navigation through various doorways in a home environment.*

## 6. MAINTENANCE

Regular maintenance ensures the longevity and safe operation of your Mobiclinic Fortuna Basic patient lift.

- **Cleaning:** Clean the lift regularly with a damp cloth and mild disinfectant. Avoid abrasive cleaners or solvents.
- **Battery Charging:** The lift operates on two batteries. Ensure batteries are fully charged before use. Connect the charger to the control unit and a power outlet. A full charge typically takes several hours. The lift can be used while charging.
- **Component Inspection:** Periodically inspect all bolts, nuts, and moving parts for tightness and wear. Check the harness for any signs of damage or fraying. Replace worn or damaged parts immediately.
- **Wheel Functionality:** Ensure the wheels rotate freely and the brakes engage properly. Keep wheels free of debris.

## 7. TROUBLESHOOTING

If you encounter any issues with your Mobiclinic Fortuna Basic patient lift, refer to the following common troubleshooting steps:

- **Lift Not Responding:**

- Check if the emergency stop button is pressed. Release it by twisting or pulling, as indicated on the button.
- Ensure all cables (actuator, power, control) are securely connected to the control unit.
- Verify that the battery is charged and the power cable is properly plugged into a working electrical outlet.

- **Unusual Noises:** If the lift makes unusual noises during operation, stop use immediately and inspect all moving parts for obstructions or damage. Contact customer support if the issue persists.

- **Harness Issues:** If the harness appears damaged or is not holding the patient securely, discontinue use and replace the harness.

For issues not covered here, please contact Mobiclinic customer support.


## 8. SPECIFICATIONS


Feature	Specification
Brand	Mobiclinic
Model Number	QA-00561/04
Color	Silver
Item Weight	38 Kilograms
Maximum Supported Weight	135 kg
Harness Included	Yes
Batteries	Two batteries
Adjustable Height	90-174 cm
Base Width	55.5 - 80 cm
Length	120 cm
Harness Dimensions	88 cm x 85 cm
Motor	24V/6000N
Turning Radius	127 cm
Crane Arm	130 cm wide, 180° rotation
Noise Level	53 dB (A)
ASIN	B0CDL7LZF8
Date First Available	August 3, 2023



*Detailed dimensions and weight capacity of the Mobiclinic Fortuna Basic patient lift.*

## FORTUNA BASIC

 Poids du produit: 38kg  
Poids supporté: 135 kg

Non pliable 

Montage de  
la structure



Hauteur réglable  
Jusqu'à 174 cm

Abordable  
pour tous





Dimensions  
120x80x174



## FORTUNA



 Poids du produit: 45kg  
Poids supporté: 135 kg

Pliable 

Montage de la structure  
sans outils

Hauteur réglable  
Jusqu'à 184 cm

PREMIUM  
QUALITY

Dimensions  
134x80x184



*Comparison of key features between the Fortuna Basic and Fortuna models.*

## 9. WARRANTY AND SUPPORT

Mobiclinic S.L. is a leading manufacturer of clinical and hospital furniture, technical aids, and orthopedics, with a history dating back to 1985. We are committed to providing high-quality and reliable products.

For warranty information, technical support, or to order spare parts, please refer to the contact information provided with your purchase documentation or visit the official Mobiclinic website.