

Pride Mobility SC7132

Pride Mobility Pursuit 2.0 SC7132 Mobility Scooter Instruction Manual

Model: SC7132

1. INTRODUCTION

This manual provides essential information for the safe and efficient operation, maintenance, and care of your Pride Mobility Pursuit 2.0 SC7132 Mobility Scooter. Please read this manual thoroughly before operating your scooter to ensure proper usage and to maximize its lifespan. The Pursuit 2.0 is designed for adults requiring mobility assistance, offering robust performance and comfort for various terrains.

2. KEY FEATURES

- **High Weight Capacity:** Supports up to 400 lbs.
- **Speed and Range:** Maximum speed of 9.3 mph and a range of up to 40 miles per charge with an upgradable dual lithium battery configuration.
- **Durable Construction:** Equipped with durable plastic shrouds for navigating various environments.
- **Premium Seating:** Memory foam seat with a 6-inch slider for adjustable support and comfort.
- **Adjustable Tiller with USB Charger:** Infinitely adjustable tiller for optimal driving position, includes a built-in USB charger.
- **Enhanced LED Lighting:** Full LED lighting system for visibility, including directional turn signals, brake lights, and ultra-bright headlights.
- **Optimal Traction:** Pneumatic tires mounted on split rims for exceptional handling and traction on all terrains.
- **Color Options:** Available in Orange, Gray, and Black.

Note: This product is an FDA Class II Medical Device, designed to aid individuals with mobility impairments.

3. PRODUCT OVERVIEW



Figure 1: Pride Mobility Pursuit 2.0 SC7132 Mobility Scooter (Black model).

PURSUIT^{2.0}



Figure 2: Key features of the Pursuit 2.0, including memory foam seat, durable plastic shrouds, tiller adjustment, front basket, FDA Class II Medical Device status, pneumatic tires, LED lighting, and front/rear CTS suspension.

PURSUIT²

INSERT YOURSELF: ANYWHERE

When it comes to adventure, you never hold back. So why should your mobility unit? Meet the Pursuit 2, built for when the road ends and the real journey begins.

The upgradable dual lithium batteries will have you traveling up to 40 miles per charge at a fast cruising speed of 9.3 mph.

The Pursuit 2 is the top choice for your next excursion, whether that's right down the sidewalk or miles down a dirt road.



FDA Class II
Medical Device*



40 MILE
RANGE*



LITHIUM-ION BATTERY
LONGER LIFE



9.3 MPH
TOP SPEED

PRIDE
MOBILITY
LIVE YOUR BEST®

Figure 3: The Pursuit 2.0 offers a 40-mile range, features a lithium-ion battery for longer life, and achieves a top speed of 9.3 mph.

4. SETUP AND INITIAL USE

4.1 Unpacking

Carefully remove all packaging materials. Inspect the scooter for any signs of damage during transit. Retain packaging for future transport or storage.

4.2 Battery Charging

Before first use, ensure the lithium battery is fully charged. Connect the 8 amp battery charger to the scooter's charging port and a standard electrical outlet. Refer to the separate battery charging instructions for detailed steps and safety precautions.

4.3 Seat and Tiller Adjustment

Adjust the seat and tiller to your comfort. The memory foam seat features a 6-inch slider for depth adjustment. The tiller is infinitely adjustable to find your preferred driving position. The headrest height can also be

adjusted.

5. OPERATING INSTRUCTIONS

5.1 Starting the Scooter

The scooter can be started in two ways:

- **Remote Start:** Double-press the lightning bolt button on the ignition fob.
- **Manual Start:** Insert the key into the ignition and turn clockwise.

5.2 Driving

Once powered on, control the lever to move forward and backward. The scooter can reach a maximum speed of 9.3 mph. Always operate the scooter at a safe speed appropriate for the environment and your skill level.

5.3 Lighting and Signals

The scooter is equipped with a full LED lighting package for visibility:

- **Headlight:** Ultra-bright headlight for illumination.
- **Turn Signals:** Directional turn signals for indicating turns.
- **Brake Lights:** Rear brake lights for increased safety.
- **Running Lights:** Front and rear LED running lights.

5.4 Storage

The Pursuit 2.0 offers versatile storage solutions:

- **Front Basket:** A detachable front basket for personal items.
- **Hidden Compartment:** A hidden compartment under the tiller.
- **Rear Storage:** A large enclosed storage area behind the seat, secured with keys for added safety.

5.5 USB Charging

A built-in USB charging port is available to keep your devices powered up while on the go.

5.6 Safety Belt

Always fasten the safety belt to protect your safety during operation.

6. MAINTENANCE

Regular maintenance is crucial for the longevity and safe operation of your scooter. Perform the following checks regularly:

- **Battery Care:** Keep the battery fully charged. Avoid deep discharges.
- **Tire Pressure:** Check tire pressure regularly to ensure optimal performance and safety.
- **Cleaning:** Clean the scooter regularly with a damp cloth. Avoid using harsh chemicals or abrasive cleaners.
- **Fasteners:** Periodically check all nuts, bolts, and fasteners for tightness.

- **Professional Service:** For any complex repairs or issues, contact an authorized Pride Mobility service technician.

7. TROUBLESHOOTING

If you encounter any issues with your scooter, refer to the following basic troubleshooting steps:

- **Scooter Not Starting:** Ensure the key is fully inserted and turned, or the remote start button is pressed correctly. Check battery charge level.
- **Reduced Performance:** Verify battery charge. Check tire pressure.
- **Unusual Noises:** Stop the scooter immediately and inspect for any loose parts or obstructions. Contact service if the noise persists.

For persistent problems or issues not covered here, please contact Pride Mobility customer support or an authorized service center.

8. SPECIFICATIONS

MODEL NAME	Pursuit 2.0
MODEL NUMBER	SC7132
WEIGHT CAPACITY	400 lbs.
MAXIMUM SPEED	Up to 9.3 mph
GROUND CLEARANCE	3.8" at transaxle
TURNING RADIUS	69.3"
OVERALL LENGTH	57.5"
OVERALL WIDTH	29.7"
SEAT-TO-GROUND HEIGHT	21.9" - 24.9"
SEAT-TO-DECK HEIGHT	17.9"- 20.9"
FRONT TIRES	13" Pneumatic
REAR TIRES	13" Pneumatic
SUSPENSION	Front and Rear (CTS Suspension)
RANGE PER CHARGE	Up to 40 miles with upgradable dual lithium battery configuration
TOTAL WEIGHT	With standard battery: 242.6 lbs. Without batteries: 220.6 lbs.
HEAVIEST PIECE WHEN DISASSEMBLED	N/A (Unit does not disassemble)
STANDARD SEATING	Type: Limited reclining high back with sliders Weight: 54.6 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 18"
DRIVE SYSTEM	Rear-wheel drive, sealed transaxle, 24-volt DC motor
DUAL BRAKING SYSTEM	Regenerative and electromechanical
BATTERY REQUIREMENTS	25 volt, 50Ah lithium ion
BATTERY CHARGER	8 amp

Figure 4: Detailed specifications for the Pursuit 2.0 SC7132.

Pride Mobility Pursuit 2.0 SC7132 Specifications

Feature	Specification
Model Name	Pursuit 2.0
Model Number	SC7132
Weight Capacity	400 lbs.

Feature	Specification
Maximum Speed	Up to 9.3 mph
Ground Clearance	3.8" at transaxle
Turning Radius	69.3"
Overall Length	57.5"
Overall Width	29.7"
Seat-to-Ground Height	21.9" - 24.9"
Seat-to-Deck Height	17.9" - 20.9"
Front Tires	13" Pneumatic
Rear Tires	13" Pneumatic
Suspension	Front and Rear (CTS Suspension)
Range Per Charge	Up to 40 miles with upgradable dual lithium battery configuration
Total Weight (with standard battery)	242.6 lbs.
Total Weight (without batteries)	220.6 lbs.
Heaviest Piece When Disassembled	N/A (Unit does not disassemble)
Standard Seating	Type: Limited reclining high back with sliders Weight: 54.6 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 18"
Drive System	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements	25 volt, 50Ah lithium ion
Battery Charger	8 amp

9. SAFETY INFORMATION

Always prioritize safety when operating your mobility scooter. Adhere to all local traffic laws and regulations. Ensure you are familiar with the controls and capabilities of your scooter before operating it in public areas.

- Always wear appropriate safety gear, such as a helmet, if recommended by local regulations or personal preference.
- Avoid operating the scooter under the influence of alcohol or medication that may impair your judgment.
- Be aware of your surroundings and potential hazards, especially when navigating uneven terrain or inclines.
- Do not exceed the maximum weight capacity of 400 lbs.
- Ensure the scooter is stable before mounting or dismounting.

10. WARRANTY AND SUPPORT

For warranty information, please refer to the separate warranty document included with your purchase or visit the official Pride Mobility website. For technical support, parts, or service inquiries, please contact Pride Mobility customer service directly. Keep your purchase receipt and model number handy when contacting support.

Pride Mobility Products Corporation

Manufacturer of the Pursuit 2.0 SC7132 Mobility Scooter.