

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [Sonmol](#) /

› [Sonmol Spirometer User Manual](#)

Sonmol SMPF-2S

Sonmol Spirometer User Manual

Model: SMPF-2S

For monitoring PEF & FEV1 in adults and children

1. INTRODUCTION

The Sonmol Spirometer is an advanced electronic device designed for personal lung function monitoring. It accurately measures Peak Expiratory Flow (PEF) and Forced Expiratory Volume in 1 second (FEV1), providing essential data for managing respiratory conditions such as asthma and COPD. This portable device features Bluetooth connectivity for seamless data synchronization with a tracking application, allowing users to monitor their lung health trends over time and share results with healthcare providers.

The device incorporates a differential-pressure-type sensor for ultra-accurate measurements, meeting the U.S. Department of Health and Human Services' Guidelines for Asthma Diagnosis and Management. It also includes a risk indicator and allows for customization of normal expected values, providing immediate feedback on lung function control conditions.




Image: Sonmol Spirometer shown alongside its ePEF tracking application on a smartphone, illustrating the device's connectivity and data display capabilities.

2. KEY FEATURES


- **Ultra-Accurate Measurement:** Equipped with a differential-pressure-type sensor, the spirometer precisely measures PEF and FEV1, adhering to established medical guidelines.
- **Smart Risk Indicator:** The large green backlit LCD screen displays test values and provides a risk indicator (red, yellow, green zones) to reflect lung function control. Personal expected values can be set for tailored feedback.
- **Bluetooth Connectivity & App Analysis:** Connects to the ePEF tracking app via Bluetooth for automatic data synchronization, allowing for comprehensive tracking, analysis, and sharing of respiratory health data with

caregivers.


- **Portable & Rechargeable:** Compact and lightweight, the device is rechargeable via USB, offering over 6 months of use on a single full charge. Ideal for on-the-go monitoring.
- **300-Record Memory:** Stores up to 300 test records with date and time stamps, enabling users to track long-term trends.
- **Dual Mouthpiece Design:** Includes two types of removable and washable mouthpieces for user comfort and hygiene.

Greater Precision and Higher Repeatability

Our peak flow meter uses a differential pressure sensor, which calculates the variance by comparing the pressure at two points, thus providing accurate pressure measurements.



Turbine Peak Flow Meter



Mechanical Peak Flow Meter




Image: An internal view of the spirometer highlighting the high-precision differential pressure sensor, which ensures accurate airflow measurement.

3. SETUP GUIDE

3.1 Initial Charging

Before first use, fully charge the device for approximately 2 hours using the provided USB cord and a 5V adapter. The battery icon on the display will indicate charging status. When the battery icon shows 'x', it requires charging.



Image: The product packaging, which includes the USB charging cable.

3.2 App Download and Installation

The Sonmol Spirometer works with the ePEF tracking application. Search for "Asthma PEF", "Sonmol PEF", or "ePEF" in your device's app store (Google Play for Android or App Store for iOS) and download it. Alternatively, scan the QR code provided in the product documentation or on the device packaging to download the app.

Your browser does not support the video tag.

Video: Instructions on how to connect the spirometer via Bluetooth for Android devices. This video demonstrates launching the ePEF app and initiating the Bluetooth connection process.

Your browser does not support the video tag.


Video: Instructions on how to connect the spirometer via Bluetooth for iOS devices. This video guides users through the app launch and Bluetooth pairing steps specific to iOS.

3.3 Bluetooth Connection


Launch the ePEF app on your smartphone. Ensure the Peak Flow Meter is turned on. Press the device icon within the app to initiate connection. Once successfully connected, the icon in the app will turn green, and the spirometer's display will show a Bluetooth symbol (*). Data can then be automatically synced to your mobile phone.

Sync Data Via Bluetooth Connections


Press the icon to connect







Press the icon to connect



Press to share



Launch ePEF app, press the icon  on your phone to connect with the Peak Flow Meter turned on. The icon will turn into  green and the peak flow meter will show  once the bluetooth connection is successfully established. Press the icon  again to sync data.

*** All data stored in the Peak Flow Meter can be synced to your phone at one time once Bluetooth is connected successfully.**

Image: A visual guide to syncing data via Bluetooth, showing the app interface and the device's display during connection.

3.4 Setting Expected PEF Value

Setting your personal expected PEF value is crucial for the risk indicator to function correctly. This can be done via auto-configuration through the ePEF app or manually on the device.

- **Auto-configuration:** The device will automatically configure based on your latest expected PEF values set in the ePEF app once connected via Bluetooth 4.0.
- **Manual Configuration:**
 1. With the device powered off, hold the Home button for 5 seconds until "000L" appears, indicating parameter setting mode. The flashing digit is editable.
 2. Press the Right button to increase the value or the Left button to decrease the value. Press the Home button to move to the next digit.
 3. Hold the Home button for 3 seconds to confirm and save your PEF value.

Set Expected Value

1. Auto configuration

The device will be automatically configured based on your latest expected PEF values set in the ePEF app once connected via Bluetooth 4.0.

2. Manual configuration

- 1) Enter Setting Mode: With the device powered off, hold the Home button for 5 seconds until "000L" appears, indicating parameter setting mode. The flashing digit is editable.
- 2) Adjust Value:
Press the Right button to increase the value.
Press the Left button to decrease the value.
Press the Home button to move to the next digit.
- 3) Save Settings: Hold the Home button for 3 seconds to confirm and save your PEF value.



Home Button

Left Button
(Page up, decreasing the value)

Right Button
(Page down, increasing the value)

Spiromol

Image: Instructions for manually setting the expected PEF value on the spirometer's display using its buttons.

Your browser does not support the video tag.

Video: A vertical video demonstrating how to set the expected value within the ePEF App, which then syncs to the device.

Your browser does not support the video tag.

Video: A video tutorial on how to manually set the time on the spirometer device.

4. OPERATING INSTRUCTIONS

4.1 Performing a Measurement

1. Ensure the device is powered on and ready (display shows "go").
2. Inhale deeply with your mouth.
3. Cover the mouthpiece completely with your mouth, ensuring a tight seal.
4. Exhale as fast and as hard as you can into the mouthpiece.
5. Wait for two beeps from the device before reading the result.
6. Repeat the measurement at least three consecutive times, taking the highest value as the final result.

2 Shapes of Mouthpieces Better Meet Your Needs

Please ensure your mouth completely seals the mouthpiece



Image: A user demonstrating the correct technique for using the spirometer with one of the provided mouthpieces.

Your browser does not support the video tag.

Video: A comprehensive guide on how to properly use the peak flow meter for accurate measurements.

4.2 Interpreting Results with Risk Indicator

The device's screen will automatically display your PEF and FEV1 values. A triangular symbol will point to a color-coded zone (red, yellow, or green) to indicate your lung function control condition based on your set expected value.

- **Green Zone (PEF \geq 80% of Expected):** Indicates good lung function management.
- **Yellow Zone (50% \leq PEF $<$ 80% of Expected):** Indicates caution is needed. Respiratory airways may be narrowing, and additional medication might be required.
- **Red Zone (PEF $<$ 50% of Expected):** Indicates a medical emergency. Severe airway narrowing may be occurring, requiring immediate action and contact with a doctor or hospital.

**Note: No risk indicator will be displayed on the device without properly setting the expected value.*

Smart Risk Indicator

Each test result will be accompanied by a triangular symbol "▼" pointing to the red, yellow, or green sections, which provides insights into your lung function control condition

The "▼" point to different colors indicates different control results:

Green

PEF% \geq 80%, a peak flow reading in the green zone indicates that the lung function management is under good control.

Yellow

50% \leq PEF% $<$ 80%, your peak flow reading indicates caution.

Red

PEF% $<$ 50%, your peak flow reading indicates a medical emergency.



* No red, yellow, or green Risk Indicators will be displayed on the device without properly setting the expected value.

Image: The spirometer's display illustrating the color-coded risk indicator zones (green, yellow, red) for easy interpretation of lung function.

How does the risk indicator work?

Each test result will be accompanied by a triangular symbol “▼” pointing to the red, yellow, or green sections, which provides insights into your lung function control condition.

* No arrow will be displayed on the device without properly setting the expected value.

What do the colors mean?



Image: A visual chart explaining the meaning of the green, yellow, and red zones of the risk indicator based on PEF percentage.

4.3 Getting Your Expected PEF Value

The expected PEF value is a personal standard lung function value, typically calculated based on statistical methods considering race, gender, age, and height. Since individual physiological conditions vary, there will be differences between actual tested values and expected values. You can obtain your expected value through:

- **Calculation Formula:** Use the calculation formula provided by a specialized research institution. Scan the QR code below or visit [Sonmol PEF Calculator](#) for more information.
- **Physician Recommendation:** Your physician may provide a recommended value for you.
- **Historical Data:** Use the best of your historical PEF values from the past 2 weeks.

1. What is Expected Value?

The expected value is a personal standard lung function value, which is usually calculated based on statistics methods. It varies according to race, gender, age, and height. Since everyone's physiological conditions is different, there is a variation between the actual tested value and the expected value.

2. How to Get Your Expected Value?

You can get your expected value in the following ways:

- 1 Calculation formula provided by a specialized research institution, scan to calculate your expected value. Please read the instructions for more information.
- 2 Your physician may give a recommend value for you.
- 3 The best of your historical PEF values (in 2 weeks).



Image: Information on how to determine your expected PEF value, including a QR code to access an online calculator.

5. MAINTENANCE

5.1 Cleaning the Mouthpiece

The spirometer comes with two types of removable and washable mouthpieces. To clean, simply wipe them with

alcohol. Ensure they are completely dry before reattaching to the device. Regular cleaning prevents the accumulation of contaminants.

5.2 Storage

Store the device in a clean, dry place away from direct sunlight and extreme temperatures. Keep it out of reach of children.

6. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device not turning on	Low battery; Device malfunction	Charge the device for 2 hours. If problem persists, contact support.
Bluetooth not connecting	Device not on; App not open; Bluetooth off on phone; Distance too far	Ensure device is on and app is open. Check phone's Bluetooth settings. Keep device close to phone.
No risk indicator displayed	Expected PEF value not set	Set your expected PEF value in the app or manually on the device (refer to Section 3.4).
Inaccurate readings	Improper technique; Mouthpiece not sealed; Device needs cleaning	Ensure proper technique (refer to Section 4.1). Check mouthpiece seal. Clean mouthpiece.

7. SPECIFICATIONS

- **Measures:** PEF (Peak Expiratory Flow), FEV1 (Forced Expiratory Volume in 1 second)
- **Sensor Type:** Differential-pressure-type sensor
- **Display:** Large green backlit LCD screen
- **Connectivity:** Bluetooth 4.0
- **Memory:** 300 test records with date & time setting
- **Power:** Rechargeable battery (USB cord included)
- **Battery Life:** More than 6 months on a single full charge
- **Mouthpieces:** 2 types, removable and washable
- **Model Number:** SMPF-2S
- **Manufacturer:** Sonmol
- **Package Dimensions:** 7.4 x 2.95 x 2.44 inches
- **Item Weight:** 7.05 ounces

8. WARRANTY AND SUPPORT

For detailed warranty information and customer support, please refer to the official user manual included in your product packaging or visit the Sonmol official website. You can also find additional resources and FAQs on the ePEF

app.

User Manual (PDF): [Download Here](#)

User Guide (PDF): [Download Here](#)



© 2025 Sonmol. All rights reserved.