

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [HaYiue](#) /
- › [HaYiue Ear Wax Removal Tool Camera 1296P HD - Instruction Manual](#)

## HaYiue B0CLXZ2WPZ

# HaYiue Ear Wax Removal Tool Camera 1296P HD - Instruction Manual

## INTRODUCTION

---

Welcome to your new HaYiue Ear Wax Removal Tool Camera. This device is designed to provide a clear, visual method for inspecting and cleaning the ear canal. Please read this manual thoroughly before use to ensure safe and effective operation.

## PRODUCT OVERVIEW

---

### Key Features:

- **1296P HD Camera:** Equipped with a 360° wide-angle professional-grade lens for high-definition imaging.
- **6 LED Lights:** Enhances brightness and sharpness of the image within the ear canal.
- **3.5mm Ultra-Slim Lens:** Designed to reach the eardrum and deep ear canal comfortably, suitable for various ear sizes.
- **IP67 Waterproof:** Camera tip is waterproof for easy cleaning.
- **Rechargeable Battery:** Built-in 230mAh lithium-ion battery, rechargeable via USB Type-C.
- **Wireless Connectivity:** Compatible with most iOS and Android smartphones and tablets.
- **Non-Heating Chip:** Maintains a constant temperature to prevent discomfort during use.

### Package Contents:

- HaYiue Ear Wax Removal Tool Camera
- Silicone Ear Spoons (various types)
- USB Type-C Charging Cable
- Traditional Earwax Removal Kits (6 pieces)
- Cleaning Brush



This image displays the complete HaYiue Ear Wax Removal Tool Camera kit, featuring the main camera device, a variety of ear picks, a USB-C charging cable, and several soft silicone ear spoon attachments. A smartphone screen showing the camera's view inside an ear is also visible.



This image details the various ear picks included in the kit, such as spiral, hard spiral head, and traditional ear picks, along with a cleaning brush for maintenance.

## SETUP GUIDE

---

### 1. Download the App:

Search for "EndscopeTool" on the Apple App Store (for iOS) or Google Play Store (for Android) and install the application on your smartphone or tablet.

### 2. Power On the Device:

Press and hold the power button on the ear wax removal tool for approximately 2 seconds to turn it on. A blue indicator light will illuminate, indicating the device is active.

### 3. Connect to Wi-Fi:

- Go to your smartphone's Wi-Fi settings.
- Find and connect to the Wi-Fi network named "AiSee-XXX" (where XXX is a unique identifier for your device).
- Return to the "EndscopeTool" app and click "Start" to begin viewing the live camera feed.



A visual guide detailing the four simple steps to set up the HaYiue ear camera: downloading the 'EndscopeTool' app, pressing the power button, connecting to the device's Wi-Fi network, and initiating the camera view within the app.

## OPERATING INSTRUCTIONS

---

### 1. Attach an Ear Spoon:

Gently attach a suitable silicone ear spoon to the tip of the camera. Ensure it is securely fastened and will not fall off during use.

### 2. Position the Camera:

Carefully insert the camera into the ear canal while observing the live feed on your smartphone screen. Adjust the angle and depth slowly to get a clear view of the ear canal and eardrum.

### 3. Earwax Removal:

Using the attached ear spoon, gently scoop or dislodge earwax. Always move slowly and avoid pushing wax deeper into the ear canal. Do not insert the device too deeply or apply excessive force to prevent injury.

### 4. Inspection:

The device can also be used for general inspection of the ear, nose, mouth, throat, and other areas as needed.

## HaYiue Ear Wax Removal with Camera

Best choice for your earwax removal

- 1296 HD
- WiFi Connection
- LED lights
- IP67 Waterproof
- Type-C Charging
- iOS / Android

This image depicts an adult gently using the HaYiue ear camera to inspect a sleeping child's ear, demonstrating the device's suitability for sensitive use.



A collage showcasing the diverse applications of the HaYiue ear camera, including detailed inspections of human ears, nose, mouth, and even a pet's ear, highlighting its versatility.

## MAINTENANCE

### Cleaning the Device:

- The camera tip is IP67 waterproof. After each use, detach the ear spoon and clean the camera tip with water or wipe it with alcohol.
- Ensure the device is completely dry before storage.
- Clean the silicone ear spoons and traditional ear picks with soap and water or alcohol after each use.
- Use the provided cleaning brush to remove any residual debris from the ear spoons.

### Charging:

Connect the USB Type-C cable to the charging port on the device and to a compatible power source. A full charge provides approximately 90 minutes of continuous use. The blue indicator light will show charging status.



This image highlights the Type-C fast charging capability, the six high-brightness LED lights for clear visibility, and the stable WiFi connection feature of the HaYiue ear camera.

## TROUBLESHOOTING

### Device not connecting to phone:

- Ensure the device is powered on and the blue indicator light is active.
- Verify that your phone's Wi-Fi is enabled and connected to the "AiSee-XXX" network.
- Restart both the ear wax removal tool and your smartphone, then attempt to connect again.
- Ensure the "EndscopeTool" app has necessary permissions (e.g., Wi-Fi, camera access) on your phone.

### Poor image quality:

- Check if the camera lens is clean. Wipe it gently with a soft, dry cloth.
- Ensure there is sufficient lighting from the built-in LEDs.
- Adjust the focus within the app if available.

### Device not charging:

- Verify that the USB Type-C cable is securely connected to both the device and the power source.
- Try using a different USB Type-C cable or power adapter.
- Ensure the charging port is free of debris.

## SPECIFICATIONS

---

- **Model:** B0CLXZ2WPZ
- **Camera Resolution:** 1296P HD
- **Lens Diameter:** 3.5mm
- **LED Lights:** 6
- **Waterproof Rating:** IP67
- **Battery Capacity:** 230mAh Lithium-Ion
- **Charging Port:** USB Type-C
- **Compatibility:** iOS and Android
- **Manufacturer:** SHENZHEN JIDING ELECTRONIC CO.,LTD HaYiue
- **Package Dimensions:** 6.69 x 2.56 x 0.75 inches
- **Item Weight:** 2.08 ounces

## SAFETY INFORMATION

---

**Legal Disclaimer:** Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. *(Note: This disclaimer is provided as per source data, though it may not directly apply to this product.)*

### General Safety Precautions:

- Do not use the device if you have an ear infection, ear pain, discharge, or any other ear condition without consulting a medical professional.
- Always use the device gently and avoid inserting it too deeply into the ear canal to prevent injury to the eardrum or ear canal.
- Keep out of reach of children. Adult supervision is recommended when the device is used on children.
- Do not share ear spoons without proper sterilization to prevent cross-contamination.
- If you experience any discomfort, pain, or unusual sensations during use, stop immediately and consult a doctor.

## WARRANTY AND SUPPORT

---

HaYiue products are designed for quality and reliability. If you encounter any issues or have questions regarding your device, please contact the HaYiue support team. We are committed to resolving any concerns you may have.

For support, please refer to the contact information provided with your product packaging or visit the official HaYiue website for further assistance.

