

HOMCOM 713-042V80BK

HOMCOM Power Lift Recliner Chair for Elderly

Model: 713-042V80BK

INTRODUCTION

This manual provides essential instructions for the safe assembly, operation, and maintenance of your HOMCOM Power Lift Recliner Chair. Please read this manual thoroughly before use and retain it for future reference. This power lift recliner is designed to assist individuals, particularly the elderly or those recovering from injuries, in transitioning between sitting and standing positions with ease, and to provide comfortable reclining for relaxation.



Image: The HOMCOM Power Lift Recliner Chair in an upright position, showcasing its design and fabric.

SAFETY INFORMATION

Important Safety Warnings:

- Do not exceed the maximum weight limit of **265 lbs (120 kg)** for safe use. Exceeding this limit may cause damage to the chair and potential injury.
- Ensure all connections are secure before operating the chair.
- Keep children and pets away from the chair's moving parts during operation.
- Avoid placing hands or feet near the mechanism under the chair to prevent injury.
- This product is certified to meet safety standards. Any modifications to the chair or its components may void the warranty and create safety hazards.
- Unplug the chair from the power outlet before cleaning or maintenance.

PACKAGE CONTENTS

Before beginning assembly, verify that all parts are present and undamaged. The package should contain:

- Recliner Base
- Backrest
- Arm Pads (if separate)
- Remote Controller
- Power Adapter and Cable
- Instruction Manual

If any parts are missing or damaged, do not proceed with assembly. Contact customer support for assistance.

ASSEMBLY INSTRUCTIONS

Assembly of this recliner chair is designed to be straightforward. No special tools are required.

1. **Unpack Components:** Carefully remove all parts from the packaging. Place them on a soft, clean surface to prevent scratches.
2. **Attach Backrest:** Align the metal brackets on the backrest with the corresponding slots on the recliner base. Slide the backrest down firmly until it clicks into place. Ensure it is securely fastened.
3. **Connect Power Cables:** Locate the power cables extending from the backrest and the base. Connect them together, ensuring a snug fit.
4. **Connect Remote Control:** Plug the remote control cable into the designated port on the chair's power system.
5. **Connect Power Adapter:** Plug the power adapter into the chair's power input, then plug the adapter into a wall outlet.
6. **Final Check:** Gently test the chair's stability and functions before full use.

For visual guidance, refer to the assembly diagrams provided in the included printed manual.

OPERATING INSTRUCTIONS

The HOMCOM Power Lift Recliner Chair is operated using the included remote control.

Remote Control Functions:

- **Up Button:** Press and hold to activate the lift mechanism, gently raising the chair to an almost standing position (up to 45 degrees). Release the button to stop at any desired height.
- **Down Button:** Press and hold to recline the chair backward and raise the footrest. The chair can recline up to 135 degrees. Release the button to stop at any desired angle.

ELECTRIC RECLINER



Image: A person demonstrating the lift assistance feature, where the chair gently tilts forward to help the user stand up.

PADDED COMFORT



Image: A visual representation of the chair's reclining capabilities, illustrating positions from upright to a 135-degree recline.

Video Demonstration:

Your browser does not support the video tag.

Video: Official demonstration of the HOMCOM Power Lift Chair's lift and recline functions, showing smooth transitions and ease of use.

The chair features convenient side pockets for storing the remote control and other small items.



Image: Detailed view highlighting the chair's side pockets, the retractable footrest mechanism, and the robust steel frame.

MAINTENANCE

Proper maintenance ensures the longevity and optimal performance of your recliner chair.

- **Cleaning:** Wipe clean with a damp cloth. Avoid harsh chemicals or abrasive cleaners, as these can damage the fabric. For spills, blot immediately with a clean, dry cloth.
- **Mechanism Check:** Periodically inspect the lift and recline mechanisms for any signs of wear or loose connections. Do not lubricate the mechanism unless specifically instructed by customer support.
- **Power Cord:** Ensure the power cord is not pinched, twisted, or placed under heavy objects.
- **Fabric Care:** The chair features breathable linen-feel fabric. Avoid direct sunlight exposure for prolonged periods to prevent fading.

TROUBLESHOOTING

If you encounter issues with your recliner, refer to the following common solutions:

- **Chair Not Responding:**

- Check if the power cord is securely plugged into both the chair and the wall outlet.
 - Ensure the remote control cable is firmly connected to the chair.
 - Verify that the wall outlet is functional by plugging in another appliance.
- **Unusual Noises:**
 - Minor operational noises are normal. If loud or grinding noises occur, discontinue use and contact customer support.
 - **Incomplete Movement:**
 - Ensure no obstructions are blocking the chair's movement or the footrest.
 - Check that the weight limit is not exceeded.

For issues not listed here, or if solutions do not resolve the problem, please contact HOMCOM customer support.

SPECIFICATIONS

Detailed specifications for the HOMCOM Power Lift Recliner Chair (Model: 713-042V80BK):

Feature	Detail
Brand	HOMCOM
Model Name	Power Lift Chair
Overall Dimensions (Upright)	26" W x 35.5" D x 37.75" H
Overall Dimensions (Reclined)	61.75" L x 26" W x 32" H
Seat Size	18.5" W x 21.5" D x 17.75" H
Weight Capacity	265 lbs (120 kg)
Material	100% Polyester Fabric, Foam Fill, Steel Frame
Special Features	Power Lift, Power Reclining, Retractable Foot Rest, Side Pocket, Spring Pack Seat
Reclining Positions	Up to 135 degrees
Lift Angle	Up to 45 degrees
Item Weight	75 pounds

HELP TO STAND

Especially for the elderly and the injured in rehabilitation



Image: Detailed dimension diagram of the recliner chair, showing height, width, depth, and seat measurements.

WARRANTY AND SUPPORT

For warranty information, please refer to the documentation included with your purchase or contact HOMCOM customer support directly. If you require assistance with assembly, operation, or troubleshooting, please reach out to our customer service team. Contact details can typically be found on the manufacturer's website or on your purchase receipt.