

## SuperHandy GBTS023

# SuperHandy GoRise FC Electric Floor to Chair Lift

## USER MANUAL

Brand: SuperHandy | Model: GBTS023

### Introduction

---

The SuperHandy GoRise FC Electric Floor to Chair Lift (Model: GBTS023) is designed to provide safe and effortless assistance for seniors and individuals with mobility challenges, helping them transition from a seated position on the floor to a chair height. This manual provides essential information for the proper setup, operation, and maintenance of your lift.

Key features of this product include:

- High capacity floor lifting, supporting up to 400 lbs, effortlessly lifting from 1.2 inches to 21 inches.
- Versatile dual power mode with DC battery and AC power options for uninterrupted assistance, including an LED recharge indicator.
- Ergonomic and safe design featuring foldable armrests, a sloped seating platform, and a secure safety belt.
- User-friendly operation with quiet, smooth lifting at 10mm/s speed and a magnetically attachable wired remote for easy control.
- Effortless mobility and storage, equipped with 5-inch silent caster wheels and a convenient top handle.

### Safety Information

## Important Safety Warnings

**This product has a maximum weight capacity of 400 lbs. Improper use may result in injury. Ensure the lift mechanism is used safely and handle with care. This product is designed for users with mobility issues; consult the user manual for specific precautions related to this user group. Always read and understand all instructions before operating the lift.**

## Components and Assembly

### Components

Your SuperHandy GoRise FC Electric Floor to Chair Lift package should include the following main components:

- Main Lift Unit with integrated seat and backrest
- Foldable Armrests (pre-attached or requiring attachment)
- Wired Remote Control with magnetic attachment
- Safety Belt
- Power Adapter/Charger
- User Manual (this document)

### Assembly

The SuperHandy GoRise FC Electric Floor to Chair Lift requires some assembly upon arrival. While the process is designed to be straightforward, it is recommended to follow the detailed instructions and diagrams provided in the physical manual that accompanies your product. Users have reported assembly taking approximately 15 minutes.





Figure 1: The SuperHandy GoRise FC Electric Floor to Chair Lift main unit.

## Setup

---

### Initial Setup Steps

1. Unpack the electric floor lift and ensure all parts are present.
2. Charge the battery fully before first use.
3. Place the lift on a flat, stable surface.

4. Ensure the person using the lift is seated comfortably with the safety belt fastened.
5. Use the wired remote control to lift the person smoothly to the desired height.
6. Move the lift (if necessary) to the target chair or seat.
7. Lower the person gently onto the chair or seat.
8. Turn off the lift when not in use.



## Need a Lift?

Effortlessly raising individuals from ground to chair height. A height range of **1.2 inches** to a comfortable **21 inches** above the ground.



**Lifting Capacity: 400lbs**

Figure 2: User on the floor, preparing for the lift. Ensure the lift is positioned correctly.

### Operating Instructions

---

The SuperHandy GoRise FC Electric Floor to Chair Lift is designed for intuitive operation, allowing both users and caregivers to manage transitions with ease.

Your browser does not support the video tag.

Video 1: A visual guide demonstrating the effective use of the SuperHandy GoRise FC Electric Floor to Chair Lift, highlighting its key features and operation.

### Using the Lift



## Floor to Sit

Get a help standing up or sitting down.

**Versatile Power Options:** Operates cordlessly using a rechargeable battery and can be used with a AC power (100~240V, 50/60Hz).

**Height Range: 1.2inches~21inches**

Figure 3: The lift effortlessly raises individuals from the ground to a comfortable chair height.

**Sloped Lift Seat:** The ergonomically designed sloped seat enhances comfort and posture by distributing weight evenly and reducing pressure on the lower back, making it easier for the user to position themselves for lifting.

# Comfort and Safety Features

Streamline your daily operations with the ergonomically designed sloped cushion, ensuring both comfort and convenience.

## Sloped Cushion Seat



## Foldable Armrest

Easily board onto the lift with ease.



## Quick Release Safety Belt



Figure 4: Key comfort and safety features including the sloped cushion seat, foldable armrest, and quick release safety belt.

**Quick Release Safety Belt:** For enhanced safety and security, always ensure the quick-release safety belt is properly fastened before initiating any movement. It is fully adjustable to provide a secure fit, ensuring comfort and stability during lifting.

**Foldable Armrests:** The armrests can be folded up, providing easier access for boarding the lift or for transfers. Once seated, they can be lowered to offer additional support and stability, streamlining daily operations.



## Hand Control

Magnetic remote, useable for more than 20 times on a single charge.

 Up  Down

Smooth and quiet lifting at **10mm/s.**

Figure 5: The magnetic remote control for easy and user-friendly operation.

**User-Friendly Control:** The lift is operated via a magnetically attachable wired remote control. This remote allows for smooth and quiet lifting and lowering at a speed of 10mm/s. The magnetic feature ensures the remote is always within reach and securely stored when not in use, usable for more than 20 times on a single charge.

**Dual Power Mode:** The lift offers versatile power options, operating on a DC battery for cordless use or by plugging directly into a wall outlet using AC power for uninterrupted assistance. An LED indicator shows the battery charge status.



## Effortless Mobility and Storage

Equipped with silent caster wheels and a top handle for effortless movement and storage.

Figure 6: The lift's design allows for effortless mobility and storage.

**Effortless Mobility and Storage:** Equipped with 5-inch silent caster wheels and a convenient top handle, the lift can be easily moved and stored when not in use, offering practical portability.



## Wide Spectrum of Uses

### Elderly Individuals

Assisting you in rising independently from the ground.

### Care Givers

Enabling effective support for patients in standing up.

### Leg Injuries

Aiding independent standing

Figure 7: The lift is suitable for elderly individuals, caregivers, and those recovering from leg injuries.

## Maintenance

Regular maintenance ensures the longevity and safe operation of your SuperHandy GoRise FC Electric Floor to Chair Lift.

- **Cleaning:** Wipe down all surfaces with a soft, damp cloth. Avoid using abrasive cleaners, solvents, or harsh chemicals, as these can damage the finish and components.
- **Battery Care:** Ensure the battery is fully charged before first use and regularly recharged to maintain optimal performance. Store the lift in a cool, dry place when not in use for extended periods.
- **Inspection:** Periodically check all bolts, connections, and moving parts for any signs of looseness, wear, or damage. Ensure the safety belt is in good condition and functions properly. If any issues are found, discontinue use and contact customer support.

## Troubleshooting

---

### Common Issues and Solutions

| Issue                              | Possible Solution  |
|------------------------------------|--|
| Lift not operating.                | Check if the power adapter is securely plugged into a working outlet or if the battery is sufficiently charged. Ensure the remote control is properly connected. |
| Lift moving slowly or erratically. | Recharge the battery fully. Check for any obstructions that might be impeding the lift's movement. Ensure the weight limit is not exceeded.                      |
| Unusual noises during operation.   | Discontinue use immediately. Inspect the lift for any visible damage or loose parts. Contact SuperHandy customer support for assistance.                         |

### Specifications

---

| Specification                 | Value  |
|-------------------------------|--|
| Brand                         | SuperHandy   |
| Model Name                    | SuperHandy Floor Lift 400Lbs (GBTS023)   |
| Product Dimensions            | 28"D x 22"W x 32"H   |
| Item Weight                   | 61.6 pounds  |
| Maximum Weight Recommendation | 400 Pounds   |
| Seat Height Range             | 1.2 inches (min) to 21 inches (max)  |
| Lifting Speed                 | 10mm/s   |
| Power Source                  | DC Battery (Lithium Ion, included) / AC Power                                      |
| Special Features              | Adjustable Height, Foldable, Magnetic Remote, Sloped Seating Platform, Safety Belt |

### Warranty and Support

---

For any questions, technical assistance, or warranty claims regarding your SuperHandy GoRise FC Electric Floor to Chair Lift, please contact SuperHandy customer support through the retailer where you purchased the product or visit the official SuperHandy website.

Optional 2-Year and 3-Year Protection Plans are available for extended coverage. Please refer to your purchase documentation for details on these plans and how to activate them.

Visit the official SuperHandy Store for more products and support:[SuperHandy Store](#)

© 2024 SuperHandy. All rights reserved.