

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [MAIDeSITE](#) /
- › [MAIDeSITE Electric Chair Lift - Base User Manual](#)

MAIDeSITE ce0778a9-fed6-4418-ad3d-35dc909f1231

MAIDeSITE Electric Chair Lift - Base User Manual

Model: ce0778a9-fed6-4418-ad3d-35dc909f1231

1. INTRODUCTION

This manual provides essential information for the safe and effective use of your MAIDeSITE Electric Chair Lift - Base. Please read these instructions thoroughly before assembly, operation, or maintenance. Keep this manual for future reference.

The MAIDeSITE Electric Chair Lift - Base is designed to assist individuals with limited mobility in sitting down and standing up from a seated position. It features an electric lifting mechanism for smooth and controlled movement.

2. SAFETY INFORMATION

- Always ensure the chair lift is placed on a stable, level, and hard floor surface.
- Do not exceed the maximum weight capacity specified in the 'Specifications' section.
- Keep hands and loose clothing clear of moving parts during operation.
- Ensure the power cord is not pinched or damaged. Do not operate with a damaged cord.
- Do not immerse the control unit or any electrical components in water.
- Supervise children and individuals who may not understand the safe operation of the device.
- Disconnect from power before cleaning or performing any maintenance.
- If the unit malfunctions, discontinue use immediately and contact customer support.

3. PACKAGE CONTENTS

Carefully unpack all components and verify that you have received the following items:

- Electric Chair Lift Base (main unit)
- Seat and Seat Padding
- Side Wings (2)
- Rails
- Cross Braces
- Locking Bar Rails
- Suction Cups (pre-attached or separate)

- Handheld Remote Control
- Power Adapter
- User Manual

If any items are missing or damaged, please contact MAIDeSITE customer support.

4. SETUP AND ASSEMBLY

The MAIDeSITE Electric Chair Lift - Base requires minimal assembly. Follow these steps to prepare your unit for use:

1. **Unpack the Unit:** Remove all packaging materials. Place the folded chair lift on a clean, flat surface.
2. **Unfold the Base:** Gently unfold the main base unit until it locks into its open position. Ensure all locking mechanisms are securely engaged.
3. **Attach Side Wings:** If the side wings are not pre-attached, slide them into the designated slots on either side of the seat until they click into place.
4. **Position the Unit:** Place the chair lift on the desired hard floor surface. Ensure the suction cups (if present) are clean and dry for optimal grip. Press down firmly to secure the unit.
5. **Connect Power:** Plug the power adapter into the designated port on the chair lift, then plug the adapter into a standard electrical outlet.



Figure 4.1: The MAIDeSITE Electric Chair Lift Base in its compact, folded position. This view shows the blue seat and white frame, ready for unfolding or storage.



Figure 4.2: An elderly man seated on the MAIDeSITe Electric Chair Lift Base, positioned next to a sofa, demonstrating its use in a home environment.

5. OPERATING INSTRUCTIONS

The chair lift is operated using the handheld remote control.

1. **Sitting Down:** Ensure the chair lift is in its lowest position. Carefully approach the unit and sit down on the seat. Use the armrests (if applicable) for support.
2. **Standing Up:** While seated, press and hold the 'UP' button on the remote control. The seat will slowly and smoothly elevate, assisting you to a standing position. Release the button once you are comfortably standing.
3. **Lowering Down:** While standing, carefully position yourself over the seat. Press and hold the 'DOWN' button on the remote control. The seat will slowly lower, assisting you to a seated position. Release the button once you are comfortably seated.
4. **Emergency Stop:** In case of an emergency, release the button on the remote control. The movement will stop immediately.

Operational Video

Video 5.1: Demonstration of the MAIDeSITe Electric Chair Lift in operation, showing the lifting and lowering functions.

Video 5.2: Another perspective of the MAIDeSITe Electric Chair Lift's functionality, highlighting its smooth movement.

6. MAINTENANCE

- **Cleaning:** Wipe down the seat and frame with a damp cloth and mild detergent. Do not use abrasive cleaners or solvents. Ensure the unit is dry before reconnecting power.
- **Inspection:** Regularly inspect the power cord for any signs of damage. Check all connections and moving parts for wear or looseness.
- **Storage:** When not in use for extended periods, disconnect the power adapter and store the unit in a dry, cool place. The unit can be folded for compact storage.

7. TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Unit does not power on.	Power cord not connected. Outlet not functioning.	Ensure power cord is securely plugged into the unit and a working outlet. Test the outlet with another device.
Unit does not move when buttons are pressed.	Remote control not connected. Overload protection activated.	Check remote control connection. Ensure the user's weight does not exceed the maximum capacity. Disconnect power for 5 minutes and try again.
Unusual noise during operation.	Loose component. Obstruction.	Inspect the unit for any loose parts or obstructions. If noise persists, discontinue use and contact support.

For issues not listed here, please contact MAIDeSITE customer support.

8. SPECIFICATIONS

Feature	Detail
Brand	MAIDeSITE
Model Number	ce0778a9-fed6-4418-ad3d-35dc909f1231
Product Dimensions (Folded)	28"D x 15"W x 3.9"H
Color	Blue
Material	Plastic (Frame), Foam (Seat)
Special Feature	Arm Rest
Unit Count	1.0 Count
Pattern	Solid
Arm Style	With arms
Surface Recommendation	Hard Floor
Indoor/Outdoor Usage	Indoor
Form Factor	Recliner
Is Foldable	Yes

9. WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the documentation included with your purchase or visit the official

MAIDeSITE website. If you encounter any issues or have questions regarding your Electric Chair Lift - Base, please contact MAIDeSITE customer service for assistance.

© 2024 MAIDeSITE. All rights reserved.