

CLEARCHOICE PRODUCTS 85db

CLEARCHOICE PRODUCTS Oticon Compatible Sport Lock 85db Hearing Aid Retention Tail User Manual

Brand: CLEARCHOICE PRODUCTS | **Model:** 85db Sport Lock

1. INTRODUCTION

This manual provides essential instructions for the proper use, installation, and maintenance of your CLEARCHOICE PRODUCTS Oticon Compatible Sport Lock 85db Hearing Aid Retention Tails. Designed for secure and comfortable wear, these retention tails help keep your Oticon hearing aids in place during various activities. Please read this manual thoroughly before use.

2. PRODUCT OVERVIEW

The Oticon Compatible Sport Lock is a retention accessory specifically designed for Oticon Hearing Aids utilizing 85 dB receivers. It provides an additional layer of security to prevent accidental dislodgement of your hearing aids.

Key Features:

- **Custom Fit for Oticon 85 dB Hearing Aids:** Meticulously designed for seamless compatibility and optimal performance.
- **Secure Retention:** Provides reliable retention to keep hearing aids in place during physical activities.
- **Hypoallergenic Material:** Crafted from a proprietary material gentle on sensitive skin.
- **Seamless Integration:** Engineered to integrate smoothly with your device without loose wires or awkward adjustments.
- **Multipack Convenience:** Each purchase includes a pack of 4 sport locks for replacements or backups.



Figure 2.1: Two individual Oticon Compatible Sport Locks. These clear, flexible retention tails are designed to provide a secure fit for hearing aids.

Key Features



Reliable
Performance



Secure
Positioning



Seamless
Integration



Perfect Fit for
Optimal Performance



Figure 2.2: Visual representation of the key benefits of the sport lock, including reliable performance and secure positioning.

3. INSTALLATION AND SETUP

Follow these steps to properly attach the sport lock to your Oticon hearing aid.

1. **Identify the Receiver:** Locate the receiver unit of your Oticon hearing aid. The sport lock is designed to attach to this component.
2. **Attach the Sport Lock:** Gently slide the small loop of the sport lock onto the receiver unit. Ensure it is securely seated and does not obstruct any part of the hearing aid.
3. **Positioning:** Once attached, the flexible tail of the sport lock should curve naturally along the contour of your ear, providing gentle pressure to keep the hearing aid in place.



Figure 3.1: An Oticon Compatible Sport Lock shown alongside a hearing aid receiver, demonstrating the attachment point.



Secure Hold For Your Hearing Aid's

These sport locks are engineered to be compatible with Oticon Hearing Aids, ensuring seamless integration and reliable performance.

Figure 3.2: A user demonstrating the insertion of a hearing aid with the sport lock into the ear, ensuring a secure fit.

Hypoallergenic Materials

Made from a proprietary hypoallergenic material



Quality Assurance



Easy Maintenance



Figure 3.3: The sport lock comfortably positioned in the ear, showcasing its discreet and hypoallergenic design.

4. USAGE AND WEARING

Once installed, the sport lock provides enhanced stability for your hearing aids. Ensure the tail rests comfortably against the concha of your ear.

- The sport lock is designed to be discreet and should not cause discomfort when properly fitted.
- Periodically check the positioning of the sport lock throughout the day, especially after strenuous activity, to ensure it remains secure.

Support Your Hearing Aid



These sport locks are engineered to be compatible with Oticon Hearing Aids.

These retention tails are crafted to keep your hearing aid's securely in place.

Ensuring a perfect fit for optimal performance.

These sport locks are gentle on the skin, suitable for users with allergies or sensitive skin.



Figure 4.1: A single sport lock held by hand, illustrating its form and flexibility for ear support.

5. MAINTENANCE

To ensure the longevity and hygiene of your sport locks, follow these maintenance guidelines:

- **Cleaning:** Gently wipe the sport locks with a soft, damp cloth. Avoid harsh chemicals or abrasive materials.
- **Inspection:** Regularly inspect the sport locks for any signs of wear, tear, or damage. Replace them if any damage is observed to maintain optimal retention.
- **Storage:** Store unused sport locks in a clean, dry place away from direct sunlight and extreme temperatures.

6. TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Sport lock feels loose or falls out.	Incorrect installation or wear, or the sport lock is worn out.	Re-check installation as per Section 3. If still loose, replace with a new sport lock from the multipack.
Discomfort or irritation in the ear.	Improper positioning or individual sensitivity.	Adjust the sport lock's position. If discomfort persists, discontinue use and consult a hearing care professional.

7. SPECIFICATIONS

- **Compatibility:** Oticon Hearing Aids with 85 dB receivers.
- **Material:** Proprietary Hypoallergenic Material.
- **Quantity:** Pack of 4 retention tails.
- **Item Weight:** 14.4 ounces (for the package).

8. ADDITIONAL RESOURCES

Instructional Video:

Your browser does not support the video tag.

Video 8.1: This video demonstrates the general concept of hearing aid ear grips, which can be helpful for understanding how retention tails function.

9. WARRANTY AND SUPPORT

CLEARCHOICE PRODUCTS stands by the quality of its products. For specific warranty information or technical support, please refer to the product packaging or contact CLEARCHOICE PRODUCTS customer service directly. Contact details can typically be found on the manufacturer's website or through your point of purchase.

© 2026 CLEARCHOICE PRODUCTS. All rights reserved.

This manual is for informational purposes only. Always consult with a qualified healthcare professional for any medical concerns or before making any decisions about your health or treatment.