

Eumspo JPD-HA200

Eumspo Automatic Blood Pressure Monitor User Manual

Model: JPD-HA200

1. INTRODUCTION

Thank you for choosing the Eumspo Automatic Blood Pressure Monitor. This device is designed for accurate and convenient measurement of blood pressure and pulse rate at home. It features a large LCD display, 2x99 memory storage for two users, and an arrhythmia detection function.

Please read this manual thoroughly before using the device to ensure correct operation and to obtain the most accurate results. Keep this manual for future reference.

2. IMPORTANT SAFETY INFORMATION

- Consult your physician before using this device if you have any medical conditions.
- Do not use this device for diagnosis or treatment of any health problem or disease. Always consult a medical professional for diagnosis and treatment.
- Do not adjust medication based on measurement results from this blood pressure monitor. Always follow the advice of your physician.
- This device is intended for adult use only. Keep out of reach of children.
- Avoid strong electromagnetic fields (e.g., mobile phones) near the device during measurement.
- Do not use the device in a moving vehicle.
- Ensure the cuff size is appropriate for your arm circumference (22-42cm / 8.66-16.53 inches).

3. PACKAGE CONTENTS

Upon opening the package, please verify that all items listed below are present and undamaged:

- 1x Eumspo Automatic Arm Blood Pressure Monitor
- 1x Upper Arm Cuff (22-42cm / 8.66-16.53 inches)
- 1x Storage Bag
- 1x Instruction Manual
- 1x Quick Guide
- 4x AAA Batteries (non-rechargeable)



Figure 3.1: Contents of the Eumspo Blood Pressure Monitor package. The image displays the blood pressure monitor, the large arm cuff, a black storage pouch, the instruction manual, and four AAA batteries.

4. PRODUCT OVERVIEW AND SETUP

4.1 Device Components



Figure 4.1: The Eumspo Blood Pressure Monitor and its accompanying arm cuff. The monitor features a large LCD display showing SYS, DIA, PUL, and time, along with control buttons. The cuff is designed for upper arm measurements.

The blood pressure monitor consists of the main unit with an LCD display and control buttons, and an upper arm cuff connected via an air tube.

4.2 Power Supply

The monitor supports two power supply modes:

- **Battery Power:** Insert the included 4 AAA batteries into the battery compartment located on the back of the device. Ensure correct polarity (+/-).
- **DC 5V Power:** The device can also be powered using a Micro USB B cable (not included) plugged into the DC 5V charging port. This is an alternative to battery power.

Note: The included AAA batteries are not rechargeable. Do not attempt to charge them.

Dual Power Supply Mode

4AAA battery (included, not rechargeable)+DC 5V mode
(charger not included) Suitable for use in a variety of scenarios



Figure 4.2: Illustration of the dual power supply options. Method 1 shows the monitor connected via a Micro USB B cable. Method 2 shows the battery compartment on the back of the monitor with four AAA batteries inserted.

4.3 Initial Setup (Time and User Settings)

Before first use, or after replacing batteries, you may need to set the date and time. Refer to the quick guide or the detailed instruction manual for specific steps on setting the year, month, and time. This ensures accurate recording of measurement nodes.

The device supports 2 users (User A and User B), each with 99 memory slots. This allows for individual tracking of readings. Instructions for selecting users and managing memory are detailed in the full manual.

5. OPERATING INSTRUCTIONS

5.1 Measurement Preparation

For accurate readings, follow these guidelines:

- Rest for at least 5 minutes before measurement.
- Avoid eating, drinking alcohol, smoking, or exercising for 30 minutes prior.
- Sit upright with your back supported, feet flat on the floor, and legs uncrossed.
- Place your arm on a flat surface so the cuff is at heart level.
- Ensure the room temperature is comfortable.

5.2 Applying the Cuff

1. Take out the cuff and wear it on your left upper arm. Ensure the air tube is securely inserted into the main device.
2. The cuff should wrap around your arm completely.
3. Position the cuff 2-3 cm (approximately 1 inch) above your elbow joint.
4. Ensure there is enough space for one or two fingers between the cuff and your arm. The cuff should be snug but not too tight.
5. The air tube should run along the inside of your arm.


BLOOD PRESSURE MEASUREMENT METHOD

1. Applying the cuff on the arm

- 1.1 Make sure the air plug is securely inserted in the main device.
- 1.2 Remove tight-fitting clothing from your upper left arm.
- 1.3 Push the cuff into a tubular shape, until it is suitable for the arm to get in.



- 1.4 Wrap the cuff around the upper arm, and also make the air tube located on the forearm.
- 1.5 The bottom of the cuff shall be located at the inner side of the upper arm elbow joint, going upward for 2cm~3cm.
- 1.6 Tightly wrap the cuff, and fix it with cloth buckle. Leave a finger gap between the arm and the cuff.
- 1.7 Align the air tube plug of the cuff with the air tube interface on the left side of the device body, and insert it to the end.



BLOOD PRESSURE MEASUREMENT METHOD

- User's position in normal use, including
 - 1) comfortably seated
 - 2) legs uncrossed
 - 3) feet flat on the floor
 - 4) back and arm supported
 - 5) middle the cuff at the level of the right atrium of the heart.
- a recommendation that the user relax as much as possible and do not talk during the measurement procedure;
- a recommendation that 5 min should elapse before the first reading is taken;
- a recommendation that measuring at same arm each time.



WRONG MEASUREMENT POSTURE

The cuff is wrapped too loosely



Move your body or speak during the measurement



leaning forward



⊗ If the cuff is positioned below the heart, place a cushion or towel underneath.

MÉTHODE DE MESURE DE LA TENSION ARTÉRIELLE

1. Application du brassard sur le bras

- 1.1 Assurez-vous que la prise d'air est bien insérée dans l'appareil principal.
- 1.2 Retirez les vêtements serrés de la partie supérieure de votre bras gauche.
- 1.3 Poussez le brassard en forme de tube jusqu'à ce que le bras puisse y entrer.



- 1.4 Enroulez le brassard autour de la partie supérieure du bras et placez le tube d'air sur l'avant-bras.
- 1.5 Le bas du brassard doit être situé à l'intérieur de l'articulation du coude du bras, en remontant de 2cm~3cm.
- 1.6 Enveloppez fermement le brassard et fixez-le à l'aide d'une boucle en tissu. Laissez un espace d'un doigt entre le bras et le brassard.
- 1.7 Alignez la prise du tube d'air du brassard avec l'interface du tube d'air sur le côté gauche du corps de l'appareil, et insérez-la jusqu'à l'extrémité.



MÉTHODE DE MESURE DE LA TENSION ARTÉRIELLE

- Position de l'utilisateur dans des conditions normales d'utilisation, y compris
 - 1) confortablement assis
 - 2) les jambes ne sont pas croisées
 - 3) les pieds à plat sur le sol
 - 4) le dos et le bras soutenus
 - 5) centrer le brassard au niveau de l'oreille droite du cœur.
- une recommandation selon laquelle l'utilisateur doit se détendre autant que possible et ne pas parler pendant la procédure de mesure ;
- une recommandation selon laquelle 5 minutes doivent s'écouler avant la première lecture ;
- une recommandation selon laquelle la mesure doit être effectuée sur le même bras à chaque fois.



MAUVAISE POSITION DE MESURE

La manchette est enroulée trop lâche



Bouger le corps ou parler pendant la mesure



se pencher en avant



⊗ Si le brassard est placé sous le cœur, placer un coussin ou une serviette sous le brassard.

Figure 5.1: Detailed instructions for cuff application and measurement posture. Step 1 shows wearing the cuff, Step 2 shows checking the 2-3cm distance from the elbow, Step 3 shows proper sitting posture with arm at heart level, and Step 4 shows pressing the measurement button.

5.3 Taking a Measurement

1. Ensure the cuff is correctly applied and you are in the proper measurement posture.
2. Press the 'START/STOP' button (usually the large central button) to begin the measurement.
3. The cuff will inflate, and the measurement will proceed automatically. Remain still and silent during this process.
4. The results (Systolic, Diastolic, and Pulse Rate) will be displayed on the LCD screen within approximately 45 seconds.
5. The device will automatically store the reading in the selected user's memory.

5.4 Understanding Your Readings

The display shows three main values:

- **SYS (Systolic Pressure):** The top number, indicating the pressure in your arteries when your heart beats.
- **DIA (Diastolic Pressure):** The bottom number, indicating the pressure in your arteries between beats.
- **PUL (Pulse Rate):** Your heart rate, measured in beats per minute.

The monitor also features an arrhythmia detection function, which will display a specific icon if an irregular heartbeat is detected during measurement. Consult your doctor if this icon appears frequently.

5.5 Memory Function

The device stores up to 99 readings for each of two users (User A and User B). To review past measurements:

- Press the 'Memory' button (often labeled 'M' or with a book icon) to cycle through stored readings.
- You can switch between User A and User B memory by pressing the user selection button (often with a person icon).

(2*99) Groups Of Stored Measurement Values, A More Than Two Sets Of Use




Figure 5.2: Representation of the dual-user memory storage. It shows two distinct memory banks, User A and User B, each capable of storing 99 measurement values (SBP, DBP, Pr) with corresponding dates, suitable for group data management.

6. MAINTENANCE AND STORAGE

- **Cleaning:** Wipe the monitor and cuff with a soft, dry cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in its storage bag in a cool, dry place, away from direct sunlight, extreme temperatures, and humidity.
- **Batteries:** Remove batteries if the device will not be used for an extended period to prevent leakage.
- **Cuff Care:** Do not fold the cuff too tightly or apply excessive force to the air tube.

7. TROUBLESHOOTING

Problem	Possible Cause	Solution
No display/Device won't turn on	Batteries are dead or incorrectly inserted; Power adapter not connected (if using DC 5V).	Replace batteries, ensuring correct polarity. Check DC 5V connection.
Inaccurate readings	Incorrect cuff placement; Movement or talking during measurement; Insufficient rest before measurement; Cuff size incorrect.	Re-apply cuff correctly. Remain still and silent. Rest for 5 minutes. Ensure cuff fits arm circumference (22-42cm).
Cuff does not inflate or inflates too little	Air tube disconnected or kinked; Cuff not properly sealed; Air leakage.	Check air tube connection to monitor and cuff. Ensure cuff is wrapped snugly. Inspect cuff for damage.
Error message on display	Various internal errors or measurement issues.	Refer to the full instruction manual for specific error codes. Often, turning the device off and on, or re-taking the measurement after resting, resolves the issue.

8. SPECIFICATIONS

Feature	Detail
Model Number	JPD-HA200
Measurement Method	Oscillometric
Cuff Size	22-42 cm (8.66-16.53 inches)
Memory Storage	2 users x 99 sets of readings
Power Supply	4x AAA Batteries (included) or DC 5V Micro USB B (cable not included)
Arrhythmia Detection	Yes
Measurement Time	Approx. 45 seconds
Manufacturer	LIKELINKS TRADING LIMITED
Parcel Dimensions	16.1 x 12.2 x 8.4 cm; 430 g

9. WARRANTY AND SUPPORT

Eumspo products are manufactured with care and attention to quality. For warranty information, technical support, or service inquiries, please refer to the contact details provided in the quick guide or on the product packaging. Keep your purchase receipt as proof of purchase for warranty claims.