

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

> [Semme](#) /

> [Semme Thumb Wearable LED Ear Pick Lamp Instruction Manual](#)

Semme SemmM8UhBKtxge-11

Semme Thumb Wearable LED Ear Pick Lamp Instruction Manual

Model: SemmM8UhBKtxge-11

1. PRODUCT OVERVIEW

The Semme Thumb Wearable LED Ear Pick Lamp is designed to assist in ear cleaning by providing illumination and a comfortable, hands-free experience. Its compact and lightweight design makes it convenient for personal use.

Key Features:

- **Thumb Wearable Design:** Frees up hands for easier ear picking.
- **Adjustable Brightness:** Features three white light brightness levels for optimal visibility.
- **Adjustable Illumination Angle:** Lamp beads can be adjusted to direct light into the ear canal.
- **USB Rechargeable:** Built-in lithium-ion battery for convenient charging.
- **Lightweight and Comfortable:** Designed for extended wear without discomfort.



Image 1: Semme Thumb Wearable LED Ear Pick Lamp.

2. PACKAGE CONTENTS

- 1 x LED Ear Pick Lamp
- 1 x Charging Cable

3. SETUP

Initial Charge:

Before first use, ensure the device is fully charged. Connect the provided USB charging cable to the lamp's charging port and a suitable USB power source. The indicator light will show charging status.



Image 2: Ear Pick Lamp with USB Charging Cable.

Wearing the Device:

Place the product on your thumb. The finger circle is designed to be non-scratchy and adjustable in size for a comfortable fit.



Image 3: Semme Ear Pick Lamp with adjustable thumb ring.

4. OPERATING INSTRUCTIONS

Turning On/Off and Adjusting Brightness:

1. Press the button once to turn on the light at the lowest brightness setting.
2. Press the button again to cycle through the three available brightness levels (low, medium, strong).
3. Press the button a fourth time to turn off the light.

The concentrated white light ensures clear visibility, making ear cleaning easier and safer.



Image 4: Ear Pick Lamp with light activated.

Adjusting Illumination Angle:

The lamp's illumination angle can be adjusted as needed to direct light precisely into the ear canal, enhancing visibility during use.



Image 5: Ear Pick Lamp showing flexible light arm.

5. MAINTENANCE

Cleaning:

After each use, gently wipe the ear pick lamp with a soft, dry cloth. Avoid using harsh chemicals or abrasive materials that could damage the device.

Charging:

Recharge the device when the light dims or fails to turn on. The built-in lithium-ion battery ensures long-lasting use. Always use the provided charging cable.

6. TROUBLESHOOTING

Light Not Turning On:

- Ensure the device is fully charged. Connect it to the charging cable and allow sufficient time to charge.
- Verify that the power button is being pressed correctly.

Light Dims Quickly:

- The battery may be low. Recharge the device.
- Ensure the charging cable and power source are functioning correctly.

7. SPECIFICATIONS

Feature	Detail
Item Type	LED Ear Pick Lamp
Material	Plastic, Alloy
Brightness	Three Speeds (White Light)
Built-in Battery	Lithium Ion Battery, 150mAh
Charging	USB Charging
Package Dimensions	5.39 x 3.94 x 0.91 inches
Item Weight	0.986 ounces
Model Number	SemmM8UhBKtxge-11

8. WARRANTY AND SUPPORT

For warranty information or customer support, please refer to the product packaging or contact Semme customer service directly. Contact details are typically provided with your purchase or on the official brand website.