

[Manuals.plus](#) /

> [Silhouette](#) /

> Silhouette TITAN MINIMAL ART 5599/NL Eyewear Frame User Manual

## Silhouette TITAN MINIMAL ART 5599/NL

# Silhouette TITAN MINIMAL ART 5599/NL Eyewear Frame User Manual

Model: 5599/NL

## PRODUCT OVERVIEW

---



The Silhouette TITAN MINIMAL ART 5599/NL eyewear frame is a lightweight, rimless design crafted for comfort and durability. It features ruthenium-colored titanium temples and clear lenses, embodying a minimalist aesthetic. This unisex frame is designed for prescription lenses and offers a sophisticated, understated look.

## SETUP

---

Your Silhouette TITAN MINIMAL ART 5599/NL frames require professional fitting and lens installation. Please follow these steps:

1. **Professional Lens Installation:** Take your frames to a qualified optician for the installation of your prescription

lenses. The optician will ensure proper lens shape and secure mounting.

2. **Initial Adjustment:** Have your optician adjust the nose pads and temple arms to ensure a comfortable and secure fit on your face. Proper adjustment prevents slippage and discomfort.
3. **Inspection:** Before first use, visually inspect the frames for any signs of damage or loose components.

## OPERATING (WEARING)

---

Proper handling and wearing of your eyewear frames will extend their lifespan and maintain comfort.

- **Wearing:** Always use both hands to put on and take off your frames. Grasp the temples near the hinges to avoid bending or stressing the frame.
- **Handling:** Avoid placing the frames face down on hard surfaces to prevent scratching the lenses. When not in use, store them in their protective case.
- **Adjustments:** If minor adjustments are needed (e.g., nose pad position), consult your optician. Attempting to adjust the frames yourself may cause damage.

## MAINTENANCE

---

Regular cleaning and care are essential for maintaining the clarity of your lenses and the integrity of your frames.

- **Cleaning Lenses:** Rinse lenses with lukewarm water to remove dust and debris. Apply a small drop of mild dish soap to each lens and gently rub with your fingertips. Rinse thoroughly and dry with a clean, soft microfiber cloth. Avoid paper towels or rough fabrics, which can scratch lenses.
- **Cleaning Frames:** Wipe the frame, including nose pads and temples, with a damp microfiber cloth. For deeper cleaning, use a mild soap solution and rinse well.
- **Storage:** Store your frames in the provided protective case when not in use to prevent damage from drops, pressure, or scratches.
- **Avoid Extreme Conditions:** Do not expose your frames to extreme heat (e.g., dashboard of a car in direct sunlight) or harsh chemicals, as this can damage the frame material and lens coatings.

## TROUBLESHOOTING

---

Here are solutions to common issues you might encounter:

- **Frames feel loose or slide down:** This indicates the frames need adjustment. Visit your optician for professional tightening and reshaping of the temple arms and nose pads.
- **Discomfort behind ears or on nose:** The frames may be too tight or improperly adjusted. Consult your optician for a comfort fitting.
- **Lenses appear smudged or dirty:** Follow the cleaning instructions in the Maintenance section. Ensure you are using a clean microfiber cloth and mild soap.
- **Small screws are loose:** Do not attempt to tighten small screws yourself unless you have the correct eyewear repair kit. Visit your optician for professional tightening to avoid stripping the screws or damaging the frame.
- **Frame appears bent or misaligned:** Avoid attempting to bend the frame back into shape. Titanium is durable but can be permanently damaged if improperly handled. Seek professional repair from your optician.

## SPECIFICATIONS

---

<b>Brand</b>	Silhouette
<b>Model Number</b>	TITAN MINIMAL ART 5599/NL
<b>Material</b>	Titanium
<b>Color</b>	Ruthenium
<b>Size</b>	55/17/0 (Lens Width / Bridge Width / Temple Length)
<b>Department</b>	Unisex-Adult
<b>Product Dimensions</b>	5.12 x 1.18 x 3.15 inches
<b>Item Weight</b>	2.82 ounces
<b>Manufacturer</b>	Silhouette

## WARRANTY AND SUPPORT

---

Your Silhouette TITAN MINIMAL ART 5599/NL eyewear frames come with a **2 Years International Warranty**. For warranty claims, product support, or further assistance, please contact your authorized Silhouette dealer or visit the official Silhouette website. Keep your proof of purchase for warranty validation.

**Manufacturer Contact:** Silhouette International Schmied AG

**Website:** [www.silhouette.com](http://www.silhouette.com)