

[Manuals.plus](#) /

› [HOMLAND](#) /

› HOMLAND Shower Chair User Manual (Model SHC10-V11)

HOMLAND SHC10-V11

HOMLAND Shower Chair User Manual

MODEL: SHC10-V11

Brand: HOMLAND

1. INTRODUCTION

Thank you for choosing the HOMLAND Shower Chair, Model SHC10-V11. This heavy-duty bath seat is designed to provide safety, stability, and comfort for individuals requiring assistance during bathing, including seniors, disabled persons, and those recovering from surgery. With a robust aluminum alloy frame and thoughtful design features, this shower chair enhances the bathing experience. Please read this manual thoroughly before assembly and use to ensure proper operation and safety.



Figure 1: HOMLAND Shower Chair, Model SHC10-V11.

2. SAFETY INFORMATION

Your safety is paramount. Adhere to the following guidelines:

- **Weight Capacity:** The shower chair is designed to support a maximum weight of 400 lbs. The recommended user weight is 300 lbs. Do not exceed these limits.
- **Stability:** Always ensure all four legs are evenly placed on a flat, non-slip surface within the shower or bathtub. Verify the non-slip rubber feet are securely attached and free from damage before each use.
- **Assembly:** Follow all assembly instructions carefully. Ensure all screws and fasteners are tightened securely. Do not use the chair if any part is loose or damaged.
- **Adjustments:** When adjusting the height, ensure the locking pins are fully engaged in the adjustment holes on all legs. Test stability after any adjustments.
- **Water Drainage:** The seat features drain holes to prevent water accumulation. Ensure these holes remain clear.
- **Supervision:** Individuals with severe mobility limitations should be supervised during use.
- **Cleaning:** Clean the chair regularly to prevent soap scum buildup, which can reduce the effectiveness of non-slip surfaces.

3. PRODUCT FEATURES

The HOMLAND Shower Chair offers several features designed for user convenience and safety:

- **Durable Construction:** Made with reinforced aluminum alloy frames, providing strength and corrosion resistance.
- **High Weight Capacity:** Supports up to 400 lbs (recommended user weight 300 lbs).
- **Non-Slip Rubber Feet:** Wide, skid-resistant feet with stainless steel washers ensure secure attachment to wet surfaces and prevent rust stains.
- **Adjustable Height:** Seat height can be adjusted from 14 inches to 19 inches in 1-inch increments, accommodating various user heights and bathroom setups.
- **Removable Backrest:** The backrest can be easily attached or removed, allowing the chair to be used as a shower seat or a backless stool.
- **Integrated Drainage:** Built-in handle and seat holes facilitate water drainage, preventing pooling and enhancing safety.
- **Tool-Free Assembly:** Designed for quick and easy setup without the need for additional tools.
- **Showerhead Clip:** Includes a convenient clip holder for easy placement of a showerhead.

Ours	Other
Safety Design 400lb   28.6mm 1.35mm ✓ Durable Frames (Φ28.6mm x 1.35mm)  ✓ Soft Non-Slip Pad  ✓ Shower Hose Clip	Not Secure 300lb   25mm 0.8mm ✗ Slim Frames (Φ25mm x 0.8mm)  ✗ Rough Panel  ✗ Without Hose Clip

Figure 2: Key safety and design features of the HOMLAND Shower Chair.

4. SETUP AND ASSEMBLY

The HOMLAND Shower Chair is designed for tool-free assembly. Follow these steps:

1. **Unpack Components:** Carefully remove all parts from the packaging. Verify that all components are present according to the packing list.
2. **Attach Frame to Seat:** Align the curved aluminum frame pieces with the underside of the seat. Insert the screws and hand-tighten the blue knobs to secure the frame to the seat.
3. **Insert Legs:** Insert the four adjustable leg extensions into the frame. Ensure the spring-loaded pins snap securely into the desired height adjustment holes.
4. **Attach Backrest (Optional):** If using the backrest, slide its support bars into the designated slots on the rear of the seat frame. Secure with the provided blue knobs.

5. **Attach Non-Slip Pad:** Peel the backing from the blue anti-slip pad and firmly press it onto the seat surface, aligning the drain holes.
6. **Attach Showerhead Clip:** Snap the showerhead clip onto one of the chair legs at a convenient height.
7. **Final Check:** Before first use, ensure all connections are tight and the chair is stable.



Figure 3: All components of the HOMLAND Shower Chair laid out for assembly.

Your browser does not support the video tag.

Video 1: Official HOMLAND Shower Chair assembly and feature overview. This video demonstrates the simple assembly process and highlights key features of the shower chair.

5. OPERATING INSTRUCTIONS

Using your HOMLAND Shower Chair safely and effectively:

- **Placement:** Position the shower chair in your bathtub or shower stall. Ensure all four non-slip feet are firmly on the surface and the chair does not wobble.
- **Height Adjustment:** To adjust the seat height, press the spring-loaded pin on each leg and slide the leg extension up or down to the desired height. Ensure the pin locks into place in one of the six adjustment holes. Adjust all four legs to the same height for stability. The seat height can range from 14 to 19 inches.
- **Using the Backrest:** The backrest provides additional support. If you prefer a backless stool, simply remove the backrest by loosening its knobs and sliding it out.
- **Showerhead Holder:** Utilize the integrated showerhead clip to keep your handheld shower conveniently within reach.
- **Entering/Exiting:** Always use caution when sitting down on or standing up from the chair. Use grab bars if available for added support.

14" - 19"

⑥ Levels of Adjustable Height



Figure 4: The shower chair offers 6 levels of adjustable height, from 14 to 19 inches.

Removable Backrest



Figure 5: The backrest can be easily removed to convert the shower chair into a shower stool.

6. MAINTENANCE

Proper maintenance ensures the longevity and safety of your shower chair:

- **Cleaning:** After each use, rinse the chair with clean water to remove soap residue. Wipe down with a soft cloth and mild detergent as needed. Avoid abrasive cleaners that could damage the surface.
- **Drying:** Allow the chair to air dry completely or wipe it dry to prevent mildew and water spots.
- **Inspection:** Regularly inspect all components, especially the leg extensions, locking pins, and rubber feet, for any signs of wear, cracks, or damage. Ensure all screws and knobs remain tight.
- **Storage:** If not in regular use, store the chair in a dry, clean area.

7. TROUBLESHOOTING

If you encounter any issues with your HOMLAND Shower Chair, refer to the following:

- **Chair Feels Unstable/Wobbly:**

- Ensure all four legs are adjusted to the same height and the locking pins are fully engaged.
- Check that the non-slip rubber feet are clean and making full contact with the shower/tub surface.
- Verify all assembly screws and knobs are securely tightened.

- **Water Pooling on Seat:**

- Ensure the drain holes on the seat and the anti-slip pad are clear of any debris or soap residue.

- **Difficulty Adjusting Height:**

- Ensure the spring-loaded pin is fully depressed before attempting to slide the leg extension.
- Check for any obstructions or damage to the leg extensions or adjustment holes.

If problems persist, please contact HOMLAND customer support.

8. SPECIFICATIONS

Feature	Detail
Model Number	SHC10-V11
Max Weight Capacity	400 lbs (Recommended User Weight: 300 lbs)
Overall Dimensions (D x W x H)	16"D x 16.3"W x 24.6"-29.5"H
Seat Dimensions (W x D x H)	16.3"W x 10.5"D x 14"-19"H
Minimum Height	13 Inches
Seat Depth	16 inches
Seat Width	17 Inches
Item Weight	5.07 Pounds
Materials of Frame	Aluminum Alloy
Materials of Seat	HDPE
Color	White
Application Scenarios	Bathroom and Bathtub

Ideal for Small Spaces

Perfect for bathtubs and showers



Figure 6: Detailed dimensions of the HOMLAND Shower Chair.

9. WARRANTY AND SUPPORT

HOMLAND is committed to your satisfaction. The HOMLAND Shower Chair comes with a 1-year warranty plus a 1-year extension service.

For any questions, concerns, or technical support, please contact our USA-based customer support team. We offer lifetime free technical support to ensure a worry-free experience.

Contact Information: Please refer to the contact details provided on your product packaging or the official HOMLAND website for the most up-to-date support information.

