

[Manuals.plus](#) /

> [Wokyo](#) /

> Wokyo P13 Hearing Aid Battery Charger Station User Manual

Wokyo Charger + 2pcs 13 Batteries

Wokyo P13 Hearing Aid Battery Charger Station User Manual

Model: Charger + 2pcs 13 Batteries

1. INTRODUCTION

Thank you for choosing the Wokyo P13 Hearing Aid Battery Charger Station. This compact and convenient charger is designed to provide reliable power for your size 13 rechargeable hearing aid batteries, offering an environmentally friendly and economical alternative to disposable batteries. This manual provides essential information for the proper setup, operation, maintenance, and troubleshooting of your device.



Figure 1: Wokyo P13 Hearing Aid Battery Charger Station with two size 13 rechargeable batteries.

2. SAFETY INFORMATION

- Keep the charger and batteries away from water and moisture.
- Do not attempt to disassemble or modify the charger or batteries.
- Use only the specified Type-C charging cable and a compatible power source (5V).
- Keep out of reach of children and pets.
- Dispose of old batteries according to local regulations.
- Avoid extreme temperatures.

3. PACKAGE CONTENTS

-
- 1 x Wokyo P13 Hearing Aid Battery Charger Station
 - 2 x 28maH Rechargeable Silver-Zinc Batteries (Size 13)
 - 1 x USB Type-C Charging Cable

4. PRODUCT OVERVIEW

The Wokyo P13 charger station is designed for ease of use and portability. It features two charging slots for size 13 batteries and indicator lights to show charging status.



Figure 2: Key features of the Wokyo P13 Rechargeable Hearing Aid Batteries and Charger.

5. SETUP AND INITIAL CHARGE

5.1 Connecting the Charger

To begin, connect the provided USB Type-C charging cable to the charging port on the back of the charger station. Plug the other end of the USB cable into a compatible 5V power source (e.g., a USB wall adapter, computer USB port, or power bank).



Figure 3: Location of the Type-C charging port and indicator light states.

5.2 Inserting Batteries

Carefully place the size 13 rechargeable batteries into the charging slots. Ensure the positive pole of the battery (marked with a '+' or flat side) faces the front of the charger, and the negative pole faces the rear. The two indicator lights on the front of the charging box will illuminate.

6. OPERATION

6.1 Charging Process

When the batteries are correctly inserted and the charger is connected to power, the indicator light for each battery will turn **red**, indicating that the battery is charging. A full charge typically takes approximately 6 hours. After 8 hours of full

charge, the indicator lights will become **green**, signifying that the batteries are fully charged and ready for use.

Easily

Each charging case comes with **2** rechargeable batteries, so you don't have to buy disposable batteries



Environmentally

The use of rechargeable hearing aid batteries can reduce 50-200 discarded disposable batteries a year and reduce environmental pollution



Economically



VS



One rechargeable battery can be used for **300-500** times=**2-3** years disposable battery usage

Figure 4: Battery insertion and indicator light status.

6.2 Battery Life

Each 28mAh battery provides up to 25 hours of use for non-Bluetooth hearing aids on a 6-hour charge. For Bluetooth hearing aids, the battery life may be reduced to approximately 4 hours after a full charge.

7. MAINTENANCE AND CARE

7.1 Cleaning

Wipe the charger station and batteries with a soft, dry cloth. Do not use liquid cleaners or abrasive materials.

7.2 Storage

Store the charger and batteries in a cool, dry place when not in use. Avoid prolonged exposure to direct sunlight or extreme temperatures.

7.3 Environmental Benefits and Longevity

One Wokyo rechargeable battery can be used for 300-500 charging cycles, which is equivalent to 2-3 years of disposable battery usage. By using rechargeable batteries, you can significantly reduce environmental pollution by avoiding the disposal of 50-200 single-use batteries per year.


<h1>SPECIFICATION</h1> 
Charging Time: 8H
Battery last time per charge:
25H (Non-Bluetooth hearing aids)
8H (Bluetooth hearing aids)
Battery longevity: >500 times charging-discharging
input Voltage/Current: 5VMA
Battery Capacity: 28maH/pack
Matching Batteries: 2 x 13 batteries (comes with)
Material: ABS & Metal
Product Dimensions: 12cm*10cm*10cm/4.72in*3.94in*3.94in

Figure 5: Environmental and economic advantages of Wokyo rechargeable batteries.

8. TROUBLESHOOTING

8.1 Battery Not Charging / Indicator Light Issues

If the battery does not appear to be charging (indicator light remains green when it should be red, or no light), it may be due to poor contact between the battery and the charger interface. In this case, please rotate or gently adjust the battery within the charging slot until the indicator light turns **red**.

Your browser does not support the video tag.

Video 1: This video demonstrates how to properly insert batteries into the Wokyo P13 charger and troubleshoot common charging issues related to battery contact.

If the indicator light is **green** when you just put a partially charged battery on the charger, it means that the battery may be poorly exposed or not making proper contact. Re-adjust the battery.

8.2 Reduced Battery Life

If you experience significantly reduced battery life, ensure the batteries are fully charged (green indicator light). Battery performance can vary based on hearing aid model, usage patterns (e.g., Bluetooth streaming), and battery age. If issues persist, consider replacing the batteries after 300-500 cycles.

9. SPECIFICATIONS

Feature	Detail
Brand	Wokyo
Model Number	Charger + 2pcs 13 Batteries
Battery Capacity	28maH/pack
Matching Batteries	2 x Size 13 (included)
Battery Cell Composition	Silver-Zinc
Charging Time	Approx. 6-8 hours for full charge
Battery Life (Non-Bluetooth HA)	Up to 25 hours per charge
Battery Life (Bluetooth HA)	Approx. 4 hours per charge
Battery Longevity	>500 charging-discharging cycles
Input Voltage/Current	5V
Material	ABS & Metal
Product Dimensions	12cm x 10cm x 10cm (4.72in x 3.94in x 3.94in)
Item Weight	1.9 ounces
Recommended Use	Hearing Aid

Rechargeable 13 Hearing Aid Batteries and Charger



- Battery Capacity 28maH 
- Comes with 13 Battery 
- Durable and long lasting 

Battery last time per charge:
25H (Non-Bluetooth hearing aids)
8H (Bluetooth hearing aids)



Figure 6: Detailed technical specifications.

10. WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the product packaging or contact Wokyo customer service through the retailer's platform where the product was purchased. Please retain your proof of purchase for any warranty claims.

