



Manuals.plus /

› **Vive Mobility** /

› Vive Mobility Bariatric Transport Wheelchair (Model B0DBGSJYZL) User Manual

Vive Mobility B0DBGSJYZL

Vive Mobility Bariatric Transport Wheelchair

Model: B0DBGSJYZL

1. INTRODUCTION

This manual provides essential instructions for the safe and effective use, setup, operation, and maintenance of your Vive Mobility Bariatric Transport Wheelchair. Designed for adults, this heavy-duty, foldable transport chair supports individuals up to 450 pounds, offering portability and comfort with its 22-inch wide seat and swing-away footrests.

2. SAFETY INFORMATION

Please read all safety warnings and instructions carefully before using the transport wheelchair. Failure to do so may result in injury or damage to the product.

- **Doorway Clearance:** The chair features an extra-wide frame measuring 27 inches wide. Be aware that this increased width may not fit through narrow doorways in some homes or facilities. Always measure your entryways to ensure safe passage.
- **Weight Capacity:** Do not exceed the maximum weight capacity of 450 pounds.
- **Brakes:** Always engage the handbrakes when transferring in or out of the chair, or when the chair is stationary.
- **Supervision:** Ensure the user is properly seated and secured before transport. Always have a caregiver present to operate the chair.
- **Uneven Surfaces:** Exercise caution on uneven terrain, inclines, or declines. Avoid sudden movements.
- **Folding Mechanism:** Keep hands and fingers clear of moving parts when folding or unfolding the chair to prevent pinching.

3. SETUP

Your Vive Mobility Bariatric Transport Wheelchair arrives largely pre-assembled. Follow these steps to prepare it for use:

1. **Unfolding the Chair:** Carefully unfold the chair by pushing down on the seat rails until the chair is fully open and the seat is flat. Ensure all locking mechanisms are engaged.
2. **Attaching Footrests:** Insert the swing-away footrests into the designated slots on the front frame of the

wheelchair. Secure them by pushing down until they click into place. Adjust the height of the footrests as needed for user comfort.

Maximum Stability

Supports up to 450 lbs



Figure 3.1: The Vive Mobility Bariatric Transport Wheelchair with its key dimensions, including a 22-inch wide seat and overall frame width of 27 inches.

Setup Video Guide

Your browser does not support the video tag.

Video 3.1: An overview of the Vive Mobility Bariatric Transport Chair, demonstrating its features and initial setup.

4. OPERATING INSTRUCTIONS

4.1. Transporting a User

1. Ensure the user is comfortably seated and any safety belts are fastened.
2. The caregiver should stand behind the chair, gripping the ergonomic hand grips.
3. To move, gently push the chair forward. The 12-inch smooth-rolling wheels are designed for easy navigation on various indoor and outdoor surfaces.
4. To stop, engage the dual handbrakes located on the push handles.



**Item weight:
48 Pounds**

Supports up to 450 lbs
Steel frame provides maximum support

Extra-Wide Seat
22" wide seat & padded backrest

Arrives Pre-Assembled
Use it straight out of the box

Ideal for Heights 4'11" - 6'3"
Recommend a wheelchair cushion for taller users.

Figure 4.1: A caregiver pushing the bariatric transport wheelchair, highlighting its caregiver-friendly design and ergonomic hand grips.

4.2. Using Handbrakes

The transport chair is equipped with dual handbrakes for secure stopping and parking.

- To engage the brakes, squeeze the levers on the push handles.
- For parking, push the brake levers down until they lock into place. This activates the wheel locks.
- To release, pull the brake levers up.

Comfortable Seat & Backrest

With an extra-wide 22" seat



Figure 4.2: Detail of the dual handbrakes, providing enhanced safety and control for the caregiver.

4.3. Adjusting Footrests

The swing-away footrests can be adjusted for height and swung out of the way for easier transfers.

- To adjust height: Locate the adjustment pins on the footrest hangers. Press the pins and slide the footrest up or down to the desired height. Ensure both footrests are at the same height and securely locked.
- To swing away: Lift the release lever on the side of the footrest hanger and swing the footrest outwards.

Easy to Push & Maneuver



Figure 4.3: The wide seat and adjustable footrests provide comfort and support for the user.

4.4. Folding for Storage and Transport

The chair features a simple folding design for compact storage and easy transport.

1. Ensure the footrests are either removed or swung to the side.
2. Lift the center of the seat fabric upwards. The chair frame will begin to collapse inwards.
3. Continue lifting until the chair is fully folded.

Folds Quickly & Easy

For transport or compact storage



Figure 4.4: The transport wheelchair folds quickly and easily for compact storage or transport.

5. MAINTENANCE

Regular maintenance ensures the longevity and safe operation of your transport wheelchair.

- **Cleaning:** Wipe down the frame and upholstery regularly with a damp cloth and mild soap. Avoid abrasive cleaners.
- **Wheel Inspection:** Periodically check the wheels for wear and tear. Ensure they roll smoothly and are free of debris.
- **Brake Check:** Test the handbrakes regularly to ensure they engage and disengage properly.
- **Fasteners:** Check all bolts, screws, and fasteners periodically to ensure they are tight. Do not overtighten.

6. TROUBLESHOOTING

If you encounter any issues with your Vive Mobility Bariatric Transport Wheelchair, please refer to the following general advice:

- **Difficulty Rolling:** Check wheels for debris or damage. Ensure brakes are fully disengaged.
- **Brakes Not Engaging:** Inspect brake cables and mechanisms for proper alignment and function.

- **Unstable Feel:** Verify that the chair is fully unfolded and all locking mechanisms are engaged. Check for loose fasteners.

For persistent issues or concerns not addressed here, please contact Vive Mobility customer support.

7. SPECIFICATIONS

Feature	Detail
Model Number	B0DBGSJYZL
Product Dimensions (D x W x H)	42" x 27" x 36"
Item Weight	48 Pounds
Weight Capacity	450 Pounds
Seat Width	22 Inches
Seat Height	19.5 Inches
Seat Depth	16 Inches
Wheel Size	12 Inches (Rear)
Turn Radius	27 Inches
Frame Material	Alloy Steel
Brake Style	Hand brake
Handle Type	Push Handles
Power Source	Manual Propulsion
Wheel Type	Solid

8. WARRANTY INFORMATION

The Vive Mobility Bariatric Transport Wheelchair comes with a **1-Year Guarantee**. Please retain your proof of purchase for warranty claims. For specific terms and conditions, refer to the warranty documentation included with your product or contact Vive Mobility customer support.

Pre-Assembled Frame

Ideal for Heights 4'11" – 6'3"



Figure 8.1: The Vive Health 1 Year Guarantee ensures product quality and customer satisfaction.

9. CUSTOMER SUPPORT

For assistance, questions, or to report any issues, please contact Vive Health Direct customer support. Information regarding returns and customer service can be found on the seller's platform or through the contact details provided with your purchase.

Manufacturer: Vive Health