

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Golden Technologies](#) /
- › [Golden Tech GB120 Buzzaround Carry-On Travel Scooter Instruction Manual](#)

Golden Technologies GB120-DT

Golden Tech GB120 Buzzaround Carry-On Travel Scooter Instruction Manual

Model: GB120-DT

INTRODUCTION

The Golden Tech GB120 Buzzaround Carry-On is a compact, fold-flat travel scooter designed for portability and ease of use. This manual provides essential information for the safe operation, maintenance, and care of your scooter, ensuring optimal performance and longevity. Please read this manual thoroughly before operating the device.

This portable, stylish scooter ensures comfortable travel with a 300 lb weight capacity and an infinitely adjustable tiller equipped with a handy USB port and drink holder. It boasts full front and rear lighting and flip-up armrests.

SAFETY INFORMATION

Adhering to safety guidelines is crucial for your well-being and the proper functioning of the scooter. Always operate the GB120 Buzzaround Carry-On with caution.

- **MAXIMUM USER WEIGHT:** Do not exceed the maximum user weight of 300 lbs (136 kg). Exceeding this limit can compromise the scooter's stability and performance.
- **GROUND CLEARANCE:** The scooter has a ground clearance of 4 inches (10 cm). Avoid operating on uneven, rough, or excessively sloped terrain to prevent tipping or damage.
- **BATTERY SAFETY:** Always follow the manufacturer's instructions for battery charging, handling, and storage. Use only the approved charger. Do not expose the battery to extreme temperatures or puncture it.
- **OPERATE WITH CAUTION AT NIGHT:** Utilize the integrated front and rear lights for enhanced visibility when operating the scooter in low-light conditions or at night.
- **GENERAL OPERATION:** Always ensure the scooter is fully unfolded and locked into position before use. Avoid sudden stops or turns.

SETUP AND ASSEMBLY

The Buzzaround Carry-On is designed for quick setup. Follow these steps to prepare your scooter for use.

Unfolding the Scooter

The scooter folds to a compact size for transport. To unfold:

1. Place the folded scooter on a flat, stable surface.
2. Gently lift the tiller assembly until it locks into the upright position.
3. Unfold the seat frame and ensure it locks securely into place.
4. Attach the seat by aligning it with the seat post and pressing down until it clicks.
5. Adjust the tiller to your desired height and angle using the adjustment lever.



Image: The scooter in a partially unfolded state, illustrating the folding tiller and seat mechanism.



Image: The scooter fully folded and standing upright, demonstrating its space-saving design for storage or transport.

Battery Installation and Charging

The GB120 comes with a standard 12.5 AH lithium-ion battery pack. An optional 15 AH long-range battery is available separately.

- Ensure the battery pack is fully charged before first use.
- Locate the battery compartment, typically under the seat or deck.
- Insert the battery pack, ensuring it is securely connected.
- To charge, connect the charger to the charging port on the scooter or directly to the battery pack if removable. Plug the charger into a standard wall outlet.
- The charging indicator light will show the charging status. Disconnect once fully charged.



Image: The scooter disassembled into its primary components, including the battery pack, for easier transport or storage.

OPERATING INSTRUCTIONS

Familiarize yourself with the controls and operating procedures before driving your scooter.

Controls and Features

- **Tiller:** The steering column, which also houses the controls. It is infinitely adjustable for user comfort.
- **Speed Control Dial:** Adjusts the maximum speed of the scooter.
- **Forward/Reverse Lever:** Push forward to move forward, pull back to move in reverse.
- **Battery Indicator:** Displays the current battery charge level.
- **Headlights/Taillights Switch:** Activates the front and rear lighting system for visibility.
- **USB Port:** Located on the tiller for charging personal electronic devices.

- **Drink Holder:** Integrated into the tiller for convenience.
- **Flip-up Armrests:** Provide comfort and ease of access to and from the seat.



Image: Front view of the scooter, highlighting the tiller with controls, front basket, and LED headlights.



Image: Rear view of the scooter, showing the comfortable seat, rear storage basket, and taillights.

Driving the Scooter

1. Ensure the scooter is fully assembled and the battery is charged.
2. Sit comfortably on the seat with both feet on the footboard.
3. Turn the key to the "ON" position (if applicable) or press the power button.
4. Set the speed control dial to a comfortable speed, starting slow.
5. Gently push the forward/reverse lever to begin moving. Release the lever to engage the electromagnetic brake and stop.
6. Steer by turning the tiller in the desired direction.
7. Be mindful of the scooter's operating range, which is approximately 14 miles with the standard battery.



Image: The Golden Tech GB120 Buzzaround Carry-On scooter ready for operation, showcasing its design and features.

MAINTENANCE

Regular maintenance ensures the longevity and reliable performance of your Buzzaround Carry-On scooter.

- **Cleaning:** Wipe down the scooter regularly with a damp cloth. Avoid using abrasive cleaners or excessive water, especially near electrical components.
- **Battery Care:**
 - Charge the battery after each use, even if only for a short period.
 - Do not allow the battery to fully discharge.
 - If storing the scooter for an extended period, charge the battery fully and then charge it once a month to maintain its health.
- **Tire Inspection:** Periodically check the tires for wear and tear, and ensure they are free of debris.
- **Fasteners:** Check all nuts, bolts, and fasteners regularly to ensure they are tight. Do not overtighten.
- **Professional Service:** For any complex repairs or issues, contact an authorized Golden Technologies service center.

TROUBLESHOOTING

This section addresses common issues you might encounter with your scooter. For problems not listed here, contact customer support.

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Scooter does not turn on.	Battery not charged or not connected; power switch off; circuit breaker tripped.	Ensure battery is charged and properly connected. Turn power switch on. Check and reset circuit breaker if tripped.
Scooter moves slowly or intermittently.	Low battery charge; speed control set too low; motor/controller issue.	Charge the battery. Increase speed setting. If problem persists, contact service.
Brakes are not engaging properly.	Brake mechanism obstruction; electrical issue.	Check for obstructions around the brakes. Do not attempt to repair electrical issues; contact service.
Unusual noises during operation.	Loose components; motor issue; wheel bearing issue.	Inspect for loose parts. If noise persists or is severe, discontinue use and contact service.

SPECIFICATIONS

Key specifications and features of the Golden Tech GB120 Buzzaround Carry-On Travel Scooter.

Feature	Detail
Model Number	GB120-DT
Brand	Golden Technologies
Color	Desert Tan
Material	Aluminum, Vinyl
Maximum User Weight	300 lbs (136 kg)
Ground Clearance	4 inches (10 cm)
Folded Height (without seat)	12 inches (30.5 cm)
Operating Range (Standard Battery)	Up to 14 miles (22.5 km)
Battery Type	12.5 AH Lithium-Ion (Airline Approved)
Seat Depth	13.75 inches (35 cm)
Tiller Features	Infinitely adjustable, USB port, drink holder
Lighting	Full front and rear lights
Portability	Folds flat, four carry handles, easily disassembles

Additional Features:

- Ideal for space-saving storage.
- Folding tiller and self-locking seat frame for a compact size.
- Flip-up armrests for convenience.
- Optional long-range battery (15 ah, 18-mile range) available for purchase separately.
- Standard 12.5 AH lithium-ion battery pack with easy access.

- Airline-approved battery included with every GB120.



Standard Colors:



Velvet Red
(Installed)



Coastal Blue
(Included in Box)

NEW Sporty

Optional Colors:



Crisp White



Canary Yellow



Stone Grey



Orchid Pink



Army Green
(Matte)



Desert Tan

Buzzaround GB120



The Buzzaround CarryOn fold-flat scooter is small in stature and big on fun! This super portable, sporty scooter will take you where you want to go in comfort and style! It features a 300 lb. weight capacity, an infinite adjustable tiller with a convenient USB port and drink holder, full front and rear light package, and flip up armrests. Enjoy long range coverage up to 14 miles on the airline friendly 12.5 AH lithium-ion battery, which can be removed and taken as a carry on while traveling by air, sea or train in a stylish carry case. A new rear basket is mounted on the deck to hold an optional spare 12.5 AH battery for even longer range and serves as a kickstand for storing upright. The heaviest piece is 47.5 lbs. (without seat and batteries) and the height when folded is only 12". It disassembles easily into two pieces and is designed to store in the smallest compartments. It comes standard with red color shrouds installed and an extra set in blue. New optional color shroud kits let you express a style all your own! Pack your dreams and go anywhere with the CarryOn Scooter!



NEW & IMPROVED
Sporty Fold-Flat Scooter!
Airline Friendly!
USB Port Standard!

Standard Features:

- 4-point carry handles for easy transport
- 12" height when folded (without seat)
- 4" ground clearance
- 300 lb. weight capacity
- Flip up armrests
- Infinite adjustable tiller
- USB port & drink holder
- Disassembles easily
- Up to 14 mile range
- Easy access, airline friendly 12.5 AH lithium ion battery pack

Image: An informational graphic detailing the scooter's features and available color options, including Desert Tan.

WARRANTY AND SUPPORT

Golden Technologies provides a warranty for the GB120 Buzzaround Carry-On scooter. Specific warranty terms and conditions may vary. Please refer to the warranty card included with your product or visit the official Golden Technologies website for detailed information.

For technical support, spare parts, or warranty claims, please contact Golden Technologies customer service directly. It is recommended to have your model number (GB120-DT) and purchase date available when contacting support.

You can find contact information for Golden Technologies on their official website: www.goldentech.com

© 2024 Golden Technologies. All rights reserved.

This manual is for informational purposes only. Specifications are subject to change without notice.

